**Prospektüs**  
Prospectus

# FASTSORB

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GB

**FASTSORB**Surgical, Polyglycolic Acid  
Rapid Absorbable PGA Suture,  
Synthetic, Multifilament, Braided, Undyed

## Instruction for Use

**DESCRIPTION**

FASTSORB Surgical Degraded Polyglycolic acid (PGA) suture; synthetic, sterile, absorbable, is homopolymer of glycolic acid. The empirical formula of polyglycolic acid is (-O-CH<sub>2</sub>-CO-O-CH<sub>2</sub>-CO-). The polymer chain is degraded by gamma radiation. FASTSORB is coated with the mixture of calcium stearate (C<sub>18</sub>H<sub>36</sub>O<sub>2</sub>Ca) and polycaprotone (C<sub>12</sub>H<sub>24</sub>O<sub>4</sub>n). The suture and coating material are inert, noncollagenous, nonantigenic, and non-pyrogenic. FASTSORB suture is undyed in natural beige color.

**INDICATIONS**

FASTSORB suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues. FASTSORB suture is for single use only.

**ACTIONS**

PGA suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGA suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis it's excreted from the body as carbon dioxide and water. Absorption occurs as a loss of tensile strength followed by a loss of mass. PGA suture retains approximately 75% of the original tensile strength (Break Strength Retention-BSR) in 7 days. It loses its strength in 42 days.

**CONTRAINDICATIONS**

The sutures being absorbable should not be used where extended approximation of tissue is required. The sutures should not be used for cardiovascular and neurosurgical tissue.

**WARNINGS**

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under AECT/ACC section) when selecting a suture for use in patients.

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.

Under some circumstances, notably orthopedic procedures, immobilization of joint by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Do not re-sterilize. Discard open packages and unused sutures. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

**PRECAUTIONS**

In handling FASTSORB or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point.

Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

PGA sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon.

Avoid prolonged exposure to elevated temperatures.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

**ADVERSE REACTIONS**

Adverse effects associated with the use of the device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solution such as urine and bile occurs, and transitory local irritation at the wound site.

**HOW SUPPLIED**

FASTSORB sutures are available as sterile, braided, undyed (beige) strands in sizes 8/0 thru 6 (metric sizes 0.4 – 8.0) in a variety of lengths, with a variety of needles. FASTSORB sutures are available in one, two and three dozen boxes. FASTSORB suture is supplied sterile.

**STORAGE**

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

**SYMBOLS USED ON LABELLING**

	Do not reuse!		Catalogue Number
	Do not re-sterilize		Store below 25°C
	Do not use if pack age is damaged		Keep away from sunlight
	Manufacturer		Protect from humidity
YYYY 	Date of Manufacture, Year		Recyclable pack
YYMM	Expiry Date, Year- Month		Attention, See instruction for use
	Sterile R: Irradiation		2292
	Batch Number		
	Undyed, Absorbable, Braided, Coated		

IFU-FS-rev-06-15-01-2018

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

**BOZ** SAGLIK MALZEME SANAYI VE TICARET A.Ş.  
Saglik Mah. Sıgirik Sokak No.33/5Sihhiye-Çankaya, ANKARA/TURKIYE  
Tel: +90 (312) 254 03 40 (5 hat) Faks:+90 (312) 254 03 50  
web: www.boztibbi.com e-mail: boz@boztibbi.com

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