Prospektüs

Prospectus

GLIKOLAK

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GB

GLIKOLAK**Instruction for Use**

Surgical, Poly(glycolide(90%-co-lactide)(10%))
Absorbable PGLA Suture
Synthetic, Multifilament, Braided, Violet/Undyed

DESCRIPTION

GLIKOLAK suture is synthetic absorbable sterile surgical suture which is flexible strand composed of a copolymer synthesized from 90% glycolide and 10 % L-lactide (derived from glycolic and lactic acids). The empirical formula of the copolymer is $(C_4H_6O_4)_m(C_3H_4O_3)_n$. Multifilament suture is consisting of elementary fibres which are assembled by braiding.

GLIKOLAK suture is prepared by coating suture material with a mixture composed of equal parts calcium stearate ($C_{18}H_{35}O_2$) and 30% glycolide and 70% lactide ($C_3H_4O_3$) where $m+n=3.7$. When GLIKOLAK suture is introduced into a living organism it is absorbed by that organism cause no undue tissue irritation. Copolymer and the coating with calcium stearate have been found to be nonantigenic, nonpyrogenic and elicit only a mild tissue reaction during absorption.

GLIKOLAK suture has been presented in USP 80 to 2 (metric sizes 0.4 and 5) gauge sizes in a variety of lengths, as non-needed or attached to stainless steel needles of varying types and sizes. The suture is colored with violet (D&C Violet No. 2 C.I. # 60725, FDA Part 74.3602 / Solvent Violet 13) to enhance visibility in tissue. The suture is also available in the undyed (natural beige) form.

GLIKOLAK sutures meet United States Pharmacopoeia (U.S.P.) requirements as described in the U.S.P. "Monograph for Absorbable Surgical Sutures" and European Pharmacopoeia (E.P.) requirements as described in the E.P. "Monograph for Sutures, Sterile Synthetic Absorbable Braided; 012008.0667" with the exception of an occasional slight oversize in some gauges.

INDICATIONS

GLIKOLAK suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues.

ACTIONS

GLIKOLAK suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of suture occurs by means of hydrolysis gradually and decrease the strength in the body. After hydrolysis it is excreted from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass.

Animal and in-vitro hydrolysis studies indicate that GLIKOLAK suture retains approximately 75% of the original tensile strength (Break Strength Retention-BSR) after two weeks post implantation. This suture provides tissue supporting during 21-28 days and suture absorption is essentially complete between 55 and 70 days.

Days Implantation	Approximate % Original Strength Remaining (%Break Strength Retention-BSR)
	GLIKOLAK
7 Days	%50
14 Days	%75
21 Days (B-S and larger)	%50
21 Days (I-S and smaller)	%40
28 Days (B-S and larger)	%25

CONTRAINDICATIONS

This suture, being absorbable, should not be used where extended approximation of tissues is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

As with any foreign body, prolonged contact of any suture with soft solutions, such as those found in the urinary or biliary tracts may result in calculus formation.

As GLIKOLAK is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching, or distension, or which may require additional support. Do not re-sterilize. Discard open packages and unused sutures.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed. Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur.

In handling GLIKOLAK or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

To avoid damaging needle points and suture areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the suture end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

GLIKOLAK suture, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon.

Avoid prolonged exposure to elevated temperatures.

Depose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse reactions associates with the use of GLIKOLAK suture include transitory local irritation at the wound site, acute inflammatory tissue reaction, erythema and induration during the absorption process of subcuticular sutures. Like all foreign bodies GLIKOLAK may potentiate an existing infection.

HOW SUPPLIED

GLIKOLAK suture is available sterile, as braided dyed (violet) and undyed (natural) strands in sizes 8/0 to 2 (metric sizes 0.4 - 5) in a variety of lengths, as non-needed or attached to stainless steel needles of varying types and sizes. Suture is available one, two and three dozen in boxes. GLIKOLAK suture is for single use only.

STORAGE

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

Do not reuse!

Do not resterilize

Do not use if package is damaged

Manufacturer

YYYY-M Date of Manufacture, Year

YYYY-MM Expiry Date, Year-Month

STERILE EO Sterile EO: Ethylene oxide

LOT Batch Number

Dyed/Absorbable, Braided, Coated

Undyed/Absorbable, Braided, Coated

REF Catalogue Number

Store below 25°C

Keep away from sunlight

Protect from humidity

Recyclable pack

Attention: See instruction for use

CE 2292

IFU-GL-rev-06-15-01-2018

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multi-lingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Saklık Mah. Saklık Sokak No:33/5Şişliye-Çankaya, ANKARA/TÜRKİYE
Tel: +90 (312) 254 02 40 (5 hat) Faks: +90 (312) 254 03 50
web: www.boztibbi.com e-mail: boz@boztibbi.com

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