Prospektüs

Prospectus

GLIKOSORB

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BOZ

GB

GLIKOSORBSurgical, Polyglycolic Acid
Absorbable PGA Suture
Synthetic, Multifilament, Braided, Undyed/Violet

Instruction for Use

DESCRIPTION

GLIKOSORB Surgical Polyglycolic acid (PGA) suture; synthetic, sterile, absorbable, is homopolymer of glycolic acid. The empirical formula of polyglycolic acid is $(-O-CH_2-CO-O-CH_2-CO-)_n$. GLIKOSORB is coated with the mixture of calcium stearate $(C_{18}H_{35}O_2Ca)$ and polyacrylonitrile $(C_3H_3.5O_n)$. The suture and coating material are inert, noncollagenogenic, nonantigenic, and non-pyrogenic. Coloured violet (D&C Violet No. 2 C.I. # 60725 / Solvent Violet 13) to enhance visibility in tissue.

INDICATIONS

GLIKOSORB suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological tissues. GLIKOSORB suture is for single use only.

ACTIONS

PGA suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGA suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis it's executed from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass. PGA suture retains approximately 75% of the original tensile strength (Break Strength Retention-BSR) at two weeks post implantation. At three weeks, approximately 50% of the original tensile strength (Break Strength Retention-BSR) is retained for sizes 6-0 and larger and approximately 40% of its original tensile strength (Break Strength Retention-BSR) is retained for sizes 7-0 and smaller. All of the original tensile strength (Break Strength Retention-BSR) is lost between four and five weeks post implantation. In 60-90 days it is totally absorbed.

CONTRAINDICATIONS

The sutures being absorbable should not be used where extended approximation of tissue is required. The sutures should not be used for cardiovascular and neurological tissue.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients.

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.

Under some circumstances, notably orthopedic procedures, immobilization of joint by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur.

As with any foreign body, prolonged contact of any suture with soft solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Do not resterilize. Discard open packages and unused sutures.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

PRECAUTIONS

In handling PGA or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

PGA sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon.

Avoid prolonged exposure to elevated temperatures.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of the device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solution such as urine and bile occurs, and transitory local irritation at wound site.

HOW SUPPLIED

GLIKOSORB sutures are available braided, dyed (violet) and undyed strands in sizes 8/0 thru 6 (metric sizes 0.4 – 8) in a variety of lengths, with a variety of needles. GLIKOSORB sutures are available in one, two and three dozen boxes. GLIKOSORB suture is supplied sterile.

STORAGE

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

 Do not heat	 Undyed, Absorbable, Braided, Coated
 Do not resterilize	 Catalogue Number
 Do not use if pack age is damaged	 Store below 25C
 Manufacturer	 Keep away from sunlight
YYYY  Date of Manufacture, Year	 Protect from humidity
 YYYY-MM Epiry Date, Year-Month	 Recyclable pack
 STERILE EO Sterile EO: Ethylene oxide	 Attention, See instruction for use
 LOT Batch Number	 Dyed, Absorbable, Braided, Coated
 Dyed, Absorbable, Braided, Coated	 CE 2292

IFU-GS-rev-06-15-01-2018

Issue date: 11.09.2012

"EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Saglık Mah. Sağlık Sokak No:33/35/İhyye Çarşakı ANKARA/TÜRKİYE
Tel: +90 (312) 254 03 40 (5 hat) Faks:+90 (312) 254 03 30
web: www.boztibbi.com e-mail: boz@boztibbi.com

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