Prospektüs Prospectus

MONAMID

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GB

Instructions for Use
MONAMID® NYLON SUTURE

NONABSORBABLE SURGICAL SUTURE, U.S.P/ Eur.Ph./TJS5504

DESCRIPTION

MONAMID nylon suture is a nonabsorbable, sterile surgical monofilament suture composed of the long-chain aliphatic polymers Nylon 6 or Nylon 6,6. MONAMID sutures are dyed black (C.I. color index no: 77266) or blue (C.I. color index no: 74160) to enhance visibility tissue. The suture is also available undyed (clear).

MONAMID suture meets all requirements established by the United States Pharmacopeia (U.S.P) and European Pharmacopeia for nonabsorbable surgical suture.

INDICATIONS

MONAMID suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. MONAMID suture is for single use only.

ACTIONS

MONAMID suture elicit a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While nylon is not absorbed, progressive hydrolysis of the nylon in vivo may result in gradual loss of tensile strength over time.

CONTRAINDICATIONS

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon suture should not be used where permanent retention of tensile strength is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing MONAMID suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Do not resterilize. Discard open packages and unused sutures. Acceptable surgical practices must be followed for the management of contaminated or infected wounds.

PRECAUTIONS

In handling MONAMID or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

MONAMID sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to

the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, gradual loss of tensile strength over time, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction and transitory local irritation at the wound site.

HOW SUPPLIED

MONAMID sutures are available as sterile, monofilament, black, blue or undyed (clear) strands in U.S.P sizes 10/0 thru 2 (metric sizes 0.2-5) in a variety of lengths with and without permanently attached needles. MONAMID sutures are also available attached to Polyterfluoroethylene (PTFE) pledget.

MONAMID sutures are available in one, two, and three dozen boxes.




STORAGE

Store below 25 °C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse		Undyed, Non-absorbable, Monofilament
	Do not resterilize		Catalogue Number
	Do not use if pack is damaged		Store below 25°C
	Manufacturer		Keep away from sunlight
YYYY 	Date of Manufacture, Year		Protect from humidity
	YYYY-MM Expiry Date, Year-Month		Recyclable pack
	STERILE EO Sterile EO: Ethylene oxide		Attention, See instruction for use
	LOT Batch Number		1014
	Dyed, Non-absorbable, Monofilament		

IFU-MA-rev-03-05-05-2017

Issue date: 11.09.2012

"EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
 Sağlık Sokak No:33/5 SİHİYE/ANKARA/TÜRKİYE
 Tel:+90 (312) 435 50 80-81 Faks:+90 (312) 435 41 26
 web: www.boztibbi.com e-mail: boz@boztibbi.com

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