Prospektüs
Prospectus

MONOFLORID

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DESCRIPTION

MONOFLORID suture is a nonabsorbable sterile surgical suture composed of a polyvinylidene difluoride polymer. The suture is pigmented blue (C.I. # 61568 or C.I. # 74160) or black to enhance visibility.

INDICATIONS

MONOFLORID suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures.

ACTIONS

MONOFLORID suture elicits a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. MONOFLORID suture is not absorbed and it is not degraded or weakened by the action of tissue enzymes. As a monofilament, MONOFLORID suture resists involvement in infection and has been successfully employed in contaminated and infected wounds to eliminate or minimize later sinus formation and suture extrusion. The lack of adherence to tissues has facilitated removing of MONOFLORID suture as a pull-out suture.

CONTRAINDICATIONS

Not known

WARNINGS

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing MONOFLORID suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Do not re-sterilize. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Monofilament sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of MONOFLORID suture include wound dehiscence, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site.

HOW SUPPLIED

MONOFLORID sutures pigmented are available as sterile strands in U.S.P. sizes 10/0 to 2 (metric sizes 0.2-5).

MONOFLORID sutures are available in a variety of lengths, with permanently attached needles or non-needed.

MONOFLORID sutures are available in one, two and three dozen boxes. MONOFLORID suture is for single use only.





STORAGE

Store below 25°C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse!		Catalogue Number
	Do not re-sterilize		Store below 25°C
	Do not use if pack age is damaged		Keep away from sunlight
	Manufacturer		Protect from humidity
YYYY 	Date of Manufacture, Year		Recyclable pack
	YYYY-MM Expiry Date, Year-Month		Attention, See instruction for use
	STERILE EO Sterile EO: Ethylene oxide		CE 1014
	LOT Batch Number		
	Dyed, Non-absorbable, Monofilament		

IFU-MD-rev-03-05-05-2017

Issue date: 11.09.2012

"EASSI"(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
 Sağlık Sokak No:33/5 Şişliye/ANKARA/TÜRKİYE
 Tel:+90 (312) 435 50 80-81 Faks:+90 (312) 435 41 26
 web: www.boztibbi.com e-mail: boz@boztibbi.com

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