Prospektüs

Prospectus

MONOKAPROL PLUS

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Instructions for Use

MONOKAPROL PLUS

Poly(glycolide (75%)-co-caprolactone(25%)) (PGLCL) Surgical Suture
SYNTHETIC, ABSORBABLE, MONOFILAMENT, USP/IEP

DESCRIPTION

MONOKAPROL PLUS is a short-term absorbable synthetic monofilament suture made of monofilament available in the closure of superficial wounds and in violet color for the rest of indications. MONOKAPROL PLUS offers excellent handling properties and high knot security.

MONOKAPROL PLUS suture provides an antibacterial effect by means of a broad spectrum antibacterial agent Chlorhexidine Diacetate (CHA) content no more than 60µg/cm.

INDICATIONS

MONOKAPROL PLUS is recommended for all soft tissue approximations and/or ligations including gastrointestinal surgery, gynaecology/obstetrics, urology, plastic and reconstructive surgery, skin closure (intra, sub, skin). But it shouldn't be used in cardiovascular or neurological tissues; microsurgery or ophthalmic surgery. MONOKAPROL PLUS suture is for single use only.

CONTRAINDICATIONS

This suture, being absorbable, should not be used where extended approximation of tissue under stress is required. MONOKAPROL PLUS suture should not be used in patients with known allergic reactions to Chlorhexidine Diacetate (CHA).

ACTIONS

PGLCL suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGLCL suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis is excluded from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass. PGLCL suture retains 75-70% of the original tensile strength in 7 days after implantation. At two weeks, approximately 40-50% of the original strength is retained. Complete mass is absorbed in 90 to 120 days.

MONOKAPROL PLUS suture has been shown to inhibit microbial adherence and colonization of the suture by E.coli, Klebsiella pneumoniae, Staphylococcus epidermidis, Staphylococcus aureus and its Methicillin resistant strains using in-vitro and in-vivo studies.

WARNINGS

Sutures should be familiar with surgical procedures and techniques involving absorbable sutures before employing MONOKAPROL PLUS suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.

The use of MONOKAPROL PLUS does not substitute normal observation of hygiene and/or otherwise needed antibiotic treatment.

Do not resterilize. Discard opened packages and unused sutures.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, MONOKAPROL PLUS suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching, or contraction, or which may require additional support.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration. Normally associated with absorption. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur.

In handling MONOKAPROL PLUS or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

MONOKAPROL PLUS sutures, which are treated to enhance handling characteristics, require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

Avoid prolonged exposure to elevated temperature.

To avoid damaging needle points and suture areas, grasp the needle in the area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. User should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site.

HOW SUPPLIED

MONOKAPROL PLUS sutures are available in undyed, dyed (violet) strands, in sizes USP 5/0 (1 metric) to USP 2 (5 metric), in a variety of lengths, as non-needed or attached to stainless steel needles, in various suture types and sizes.

MONOKAPROL PLUS sutures are available in one, two, and three dozen boxes.

MONOKAPROL PLUS suture is supplied sterile.

STORAGE

Store below 25°C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse		Undyed Absorbable Monofilament
	Do not resterilize		Catalogue Number
	Do not use if pack age is damaged		Store below 25°C
	Manufacturer		Keep away from sunlight
	Date of Manufacture/Year		Protect from humidity
	Expiry Date, Year-Month		Recyclable pack
	Sterile EO: Ethylene oxide		Attention, See instruction for use
	Batch Number		1014
	Dyed Absorbable Monofilament		

IU-MKP-REV-02-05-05-2017

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Saglık Sokak No:33/5 Sincihüyük/ANKARA/TÜRKİYE
Tel:+90 (312) 635 50 80-81 Faks:+90 (312) 432 41 26
web: www.boztibbi.com e-mail: boz@boztibbi.com