Prospektüs

Prospectus

MONOKAPROL

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GB

MONOKAPROL**Instruction for Use**

Surgical, Poly(glycolide)(75%)-co-caprolactone(25%)]
 Absorbable PGCL Suture
 Synthetic, Monofilament, Undyed/Violet

DESCRIPTION

MONOKAPROL is a short term absorbable synthetic monofilament suture made of poly(glycolide-co-caprolactone) available in undyed color for the closure of superficial wounds and violet color for the rest of indications. MONOKAPROL offers excellent handling properties and high knot security.

INDICATIONS

MONOKAPROL is recommended for all soft tissue approximations and/or ligations including gastrointestinal surgery, gynecology/obstetrics, urology, plastic and reconstructive surgery, skin closure (intra, sub, skin). But it shouldn't be used in cardiovascular or neurological tissues; microsurgery or ophthalmic surgery. MONOKAPROL suture is for single use only.

CONTRAINDICATIONS

This suture, being absorbable, should not be used where extended approximation of tissue under stress is required.

ACTIONS

PGCL suture elicits a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of PGCL suture occurs by means of hydrolysis gradually and decreases the strength in the body. After hydrolysis it is executed from the body as carbon dioxide and water. Absorption begins as a loss of tensile strength followed by a loss of mass. PGCL suture retains 45-70% of the original tensile strength (Break Strength Retention-BSR) in 7 days after implantation. At two weeks, approximately 20% of the original tensile strength (Break Strength Retention-BSR) is retained. Complete mass is absorbed in 90 to 120 days.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture MONOKAPROL suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds. As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distension or which may require additional support. Do not resterize. Discard opened packages and unused sutures.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off removed as indicated. Subcuticular sutures should be placed as deep as possible to minimize the erythema and induration normally associated with absorption.

Under some circumstances notably orthopedic procedures immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable sutures in tissue with poor blood supply as suture extrusion and delayed absorption may occur.

In handling MONOKAPROL or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Avoid prolonged exposure to elevated temperature.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

MONOKAPROL sutures, which are treated to enhance handling characteristics require the accepted surgical technique of flat and square lites with additional throws as warranted by surgical circumstances and experiences of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Dispose of contaminated and unused products as in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solution such as urine and bile occurs, and transitory local irritation at the wound site.

HOW SUPPLIED

MONOKAPROL sutures are available in undyed, dyed (violet) strands, in sizes USP 6/0 (0.7 metric) to USP 2 (5 metric) in a variety of lengths with a variety of needles. MONOKAPROL sutures are available in one, two and three dozen boxes. MONOKAPROL suture is supplied sterile.

STORAGE

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

IFU-MK-rev-05-15-01-2018

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
 Sığılı Mah. Sığılı Sokak No:3/35Sığılı-Garipaşa, ANKARA/TÜRKİYE
 Tel: +90 (312) 254 03 40 (5 hat) Faks:+90 (312) 254 03 50
 web: www.boztibbi.com e-mail: boz@boztibbi.com

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