Prospektüs
Prospectus

MONOPROLEN

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NONABSORBABLE SURGICAL SUTURE, U.S.P./ Eur.Ph./TSEK

DESCRIPTION

MONOPROLEN polypropylene suture (clear or pigmented) is a monofilament, non-absorbable, sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The suture is pigmented blue (C.I. #61568 or C.I. #74160) and black (C.I. #77266) to enhance visibility.

INDICATIONS

MONOPROLEN suture is indicated for use in general soft tissue approximation and /or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures.

ACTIONS

MONOPROLEN suture elicits a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. MONOPROLEN suture is not absorbed and it is not degraded or weakened by the action of tissue enzymes. As a monofilament, MONOPROLEN suture, resists involvement in infection and has been successfully employed in contaminated and infected wounds to eliminate or minimize later sinus formation and suture extrusion. The lack of adherence to tissues has facilitated removing of MONOPROLEN suture.

CONTRAINDICATIONS

Not known.

WARNINGS

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing MONOPROLEN suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Do not resterilize. Discard opened packages and unused sutures. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice must be followed for the management of infected or contaminated wounds.

PRECAUTIONS

In handling this suture material, care should be taken to avoid damage to needle and suture from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties of single suture strands. The use of additional throws is particularly appropriate when knotting polypropylene sutures.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of MONOPROLEN suture include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site.


HOW SUPPLIED

MONOPROLEN sutures pigmented or clear, are available as sterile strands in USP sizes 10/0 to 2 (metric sizes 0.2 to 5). MONOPROLEN sutures are available in a variety of lengths, with permanently attached needles. MONOPROLEN sutures, pigmented and clear, are also available attached to Polytetrafluoroethylene (PTFE) pledgets. MONOPROLEN sutures are available in one, two, and three dozen boxes. MONOPROLEN suture is for single use only.

STORAGE

Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse		Undyed, Non-absorbable, Monofilament
	Do not resterilize		Catalogue Number
	Do not use if pack age is damaged		Store below 25C
	Manufacturer		Keep away from sunlight
YYYY 	Date of Manufacture, Year		Protect from humidity
	YYYY-MM Expiry Date, Year-Month		Recyclable pack
	Sterile EO: Ethylene oxide		Attention, See instruction for use
	Lot Number		CE 1014
	Dyed, Non-absorbable, Monofilament		

IUJ-ML-rev-04-05-05-2017

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBI MALZEME SANAYI VE TICARET A.Ş.
Saglık Bostak No:33/5 Siihkiye/ANKARA/TURKIYE
Tel:+90 (312) 435 50 80-81 Faks:+90 (312) 435 41 26
web: www.boztibbi.com e-mail: boz@boztibbi.com

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