# Prospektüs

Prospectus

# MONOTEF

## Emilemeyen Cerrahi Ameliyat İpliđi

Non-Absorbable Surgical Sutures

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MONOTEF suture is a non-absorbable, monofilament, sterile, surgical suture composed of Polytetrafluoroethylene. It is manufactured from 100% high-density polytetrafluoroethylene (PTFE) polymer, extruded in such a fashion as to produce a structure with a minimal pore size and volume while maintaining integrity and tensile strength. The suture is undyed (white), contains no additives and free of PFOA.

MONOTEF sutures meet all requirements for non-absorbable surgical sutures established by the United States Pharmacopoeia (U.S.P) and European Pharmacopoeia (E.P.).

**INDICATIONS:**

MONOTEF suture is indicated for use in all types of soft tissue approximation and/or ligation, including use in cardiovascular, dental and general surgeries, as well as repair of the dura mater. The device is not indicated for use in ophthalmic surgery and peripheral neural tissue.

**ACTIONS:**

PTFE is one of the most inert materials known and MONOTEF suture elicits minimal tissue reaction. MONOTEF Suture is not absorbed or subject to weakening by the action of tissue enzymes. It does not degrade in the presence of infection.

**CONTRAINDICATIONS:**

Not known.

**WARNINGS:**

The safety and effectiveness of this suture in peripheral neural, and ophthalmic applications has not been established.

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing MONOTEF suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Do not sterilize. Discard opened packages and unused sutures. Do not use if package is opened or damaged.

Tissue invasion of MONOTEF Suture can result in attachment of the suture to the tissue it penetrates in long-term use. Such attachment may make removal of the suture difficult.

Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

**PRECAUTIONS:**

In handling MONOTEF or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle

holders. Sutures should be handled preferably with a needle holder at a third or half of the distance between the connection to the braid and the point, to avoid penetration difficulties, rupture or needle bending. Changing the bending of the needle leads to a torsion resistance loss and can cause the breakage.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. When tying knots with the MONOTEF Suture, tension should be applied by pulling each strand of the suture in opposite directions with equal force.

**ADVERSE REACTIONS:**

Adverse effects associated with the use of this device include wound dehiscence, infection, localized inflammatory reaction, and transitory local irritation at the wound site.

**HOW SUPPLIED:**

MONOTEF sutures are available as sterile, monofilament, undyed strands in U.S.P. sizes 7/0 to 3/0 (metric sizes 0,5- 2) in a variety of lengths, with and without permanently attached needles.

MONOTEF sutures are available in half, one, two, and three dozen boxes.

MONOTEF suture is for single use only.


**STORAGE:**

Store below 25°C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

**SYMBOLS USED ON LABELLING**

	Do not reuse		Catalogue Number
	Do not sterilize		Store below 25°C
	Do not use if pack age is damaged		Keep away from sunlight
	Manufacturer		Protect from humidity
YYYY 	Date of Manufacture, Year		Recyclable pack
	YYYY-MM Expiry Date, Year- Month		Attention, See instruction for use
	Sterile EO: Ethylene oxide		1014
	Batch Number		
	Undyed, Non-absorbable, Monofilament		

IFU-MT-rev-02-05-05-2017

Issue date: 19.01.2015

"EASSI (The European Association of Surgical Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

**BOZ TIBBI MALZEME SANAYI VE TICARET A.Ş.**  
Sığırcı Sokak No:33/5 Şişli/Beşiktaş/İSTANBUL  
Tel:+90 (312) 435 50 80-81 Faks:+90 (312) 435 41 26  
web: www.boztibbi.com e-mail: boz@boztibbi.com