Prospektüs
Prospectus

MONOWIRE

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Instructions for Use

MONOWIRE® SURGICAL STEEL WIRE

NONABSORBABLE, SURGICAL STAINLESS, USP/ Eur.Ph.

DESCRIPTION

Surgical stainless steel suture is a nonabsorbable, sterile surgical suture composed of AISI 316L VM grade stainless steel conforming to ISO 6832 Part 1 "Metallic Materials for surgical implants". Surgical stainless steel suture is available as a monofilament and multifilament suture. Surgical stainless steel suture meets all requirements established by the European Pharmacopoeia (E.P.) and the United States Pharmacopoeia (U.S.P.) for nonabsorbable surgical sutures. Surgical stainless steel suture may also be labelled with the Brown & Sharpe Gauge (B&S) or American Wire Gauge (AWG) classification.

Pharmacopoeia sizes of steel suture corresponding to the B&S / AWG gauges		
USP	EP	B&S / AWG Gauge
2/0	2	3/16
3/0	3	1/8
4/0	4	5/32
5/0	5	3/32
6/0	6	1/16
7/0	7	3/64
8/0	8	1/32
9/0	9	3/64
10/0	10	1/64
11/0	11	3/128
12/0	12	1/64
13/0	13	3/256
14/0	14	1/128
15/0	15	3/256
16/0	16	1/128
17/0	17	3/512
18/0	18	1/256
19/0	19	3/512
20/0	20	1/256
21/0	21	3/1024
22/0	22	1/512
23/0	23	3/1024
24/0	24	1/512
25/0	25	3/2048
26/0	26	1/1024
27/0	27	3/2048
28/0	28	1/1024
29/0	29	3/4096
30/0	30	1/2048
31/0	31	3/8192
32/0	32	1/1024
33/0	33	3/16384
34/0	34	1/512
35/0	35	3/16384
36/0	36	1/512
37/0	37	3/32768
38/0	38	1/256
39/0	39	3/32768
40/0	40	1/256
41/0	41	3/65536
42/0	42	1/128
43/0	43	3/65536
44/0	44	1/128
45/0	45	3/131072
46/0	46	1/64
47/0	47	3/131072
48/0	48	1/64
49/0	49	3/262144
50/0	50	1/32
51/0	51	3/262144
52/0	52	1/32
53/0	53	3/524288
54/0	54	1/16
55/0	55	3/524288
56/0	56	1/16
57/0	57	3/1048576
58/0	58	1/8
59/0	59	3/1048576
60/0	60	1/8
61/0	61	3/2097152
62/0	62	1/4
63/0	63	3/2097152
64/0	64	1/4
65/0	65	3/4194304
66/0	66	1/2
67/0	67	3/4194304
68/0	68	1/2
69/0	69	3/8388608
70/0	70	1/1
71/0	71	3/8388608
72/0	72	1/1
73/0	73	3/16777216
74/0	74	1/1
75/0	75	3/16777216
76/0	76	1/1
77/0	77	3/33554432
78/0	78	1/1
79/0	79	3/33554432
80/0	80	1/1

Polypropylene (PP) button, the accessory, is used for buttressing the suture thread to avoid cutting the tissue that the suture thread intended to be used on tendon repair. The set involving the button is called "Tendon Repair Kit". The curved side is the top of the button. Flat side should be used to spread suture tension by dermal touch. Another accessory is Polyethylene (PE) oval plate, which is used for buttressing/cushioning the suture thread to avoid cutting the tissue and to spread the suture tension of where the suture thread intended to be used on abdominal closure. The set involving the plate is also called "Retention Suture". The plate consists of a hard high density polyethylene (HDPE) layer at the top, a low density polyethylene (LDPE) foam in the bottom layer which will be in touch with the area of application. Both layers are connected by double sided adhesive tape.

INDICATIONS

Surgical stainless steel suture is indicated for use as nonabsorbable suture in tissue approximation such as abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.

ACTIONS

Surgical stainless steel suture elicits a minimal acute inflammatory reaction in tissue and is not absorbed.

CONTRAINDICATIONS

The use of this suture is contraindicated in patients with known sensitivities or allergies to 316L stainless steel and/or its principal metallic components metals such as chromium and nickel.

The presence of steel may interfere with certain radio diagnostics and its use is contraindicated where radio transparency of suture material is required.

WARNINGS

Users should be familiar with surgical procedures and techniques involving non-absorbable, stainless steel sutures before employing for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Stainless steel sutures are not to be used in direct contact with the heart, the central circulatory system or the central nervous system. Acceptable surgical practice must be followed for the management of contaminated or infected wounds.

Accessories are non-absorbable, short term external use, not for inner body or cavity usage, only dermal touch. Do not resterilize. Discard opened, unused sutures.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage to needle and suture, such as kinking or excessive twisting.

Acceptable surgical practice should be followed to place steel sutures sufficiently distal to wound edges to reduce the incidence of the suture cutting through tissue, when the wounds is placed under extreme stress.

Care should be taken to cut the knots/twists to allow for the bending of sharp ends away from surrounding tissue.

Grasping at the attachment area could cause breakage of the needle barrel or the suture in the attachment area. Grasping needle near the point could result in fracture of the needle.

Reshaping needles may result in decreased resistance to bending and breaking. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point.

Avoid inadvertent needle stick injury, which may result in transmission of bio- borne pathogens. Suitable containers must be used for disposal of needles. Dispose of contaminated suture and packaging materials utilizing standard hospital procedures and universal precautions for biohazardous waste.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include dehiscence, allergic response in patients with known sensitivities to 316L stainless steel, or constituent metals such as chromium and nickel infection, minimal acute inflammatory tissue reaction, pain, edema, and local irritation at the wound site. Needle breakage during surgery may result in extended anesthesia and/or surgery time, retention of foreign body in the patient, and/or additional surgery.

HOW SUPPLIED

Surgical stainless steel sutures are available in sizes 6/0 to 8 (metric sizes 0.7 – 10.0) in a variety of pre-cut lengths with and without permanently attached needles by permanent needle attachment or rotating techniques, in one, two and three dozen boxes. MONOWIRE suture is supplied as sterile. MONOWIRE® is for single use only.

STORAGE

Keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

Do not reuse

Catalog Number

Do not resterilize

Store below 25°C

Do not use if pack age is damaged

Keep away from sunlight

Manufacturer

YYYY MM Date of Manufacture, Year

Protect from humidity

YYYY-MM Expiry Date, Year- Month

Recyclable pack

STERILE EO Sterile EO: Ethylene oxide

Attention, See instruction for use

LOT Lotch Number

CE 1014

Unsterilized, Non-absorbable, Monofilament

IU-MW-rev-04-05-05-2017

Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Sığlık Sokak No:33/5 Sırtlıbaşı/ANKARA/TÜRKİYE
Tel: +90 (312) 435 50 80-81 Faks: +90 (312) 435 41 26
web: www.boztibbi.com e-mail: boz@boztibbi.com

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