Prospektüs

Prospectus

PLAIN CATGUT

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Instructions for Use

CATGUT

ABSORBABLE, SURGICAL SUTURES, USP/EP/TS3497

DESCRIPTION

Surgical gut suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from submucosal fibrous layer of sheep (ovine) intestines. Surgical gut sutures are available in plain or chromic. Chromic gut is chemically treated with chromic salt solutions to provide greater resistance to absorption. Plain catgut is naturally cream coloured. Chromic catgut is coloured brown by pyrogallol. Chromic gut suture is also available in blue (Indigo) dyed. Surgical catgut is packaged in dry or in conditioning fluid. The fluid is consist of;

- Isopropyl alcohol, % 89.50
- Diethyl ethanolamine, %1.00
- Sodium benzoate, %0.25
- Water, %9.25

Surgical gut suture meets all requirements established by the United States Pharmacopoeia (U.S.P.) monographs.

INDICATIONS

Surgical gut suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.

ACTIONS

When surgical gut suture is placed in tissue, a moderate tissue inflammation response which is characteristic response of body to a foreign substance. This is followed by loss of tensile strength and a loss of suture mass, as the proteolytic enzymatic digestive process dissolves the surgical gut. This process continues until the suture is completely absorbed. Many variable factors may affect the rate of absorption. Some of the major factors which can affect tensile strength loss and absorption rates are:

- Type of suture - plain gut generally absorbs more rapidly than chromic gut.
- infection - surgical gut is absorbed more rapidly in infected tissue than in non-infected tissue.
- Tissue sites - surgical gut will absorb more rapidly in tissue where increased levels of proteolytic enzymes are present as in the secretions exhibited in the stomach, cervix, and vagina.

CONTRAINDICATIONS

This suture, being absorbable, should not be used where extended approximation of tissue is required.

The use of this suture is contraindicated in patients with known sensitivities or allergies to collagen or chromium, as gut is a collagen based material, and chromic gut is treated with chromic salt solutions.

WARNINGS

Users should be familiar with surgical procedures and techniques involving gut suture before using surgical gut suture for wound closure, as the risk of wound dehiscence may result with the site of application and the suture material used. User should consider the in vivo performance when selecting a suture. The use of this suture may be inappropriate in elderly malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.

Do not resterilize. Discard open packages and unused sutures.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion,

stretching, or distension, or which may require additional support.

PRECAUTIONS

In handling surgical CATGUT suture or any other suture material, care should be taken to avoid to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needleholders.

Under some circumstances, notably orthopaedic procedures immobilization of joints by external support may be employed at the discretion of the surgeon.

The surgeon should avoid unnecessary tension when knotting, to reduce the occurrence of surface fraying and weakening of the strand.

Surgical gut sutures require the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Avoid prolonged exposure to elevated temperatures.

ADVERSE ACTIONS

Adverse effects associated with the use of this device include wound dehiscence, variable rates of absorption over time (depending on such factors as the type of suture used, the presence of infection, and the tissue site), failure to provide adequate wound support in closure of sites where expansion, stretching, or distension occurs, etc., unless additional support is supplied through the use of nonabsorbable suture material failure to provide adequate wound support in elderly malnourished, or debilitated patients, or in patients suffering from cancer, anaemia, obesity, diabetes, infection, or other conditions which may delay wound healing allergic response in patients with known sensitivities to collagen or chromium which may result in an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound supuration and bleeding, as well as sinus formation, infection, moderate tissue inflammatory response characteristic of foreign body response, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site.

HOW SUPPLIED

Surgical gut sutures are available in U.S.P. sizes 6/0 thru 4 (metric sizes 1 - 8) in a variety of lengths with and without permanently attached needles. The suture is supplied sterile in one, two, and three dozen boxes. Surgical gut suture is supplied sterile.

STORAGE

Store below 25°C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse		Lot Batch Number
	Do not resterilize		Catalogue Number
	Do not use if package is damaged		Store below 25°C
	Manufacturer		Keep away from sunlight
	YYYY Date of Manufacture/Year		Protect from humidity
	YYYY-MM Expiry Date, Year-Month		Recyclable pack
	Sterile EO: Ethylene oxide		Attention, See instruction for use
	Sterile R: Irradiation		

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