Prospektüs

Prospectus

POLIDIOX

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Boz

GB**POLIDIÖX**

Surgical, Polydioxanone
Absorbable PDO Suture
Synthetic, Monofilament, Violet

Instruction for Use**DESCRIPTION**

POLIDIÖX (polydioxanone) monofilament synthetic absorbable suture is prepared from the polyester, poly(p-dioxanone). The empirical molecular formula of the polymer is $(C_8H_{12}O_2)_x$. POLIDIÖX sutures are dyed violet (D&C Violet No. 2 C.I. # 60725/ Solvent Violet 13) to enhance visibility in tissue. Polydioxanone polymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

INDICATIONS

POLIDIÖX monofilament synthetic absorbable sutures are indicated for use in all types of soft tissue approximation, including use in pediatric cardiovascular tissue where growth is expected to occur and ophthalmic surgery. POLIDIÖX suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue. These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to the six week) is desirable.

ACTIONS

POLIDIÖX synthetic absorbable suture has been formulated to provide wound support through an extended healing period. The results of implantation studies of POLIDIÖX monofilament suture in animals indicate that approximately 70% of its original strength remains (Break Strength Retention-BSR) two weeks after implantation. POLIDIÖX suture provides tissue supporting during 60 days post implantation. Absorption is essentially complete within six months (180-240 days).

CONTRAINDICATIONS

These sutures, being absorbable, are not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and are not to be used in conjunction with prosthetic devices, i.e., heart valves or synthetic grafts.

WARNINGS

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients.

The safety and effectiveness of POLIDIÖX (polydioxanone) sutures have not been established in neural tissue, adult cardiovascular tissue, or for use in microsurgery. Under certain circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon. Do not resterilize.

PRECAUTIONS

In handling POLIDIÖX or any other suture materials, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

PDO sutures, which are treated to enhance handling characteristics requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and experiences of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

Conjunctival and vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation. Subcuticular sutures should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

Dispose of contaminated and unused products are in accordance with local and facility requirements. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Due to prolonged suture absorption, some local inflammatory irritation, acute inflammatory tissue reaction, wound/surgical site infection, foreign body reaction, wound dehiscence, and bleeding have been observed.

HOW SUPPLIED

POLIDIÖX sutures are available as sterile, monofilament, dyed (violet) strands in sizes 7/0 to 2 (metric sizes 0.5-5.0), and sterile, in a variety of lengths with permanently attached needles.

STORAGE

Store below 25°C and keep away from sunlight.

Protect from humidity.

Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not reuse		Catalogue Number
	Do not resterilize		Store below 25°C
	Do not use if pack age is damaged		Keep away from sunlight
	Manufacturer		Protect from humidity
	YYYY-MM Date of Manufacture, Year		Recyclable pack
	YYYY-MM Expiry Date, Year- Month		Attention, See instruction for use
	STERILE EO Sterile EO; Ethylene oxide		2292
	LOT Batch Number		
	Dyed, Absorbable, Monofilament		

IFU-PX-ref-07-15-01-2018

Issue date: 11.09.2012

"EASS(IE The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of this system is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Saklık Mah. Sağlık Sokak No:33/Sisliye Çarşısı, ANKARA/TÜRKİYE
Tel: +90 (312) 254 03 40 (5 hat) Faks: +90 (312) 254 03 50
web: www.boztibbi.com e-mail: boz@boztibbi.com

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