# Prospektüs

Prospectus

# POLYTEF

## Emilemeyen Cerrahi Ameliyat İpliđi

Non-Absorbable Surgical Sutures

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NONABSORBABLE, SURGICAL SUTURE, U.S.P./Eur. Ph./ TS-5503

## DESCRIPTION:

POLYTEF polyester suture is a non-absorbable, braided, sterile, surgical suture composed of Poly (ethylene terephthalate). It is prepared from fibers of high molecular weight, long-chain, linear polyester. POLYTEF braids are coated with a polymer of tetrafluoroethylene (PTFE). Coating compound which is nonabsorbable, acts as a lubricant and mechanically improves the physical properties of the braided uncoated suture by improving easy passage through tissues and overall handling qualities. Heavy PTFE coating and tight weave provide minimum friction and eliminates dead space. POLYTEF suture is braided for optimal handling properties and dyed green to increase visibility. The dye is D&C #6. POLYTEF sutures meets all requirement for non-absorbable surgical sutures established by the United States Pharmacopoeia (U.S.P.) and European Pharmacopoeia (E.P.).

## INDICATIONS:

POLYTEF suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, vascular, ophthalmic, and neurosurgical procedures.

## ACTIONS:

POLYTEF suture elicits a minimal acute inflammatory reaction in tissue, followed by a gradual encapsulation of the suture by fibrous connective tissue. POLYTEF sutures are not absorbed and there is no degradation or loss of strength.

## CONTRAINDICATIONS:

Not known.

## WARNINGS:

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing POLYTEF suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Do not resterilize. Discard opened packages and unused sutures. Do not use if package is opened or damaged.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Acceptable surgical practice should be followed for the management of infected or contaminated wounds.

## PRECAUTIONS:

In handling POLYTEF or any other suture material, care should be taken to avoid damage to needle and suture. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Sutures should be handled preferably with a needle holder at a third or half of the distance between the connection to the braid and the point, to avoid penetration difficulties, rupture or needle bending. Changing the bending of the needle leads to a torsion resistance loss and can cause the breakage.

As with any suture material, adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. Beware when laser or electro-cautery is used as the sutures can be thermo fused.

## ADVERSE REACTIONS:

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site.

## HOW SUPPLIED:

POLYTEF sutures are available as sterile, braided, green strands in U.S.P. sizes 7/0 to 2/0 (metric sizes 0.5-3) in a variety of lengths, with and without permanently attached needles.

POLYTEF sutures, green, are also available attached to Polytetrafluoroethylene (PTFE) pledget.




POLYTEF sutures are available in half, one, two, and three dozen boxes.

POLYTEF suture is for single use only.

## STORAGE:

Store below 25°C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

## SYMBOLS USED ON LABELLING

	Do not reuse		Catalogue Number
	Do not resterilize		Store below 25°C
	Do not use if pack age is damaged		Keep away from sunlight
	Manufacturer		Protect from humidity
	YYYY MM Date of Manufacture, Year		Recyclable pack
	YYYY-MM Expiry Date, Year- Month		Attention, See instruction for use
	STERILE EO Sterile EO: Ethylene oxide		1014
	LOT Batch Number		
	Dyed, Non-absorbable, Braided, Coated		

RU-PT-rev-02-05-05-2017

Issue date: 19.01.2015

\*EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations.\*

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