Prospektüs Prospectus

TEMPOWIRE

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NONABSORBABLE, USP/ Eur.PH.

DESCRIPTION:

TEMPOWIRE Temporary Pacing Wire is twisted or braided 316L stainless steel wire, with PE or FEP coated shielded, with double retractor.

TEMPOWIRE meets requirements established by European Pharmacopoeia (EP) and United States Pharmacopoeia (USP) for nonabsorbable surgical suture, except insulation thickness. TEMPOWIRE may also be labeled with the CE mark classification.

INDICATIONS:

TEMPOWIRE is indicated for use as temporary cardiac pacing wire to transfer electrical energy along its length between the external pacemaker and the patient's heart for temporary cardiac pacing and intracardiac ECG monitoring and after-heart-block therapy.

TEMPOWIRE has one end stripped of insulation for some length and conductively coupled to a curved lead point suture needle for facilitating insertion in the myocardium. The other end of the insulated lead is conductively attached to a curved cutting edge lead with needles for easy insertion from the anterior to the posterior.

The Keith needle may be equipped with a special break-off groove. The cutting part of the needle is broken off from the groove, and the main ends an electrode plug which is connected to conventional pacing and monitoring equipment.

ACTIONS:

TEMPOWIRE elicits a minimal acute inflammatory reaction and is not absorbed.

CONTRAINDICATIONS:

Use of this device is contraindicated in patients with known sensitivities or allergies to 316L stainless steel, or to constituent metals such as chromium and nickel. TEMPOWIRE is contraindicated if a permanent pacemaker has already been implanted. Use of the temporary cardiac pacing wire is contraindicated when permanent cardiac pacing or monitoring is required. Because it can be difficult to distinguish pace wires in the atrial myocardium from those in the ventricle after the chest is closed, TEMPOWIRE should not be used for defibrillation. The presence of steel may interfere with certain radio diagnostics and its use is contraindicated where radio transparency of tissue materials is required.

WARNINGS:

Users should be familiar with surgical procedures and techniques involving temporary pacing and sensing before insertion of the wire.

The presence of steel may cause electrolytic reactions in the presence of the other metals or alloys.

This device is single use. Known hazards with re-use or reutilization of opened / original packages of the device are cross-contamination, infection and trauma.

Avoid disconnection of the wire from heart during in use. This device is only a conductor, does not involve any power source. Pacing of skeletal muscles or nerves can arise from pacing impulses which are set too high. Use this device with a CE-Approved Temporary Pacing Wire Maker device and follow the manufacturer's instructions for use for temporary pacing and sensing.

It is important that the wires should be placed in myocardial tissue on the heart. Faty tissue cannot be stimulated electrically. Tissue that is damaged because of fatty cardiac infarction is also unsuitable for the placement of pacing leads. Occasionally there is a risk that cannot be eliminated which is referred to as "atrial cavitation". Mechanical irritation of the lead may cause further injury to the myocardium and lead to an increase in the pacing threshold. This is so great that they can cause an exit block. Exit block is the condition in cardiac pacing where the heart stimulation cannot be triggered, even with higher energy output from the pacemaker. Therefore, it is important to make sure that the wire sits unobstructed, in the chest and is fixed at the site where it is inserted into the thorax. This way it can easily follow the heart's movement without extrusion from external mechanical factors.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

Strictly follow the manufacturer's instructions for use to avoid infection.

Do not use at room temperature and avoid prolonged exposure to elevated temperatures.

Discard opened packages and unused devices. Do not use opened or damaged products.

PRECAUTIONS:

In handling this device care should be taken to avoid damage from handling such as kinking or excessive heating. This is useful for insulating, do not damage the wires through the use of electrolysis.

It is recommended that the time of pacing with temporary pace wires should be as short as possible and continuous ECG monitoring of the patient's sinus rate is essential in the pacing threshold, understanding and exit block are possible.

When wires are kept a direct low current flow to the heart. As a result, small amounts of red or current are sufficient to cause the heart to fibrillate. Therefore, users should not touch the bare conductor to come in contact with surfaces which are conductive or wet.

Users should avoid insecure connection and usage of inadequate connectors and adapters. Direct connection (i.e. without an adaptor) of a temporary pace wire to external pacemaker or generator cables is recommended. This will reduce the potential for disconnection and will increase the compatibility of the electrical connection. This will reduce the need for an external pacemaker, ensure that it is compatible with the types of temporary cardiac pacing leads that are used in your department.

All static electricity must be kept far away from the pacing system. Before handling the external pulse generator, the patient cable or insulating leads that are attached to the wires through the electrocardiogram (ECG) leads and the patient, for example by touching the patient at a site remote from the pacing lead. Measurements may be carried out with CE-Approved CE-marked equipment. This must include the use of a CE-Approved CE-marked equipment. It should also be noted that heat or your neck might be transmitted to a focal area of cardiac muscle. Use of generator heat that causes the lead or an arrhythmia. Non-optimal pacing parameters may cause Arrhythmia. The danger of fibrillation resulting from an alternative current leakage is greatly increased when SUPPLY MAINS connected pacemaker is connected to the lead system. Extreme caution has to be taken to have proper grounding of SUPPLY MAINS oriented pacemaker used in the vicinity of the patient. The patient should be connected to a ground conductor and a sufficient available in case of emergency. Under no circumstances may any part of the wire be touched during defibrillation of the patient. Cardiac tissue cells can be destroyed and/or damaged when the wires are fixed to the heart muscle (insertion of the wire), i.e. the heart polepieces are not to be destroyed.

contributes to a temporary increase in the pacing threshold. After fixation of the leads, 2-3 minutes should elapse before pacing threshold measurements so that the injury potential can dissipate and the acute pacing threshold can be accurately measured. Temporary pacemaker wire should be placed in such way without interfering since node Atrial and ventricular pacing wires could cause phrenic nerve irritation and diaphragmatic stimulation and careful choice of site of lead placement should be made to avoid these problems.

Fixation of the temporary pace wire in the heart should be sufficient to guarantee reliable heart pacing during the entire pacing period. Users should remember the method of fixation before removing the wire. Removal of the wire must take place with utmost care and be undertaken by an experienced physician. Disconnect the external pacemaker first. Wires are removed with gentle transcutaneous reaction. The patients' vital signs should be monitored following wire removal to allow early identification of the rare potential complications, watch for signs of tamponade in haemodynamically unstable, preunconscious. Arrhythmias patients are at higher risk of haemorrhagic complications. Ultimately, the fixation that is not in the operating room depend on the experience and preference of the physician.

Patients with temporary myocardial leads should not be exposed to any large electrical or magnetic fields. Therefore, machines that work based on the principle of magnetic resonance imaging (MRI) may not be used for the purpose of diagnosis or treatment of these patients.

Temporary pace wires may stay in the body for a maximum of 30 days. However, it is recommended that wires be removed after 10 days at the most in order to prevent infections and guarantee safe patients. If they are not yet stabilized should subsequently be implanted with intracardiac temporary leads and shields, if necessary, receive a permanent pacemaker.

Hold the straight (theta) needle distal from the breakaway marking (if breakaway) and run it through the thoracic cage superior to the anterior. Strip the needle distal from the breakaway marking. Disconnect the tip of the thoracic wire using the theta as an electrode plug. Connect pin of breakaway needle or needle (if not breakaway) to the patient cable (extension cable) or external pacemaker with following polarity. Verify cable to external pulse generator compatibility with use. Patient cable shall be connected to the non-implantable pulse generator before the pacing wires are connected to the patient cable. Connect the patient cable to the appropriate terminal on the pulse generator. Follow pulse generator manufacturer's instructions and polarity and chamber designation. Cut off the excess wire including the heart needle.

Do not pull down the needle points and swing arms, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the needle points of the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

Dispose of contaminated and unused products as are in accordance with local and facility requirements. Discard used devices in accordance with local regulations.

ADVERSE REACTIONS:

Adverse effects associated with the use of this device include, allergic response in patients with known sensitivities to 316L stainless steel, or constituent metals such as chromium and nickel. Local infection, minimal acute inflammatory tissue reaction, myocardial fibrillation, pacing wire displacement, and arrhythmia (myocardial) infection.


NOW SUPPLIED:

TEMPOWIREs available in sizes USP 30 to USP 6 (metric sizes 2.5-3.5) in a variety of lengths with permanently attached needles, with straight or breakaway needle in one, two, three, and three dozen boxes. TEMPOWIREs are supplied sterile.

STORAGE:

Store below 25°C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

	Do not touch!		Not Used Non-absorbable , Biodegradable, Coated
	Do not sterilize		Catalogue Number
	Do not use if package is damaged		Store below 25°C
	Manufacturer		Keep away from sunlight
YYYY 	Date of Manufacture, Year		Protect from humidity
YYYY-MM 	Expiry Date, Year-Month		Recyclable pack
	Sterile EO, Ethylene oxide		Attention, See instruction for use
	Batch Number		Not Used Non-absorbable , Twisted, Coated

RU-TW-rev-04-05-2017

Issue date: 11.09.2012

"EASSI (The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBI MALZEME SANAYE VE TICARET A.Ş.
 Bağlık Sokak No:33/5 SİHAYI/ANKARA/TURKYE
 Tel: +90 (312) 436 50 80-81 Faks: +90 (312) 433 41 26
 web: www.boztibbi.com e-mail: boztibbi@boztibbi.com

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