REOXCEL®

Absorbable Haemostat, Oxidized Regenerated Cellulose (ORC)

REOXCEL®-Sigma Knit

Absorbable Haemostat, Oxidized Regenerated Cellulose (ORC)

EN Instructions for Use

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TR Kullanma Talimatı

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Absorbable Haemostat, Oxidized Regenerated Cellulose (ORC)

EN
Instructions for Use**DESCRIPTION**

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is a sterile haemostatic preparation available in the flat form with textile structure and made of Oxidized Regenerated Cellulose (ORC) (polyanhydro glucuronic acid). REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Regenerated Cellulose". REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat does not contain any animal or collagen additives. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat can be sutured and cut without fraying. Its structure is stable and should be stored at controlled room temperature. The fabric is pale yellow and has a faint, caramel-like aroma. Although a slight discoloration may occur with time, this does not affect its performance.

ACTIONS

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat can be laid on, held against, closed exactly or sutured to a bleeding surface. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is applied as dry. Action mechanism of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is independent from blood coagulation mechanism of the body. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes when it contacts with blood. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is absorbed from the sites of implantation without tissue reaction when it is used properly in minimal amounts. Absorption depends on tissue type, haemostat amount used and saturation degree of blood. It is completely absorbed in 7-14 days. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr(-) (gram negative) microorganisms including those of:

Staphylococcus aureus	Escherichia coli	Klebsiella pneumoniae
Meticillin Resistant Staphylococcus epidermidis (MRSE)	Salmonella typhimurium	Klebsiella aerogenes
Meticillin Resistant Staphylococcus aureus (MRSA)	Proteus mirabilis	Bacillus subtilis
Staphylococcus epidermidis	Proteus vulgaris	Branhamella catarrhalis
Micrococcus luteus	Shigella boydii	Serratia Marcescens
Vancomycin Resistant Enterococcus sp (VRE)	Shigella flexnerii	Clostridium perfringens
Penicillin Resistant Streptococcus pneumoniae (PRSP)	Shigella dysenteriae	Clostridium tetani
Enterococcus faecalis	Salmonella enteritidis	Listeria monocytogenes
Enterococcus faecium	Bacteroides fragilis	Bordetella pertussis
Enterococcus sp.	Lactobacillus casei	MOT (Mycobacterium Other than Tuberculosis)
Streptococcus pyogenes Group A	Lactobacillus sp.	Mycobacterium phlei
Streptococcus pyogenes Group B	Pseudomonas aeruginosa	Streptococcus faecalis
Streptococcus viridans	Pseudomonas stutzeri	
Corynebacterium xerosis	Enterobacter cloacae	

In addition, the bactericidal effect of REOXCEL - REOXCEL SIGMA KNIT Absorbable Hemostat has been shown in vivo studies with MRSA (Methicillin-resistant Staphylococcus aureus) and E. coli (Escherichia coli).

INDICATIONS

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is designed to arrest capillary bleeding and bleeding from parenchymatous organs and resection areas at surgical interventions. It shall be used to stop bleeding when other haemostatic procedures are not applicable. It is suitable for use in general surgery and digestive surgery, neurosurgery (especially cerebral operations), plastic surgery, orthopedy, gynecology, urology, stomatology, operations to the throat or nose, traumatology and all other branches of surgery such as cardiovascular surgery, implantation of vascular prostheses, surgery to the face and jaw, lung operations, haemorrhoidectomy, biopsies, liver and gall bladder operations, gastric resection, thoracic and abdominal sympathectomies. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat can be applied into cavities (after extirpation of tumors) as well as endoscopic interventions or dental praxis.

CONTRAINDICATIONS

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should not be used for implantation in bone defects, such as fractures since it may interfere with callus formation and may cause cyst formation. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should not be used for packing or wrapping. It should be removed after haemostasis is achieved if you have to use in this manner. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should not be used to control haemorrhage from large arteries. Since body fluids, except whole blood, such as serum, don't react with REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat it should not be used to produce haemostatic effect on non-haemorrhagic serous oozing surface. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should not be used as an adhesion prevention product.

WARNINGS

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is suitable for gamma irradiation and supplied as sterile. Should not be re-sterilized. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is not compatible with autoclaving or ethylene oxide sterilization. Closing REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat in a contaminated wound without drainage may lead to complications

and should be avoided this situation. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should not be moistened with water or saline, because the haemostatic effect of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is greater when it is applied dry. Anti-infective agents, buffering materials or haemostatic substances should not be impregnated to REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat. Since thrombin activity is destroyed due to the low pH of the product, haemostatic effect of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat isn't enhanced by addition of thrombin. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat may be left in tissue when necessary, but it is advisable to remove it once bleeding has been controlled. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine the spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure, because REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat may exert pressure resulting in paralysis and/or nerve damage by swelling. Special care must be taken by physicians, regardless of the type of surgical procedure, to consider the advisability of removing REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat after haemostasis is achieved. Although REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is bactericidal against a wide range of pathogenic microorganisms, it should not be used as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections. Particular attention should be paid when haemostat stay in the body after surgical operations, in case of severely bleedings, although it is recommended to remove it after haemostasis.

PRECAUTIONS

Use only as much REOXCEL - REOXCEL SIGMA KNIT Absorbable Hemostat as is necessary for haemostasis, holding it firmly in place until bleeding stops. In order to facilitate absorption and minimize the possibility of foreign body reaction, remove any excess before closure. In urological procedures, minimal amounts of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.

Use of REOXCEL - REOXCEL SIGMA KNIT Absorbable Hemostat should not be preceded by application of silver nitrate or any other escharotic chemicals since its absorption could be prevented in chemically cauterized areas. If REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.

Precautions should be taken in otolaryngologic surgery to assure that none of the material is aspirated by the patient (examples: controlling hemorrhage after tonsillectomy and controlling epistaxis). When REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is used as a wrap during vascular surgery care should be taken not to apply it too tightly.

ADVERSE REACTIONS

"Encapsulation" of fluid, foreign body reactions, hemorrhage, leucocytosis and haematoma have been reported in some cases when the residue of hemostat stay in the body after hemostasis achieved. Abscess of organs and wound dehiscence have been observed in some cases, even the cause of these cases are not related with Haemostat. Paralysis and nerve damage have been reported when Oxidized Regenerated Cellulose (ORC) Haemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy has been reported. Occasional reports of "burning" and "stinging" sensation and sneezing when Oxidized Regenerated Cellulose (ORC) Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product. Headache, burning, stinging, and sneezing in epistaxis control and other rhinological procedures, and stinging when Haemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) also might be possible, rarely.

POSAGE AND ADMINISTRATION

Sterile technique should be observed in removing REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat from its sterile container. Minimal amount of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat in appropriate size is laid on the bleeding site or held firmly against the tissues until haemostasis is obtained. Required Haemostat amount depends on the nature and intensity of the haemorrhage to be stopped. The haemostatic effect of REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is particularly pronounced when used dry. Moistening the material with water or physiological saline solution is not recommended.

CAUTION

Since REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat cannot be re-sterilized, opened and unused REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should be discarded. Dispose of contaminated or unused products in accordance with local and facility requirements.

NOTES AND INFORMATION ON SHELF-LIFE

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat should be stored dry under controlled conditions (max. 25°C) and protected from direct sunlight, in the original packaging. Shelf life of the product is 3 years. The expiry date REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is printed on the pack. Do not use REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat after this date.

PRESENTATION AND CONTENTS

REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat is available sterile, as an original pack in different sizes. REOXCEL - REOXCEL SIGMA KNIT Absorbable Haemostat strips are sterile-packed individually.

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