EN : Do not reuse ES : No reutilizar PT : Não reutilize IT : Non riutilizzare DE : Nur zum Einmalgebrauch NL : Bestemd voor eenmalig gebruik TR : Tek Kullanımlık GR : Διευχρησιστος PL : Nie stosować ponownie DK : Må ikke genbruges FR : Pour utilisation unique AR : لا تستعمل مرة و واحد FA : تک و جا مورد استفاده قرار نگیرد	EN : Do not resterilize ES : No reesterilice PT : FAVOR não esterilizar novamente IT : Non riesterilizzare DE : Nicht nochmal sterilisieren. NL : Niet opnieuw steriliseren TR : Tekrar steril etmeyiniz GR : Μην ξανά σποστεριζώετε PL : Nie reesterilizować DK : Må ikke resteriliseres FR : Ne pas stériliser à nouveau AR : لا تمجد مرة أخرى FA : مجدداً استفاده قرار نگیرد	EN : Batch Number ES : Código de lote PT : Número de Lote IT : Numero di Lotto DE : Chargennummer NL : Lotnummer TR : Seri No GR : Αριθμός Παρτίδας PL : Numer serii DK : Batchnummer FR : No de série AR : رقم التفتة FA : شماره میل	EN : Manufacturer ES : Productor PT : Fabricante IT : Fabricante DE : Hersteller NL : Producent / fabrikant TR : Üretici GR : Παραγωγός PL : Producent DK : Producent FR : Fabricant AR : المصنع FA : تولید کننده	EN : Catalogue number ES : Numero de Catalogo PT : Referência IT : Numero di catalogo DE : Katalognummer NL : Catalogusnummer TR : Katalog numarası GR : Αριθμός καταλόγου PL : Numer katalogowy DK : Katalognummer FR : Numéro de catalogue AR : رقم القهراس / المرجع FA : شماره کاتالوگ

EN : Date of Manufacture, Year ES : Fecha de fabricación, Año PT : Data de Fabrico, Ano IT : Data di produzione, anno DE : Datum der Produktion Jahr NL : Productiedatum / jaar TR : Üretim Tarihi, Yıl GR : Ημερομηνία κατασκευής, έτος PL : Data produkcji, Rok DK : Dato / År af produktionen FR : Date de production, année AR : تاريخ الإنتاج / السنة FA : تاریخ و سال تولید	EN : Expiry Date, Year-Month ES : Fecha de Vencimiento, Año-Mes PT : Data de validade, Ano-Mês IT : Data di scadenza, anno - mese DE : Datum des Verfalls, Jahr - Monat NL : Uiterste houdbaarheidsdatum / jaar - maand TR : Son Kullanma Tarihi, Yıl-Ay GR : Ημερομηνία λήξης, έτος και μήνας PL : Data ważności, Rok - Miesiąc DK : Dato / År - Måned for Udløb FR : Date d'expiration, année - mois AR : تاريخ انتهاء الصلاحية / السنة / الشهر FA : آخرین تاریخ مصرف، سال-ماه	EN : Keep away from sunlight ES : Mantener alejado de la luz del sol PT : Mantenha longe da luz solar IT : Tenere lontano dalla luce diretta DE : Vor Sonneneinstrahlung schützen NL : Beschermen tegen directe zonnestraling TR : Güneğten uzak tutunuz GR : Νη από ηλιακή ακτινοβολία PL : Chronić przed światłem słonecznym DK : Må ikke udsættes for direkte sollys FR : Protéger du soleil AR : تحفظ بعيد عن أشعة الشمس FA : دور محظی به دور از تابش مستقیم خورشید نگهداری شود	EN : Protect from humidity ES : Proteger de la humedad PT : Proteja contra a Umidade IT : Proteggere dall'umidità DE : Vor Feuchtigkeit schützen NL : Beschermen tegen vocht TR : Kuru yerde muhafaza ediniz GR : Ξηροχώροστο εφόσον υπάρχει υγρασία PL : Chronić przed wilgocią DK : Beskyttes mod fugt FR : Conserver dans un lieu sec AR : تحفظ مكان جاف عن الرطوبة FA : در محل خشک نگهداری شود

EN : Consult instructions for use ES : Vea el manual del usuario PT : Verifique as instruções de uso IT : Vedi il Manuale d'Uso DE : Siehe Gebrauchsanweisung NL : Zie gebruiksaanwijzing TR : Kullanma kılavuzuna bakınız GR : Ανατρέξτε στο εγχειρίδιο χρήσης PL : Patrz: Instrukcja Obsługi. DK : Se brugsanvisningen FR : Se confier au Manuel d'Utilisation AR : اقرأ دليل المستخدم FA : به دستور العمل استفاده نگاه کنید	EN : Attention, See instruction for use ES : Consulte la guía de usuario PT : Consulte o manual de utilização IT : Attenzione! vedere manuale d'uso DE : Achtung! finden Sie in der Bedienungsanleitung NL : Attention! raadpleeg de gebruiksaanwijzing TR : Dikkat! Kullanma kılavuzuna bakınız GR : Προσοχή, ανατρέξτε στο εγχειρίδιο οδηγιών PL : Uwaga, zaprzyj do instrukcji używania DK : Fare! Der henvises til brugsanvisningen FR : Attention! voir le manuel d'utilisation AR : انتبه! انظر الى دليل المستخدم FA : دقت، به دستور العمل استفاده نگاه کنید	EN : Recyclable pack ES : Paquete reciclable PT : Embalagem reciclável IT : Pacco da riciclare DE : Recycelbares paket NL : Recyclebaar pakket TR : Geri dönüştürülebilir paket GR : Ανακυκλώσιμο πακέτο PL : Opackowanie recyklingowe DK : Genanvendelig pakke FR : Emballage recyclable AR : قبة لإعادة التدوير FA : بسته قابل بازیافت	EN : Store below 25°C ES : Conservar a una temperatura inferior a 25°C PT : Armazenar abaixo de 25°C IT : Conservare ad una temperatura inferiore al 25°C DE : Unter 25°C lagern NL : Bewaren bij een temperatuur onder 25°C TR : 25°C'nin altında muhafaza ediniz GR : Φυλάξτε σε θερμοκρασία κάτω των 25°C PL : Przechowywać w temperaturze poniżej 25°C DK : Opbevareres under 25°C FR : Conserver sous 25°C AR : تحفظ تحت 25 درجة مئوية FA : در دما پایین تر از 25 درجه و سستی گرد نگهداری شود

EN : Authorized representative in the European Community ES : El representante competente en el Unión Europea PT : Representante autorizado na União Europeia IT : Rappresentante autorizzato nell'Unione Europea DE : Bevollmächtigter EU-Representant NL : De bevoegde vertegenwoordiger bij de Europese Gemeenschap TR : Avrupa Birliği'ni temsil eden yetkili temsilci GR : Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Ένωση PL : Upoważniony Przedstawiciel w Unii Europejskiej DK : Autoriseret Repræsentant i Den Europæiske Union FR : Délégué autorisé dans l'union européenne AR : الممثل المعتمد في الاتحاد الأوروبي FA : نماینده مجاز در اتحادیه اروپا	EN : Do not use if package is damaged ES : No utilizar si el paquete está dañado PT : Não usar se a embalagem estiver danificada IT : Non usare se la confezione è danneggiata DE : Bitte nicht verwenden, falls die Verpackung beschädigt ist NL : Gebruik geen als de verpakking beschadigd is TR : Paket zarar görülmüşse kullanmayınız GR : Μην χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά PL : Nie używać jeśli opakowanie jest uszkodzone DK : Brug ikke, hvis emballagen er beskadiget FR : Ne pas utiliser si l'emballage est endommagé AR : لا تستخدم إذا كانت العبوة التالفة FA : اگر بسته آسیب دیده باشد استفاده نکنید	EN : Sterile R: Gamma Irradiation ES : Esterilizado mediante irradiación: Irradiación gamma PT : Estéril R: Radiação Gamma IT : Sterile R: Raggi Gamma DE : ChargennummerSteril R: Gammastrahlen NL : Sterilisatie methode: Bestraling TR : Steril R: Radyasyon GR : Αποστείρωση PL : Sterylizowany przez napromieniowanie DK : Sterile R: Gamma Radiation FR : Stérile R : Radiation AR : معقم بالإشعاع FA : استریل‌شده با تابش	EN : Sterile R: Gamma Irradiation ES : Esterilizado mediante irradiación: Irradiación gamma PT : Estéril R: Radiação Gamma IT : Sterile R: Raggi Gamma DE : ChargennummerSteril R: Gammastrahlen NL : Sterilisatie methode: Bestraling TR : Steril R: Radyasyon GR : Αποστείρωση PL : Sterylizowany przez napromieniowanie DK : Sterile R: Gamma Radiation FR : Stérile R : Radiation AR : معقم بالإشعاع FA : استریل‌شده با تابش

IU-ZC-rev11/26.05.2022
Issue Date:11.09.2012

ZEOCEL®

Absorbable Haemostat, Oxidized Cellulose (OC)

EN Instructions for Use

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TR Kullanma Talimatı

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DESCRIPTION

ZEOCEL Absorbable Haemostat is a woven sterile hemostatic agent formed of oxidized cellulose (OC) (polyanhydro glucuronic acid). It has been manufactured from alpha grade cotton. Agent is a pale yellow material and has a caramel like smell. As it is stored in controlled room temperature, there might be a slight change in color in time but this does not effect the performance of the material. ZEOCEL Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Cellulose". ZEOCEL Absorbable Haemostat does not contain any animal or collagen additives.

ACTIONS

ZEOCEL Absorbable Haemostat can be applied to bleeding surface in different forms. The haemostatic agent material can be closed exactly, pressured against or sutured to the desired usage area. ZEOCEL Absorbable Haemostat should be used as dry. Action mechanism of OC ZEOCEL Absorbable Haemostats is independent from blood coagulation mechanism of the body. ZEOCEL Absorbable Haemostats form a brownish or black gelatinous mass which aids clot formation in 3-4 minutes when it contacts with blood. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. When ZEOCEL Absorbable Haemostat is used only as needed, it will be absorbed from the site of usage without tissue reaction. Absorption of the haemostatic agent from the used site depends on the amount of haemostatic agent used, tissue type and blood saturation degree. The material will be fully absorbed in 7-14 days. ZEOCEL Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr(-) (gram negative) microorganisms including those of:

Staphylococcus aureus	Escherichia coli	Klebsiella pneumoniae
Meticillin Resistant Staphylococcus epidermidis (MRSE)	Salmonella typhimurium	Klebsiella aerogenes
Meticillin Resistant Staphylococcus aureus (MRSA)	Proteus mirabilis	Bacillus subtilis
Staphylococcus epidermidis	Proteus vulgaris	Branhamella catarrhalis
Micrococcus luteus	Shigella boydii	Serratia Marcescens
Vancomycin Resistant Enterococcus sp. (VRE)	Shigella flexnerii	Clostridium perfringens
Penicilline Resistant Streptococcus pneumoniae (PRSP)	Shigella dysenteriae	Clostridium tetani
Enterococcus faecalis	Salmonella enteritidis	Listeria monocytogenes
Enterococcus faecium	Bacteroides fragilis	Bordetella pertussis
Enterococcus sp.	Lactobacillus casei	MOT (Mycobacterium Other than Tuberculosis)
Streptococcus pyogenes Group A	Lactobacillus sp.	Mycobacterium phlei
Streptococcus pyogenes Group B	Pseudomonas aeruginosa	Streptococcus faecalis
Streptococcus salivarius	Pseudomonas stutzeri	
Corynebacterium xerosis	Enterobacter cloacae	

INDICATIONS

ZEOCEL Absorbable Haemostat is designed to control the oozing of the blood, other different types of capillary bleeding and bleeding due to resection of the tissue during any type of surgical interventions. ZEOCEL Absorbable Haemostat can be used when other methods of haemostatic procedures are not usable. ZEOCEL Absorbable Haemostat can be used in many types of surgeries including but not limited to general surgery, digestive surgery, plastic surgery, orthopedy, gynaecology, urology, stomatology, traumatology and all other branches of surgery. Haemostat can be used to fill cavities formed by removal of tumors during surgeries as well as endoscopic interventions and dental praxis.

CONTRAINDICATIONS

ZEOCEL Absorbable Haemostat should not be used for implantation in bone defects, such as fractures since it may interfere with callus formation and may cause cyst formation. ZEOCEL Absorbable Haemostat should not be used for packing or wrapping. It should be removed after hemostasis is achieved if you have to use in this manner. ZEOCEL Absorbable Haemostat should not be used to control haemorrhage from large arteries. Since body fluids, except whole blood, such as serum don't react with ZEOCEL Absorbable Haemostat, it should not be used to produce hemostatic effect on non-hemostatic serous oozing surface. ZEOCEL Absorbable Haemostats should not be used as an adhesion prevention product.

WARNINGS

ZEOCEL Absorbable Haemostat is suitable for gamma irradiation and supplied as sterile. Should not be re-sterilized. ZEOCEL Absorbable Haemostat is not compatible with autoclaving or ethylene oxide sterilization. Usage of ZEOCEL Absorbable Haemostat in a contaminated wound without drainage can cause complications and contaminated areas should be cleaned before usage. Hemostatic properties of ZEOCEL Absorbable Haemostat are at best if applied dry; therefore it is not recommended to wet or moisten ZEOCEL Absorbable Haemostat before usage with water or saline.

ZEOCEL Absorbable Haemostat should not be impregnated with anti-infective agents, buffering materials or other hemostatic agents. Since thrombin activity is destroyed due to the low pH of the oxidized cellulose haemostats, haemostatic effect of ZEOCEL Absorbable Haemostat isn't enhanced by addition of thrombin. When necessary ZEOCEL Absorbable Haemostat can be left within the tissue. It is recommended to remove ZEOCEL Absorbable Haemostat when hemostasis is achieved. It must always be removed from

the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure, because ZEOCEL Absorbable Haemostat may exert pressure resulting in paralysis and/or nerve damage by swelling. Special care must be taken by physicians, regardless of the type of surgical procedure, to consider the advisability of removing ZEOCEL Absorbable Haemostat after hemostasis is achieved. Although ZEOCEL Absorbable Haemostat is bactericidal against a wide range of pathogenic microorganisms, it should not be used as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections. Attention should be paid when hemostasis stay in the body after surgical operations, in case of severely bleedings, although it is recommended to remove it after hemostasis.

PRECAUTIONS

Use only as much ZEOCEL Absorbable Haemostat as is necessary for haemostasis, holding it firmly in place until bleeding stops. Any excess ZEOCEL Absorbable Haemostat should be removed before the closure of the operated area in order to provide efficient absorption of the material and minimize the risk of foreign body reaction. In urological procedures, minimal amounts of ZEOCEL Absorbable Haemostat should be used and care must be exercised to prevent the plugging of the urethra, ureter, or a catheter by dislodged portions of the product. Since absorption of ZEOCEL Absorbable Haemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals. If ZEOCEL Absorbable Haemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped. Precautions should be taken in otolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling haemorrhage after tonsillectomy and controlling epistaxis). Care should be taken not to apply ZEOCEL Absorbable Haemostat too tightly when used as a wrap during vascular surgery.

ADVERSE REACTIONS

"Encapsulation" of fluid and foreign body reactions, hemorrhage, leucocytosis and haematoma have been reported in some cases when the residue of hemostat stay in the body after hemostasis achieved. Abscess of organs and wound dehiscence have been observed in some cases, even the cause of these cases are not related with Haemostat. Paralysis and nerve damage have been reported when Oxidized Cellulose (OC) Haemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures.

Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy has been reported. Occasional reports of "burning" and "stinging" sensation and sneezing when Oxidized Cellulose (OC) Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product. Headache, burning, stinging and sneezing in epistaxis control and other rhinological procedures, and stinging when Haemostat was applied on surface wounds (varicose ulcerations, demabrasions, and donor sites) also might be possible, rarely.

DOSAGE AND ADMINISTRATION

Sterile technique should be observed in removing ZEOCEL Absorbable Haemostat from its sterile container. Minimal amount of ZEOCEL Absorbable Haemostat in appropriate size is laid on the bleeding site or held firmly against the tissues until hemostasis is obtained. Required Haemostat amount depends on the nature and intensity of the hemorrhage to be stopped. The haemostatic effect of ZEOCEL Absorbable Haemostat is particularly pronounced when used dry. Moistening the material with water or physiological saline solution is not recommended.

CAUTION

Since ZEOCEL Absorbable Haemostat can not be re-sterilized, opened and unused ZEOCEL Absorbable Haemostat should be discarded. Dispose of contaminated or unused products in accordance with local and facility requirements.

NOTES AND INFORMATION ON SHELF-LIFE

ZEOCEL Absorbable Haemostat should be stored dry under controlled conditions (max. 25°C) and protected from direct sunlight, in the original packaging. Shelf life of the product is 3 years. The expiry date ZEOCEL Absorbable Haemostat is printed on the pack. Do not use ZEOCEL Absorbable Haemostat after this date.

PRESENTATION AND CONTENTS

ZEOCEL Absorbable Haemostat is available sterile, as an original pack in different sizes. ZEOCEL Absorbable Haemostat strips are sterile-packed individually.

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