



www.boztibbi.com
boz@boztibbi.com

Product Catalogue



Boz Medical Device Industry and Trade Inc.

Product
Catalog Link



www.boztibbi.com/bozkatalog.pdf

About our Company

Boz Medical Device Industry and Trade Inc. is founded by Adil Boz and Ali Aykut Boz in August of 2009 in Ankara. Common grandfather Adil Boz of the founders are one of the oldest people in health sector in Turkey.

Adil Boz (Grandfather) began his trade life in 1940s. He founded a small family company, began to produce the gauze in 1960s, plaster of paris bandage and absorbent cotton in 1970's to the health sector. He is the one of oldest medical device producer of Türkiye. The company founded by Adil Boz located in Uşak, produced integratedly from raw cotton to the yarn; from yarns to the woven; woven textile raw materials, home and hospital textile fabrics and provided job oppornuties to 1500 people in its facilities.

Adil Boz founded a company that gives his own name in 1983, and has begun to his trade life with selling product of gauze, plaster of paris bandage and absorbent cotton that is produced by Adil Boz in Uşak. After noticing lack of the surgical suture goods, in 1993 he has begun to trial production after set up factory for the purpose of producing surgical sutures with the name of BOZ MEDICAL DEVICE CO. In year of 1998 he manufactured the products of gauze compress, sponge, cotton pad, gauze roll that are produced from gauze. Boz family, since they have gone into the bussine life in Republic of Turkiye, created a competitive enviroment against the a lot of international companies and has played a important role in our natioal interests.

In 2009, Boz Medical Device Industry and Trade Inc. founded by Ali Aykut Boz and Adil Boz, has aimed in accordance with emerging market needs and conditions to be in position of manufacturer and has achived its goal in month of April of 2012.

The factory established in Ankara, has 8600 m2 closed areas and 1200 m2 clean rooms. In the facility has equipped with modern climate system, the white areas are the level of class III. The minimum requirements in production, is executed prepared towards the World Health Organization's (WHO) proposed principles and published "Good Production Applications" by T.C. Ministry of Health. Rigor and Hygenic are essential at all level of Production.



Needle Information & Use of Surgical Needles **1**

Tip Shapes According to Lengths	2
Needle Point Profiles	3
Needle Shapes According to Lengths	4 - 8
Angle of Needles & Needle Shape	9
Comparision of Suture Size	10
Use of Needle Holders	11

Packaging & Labelling Information **12-13**

Suture Labelling Guide	14
Packaging & Labelling Information	15
Multi - Lingual Symbols	16
Package Types	17
Comparision of Suture Size	18

Absorbable Sutures **19**

Glikosorb Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided	20-33
Glikolak Surgical Suture, Poly[glycolide(90%-co-lactide(10%)](PGLA) Absorbable, Synthetic, Braided	34-43
Fastsorb Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed	44-49
Fastlak Surgical Suture, Poly[glycolide(90%-co-lactide(10%)](PGLA) Fast Absorbable, Synthetic, Braided, Undyed	50-55
Monokaprol Surgical Suture, Poly[glycolide(75%-co-caprolactone (25%)](PGCL) Absorbable, Synthetic, Monofilament	56-61
Polidiox Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet	62-67

Antibacterial Absorbable Sutures **68**

Glikosorb Plus Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided	69-77
Glikolak Plus Surgical Suture, Poly[glycolide(90%-co-lactide(10%)](PGLA) Absorbable, Synthetic, Braided	78-86
Monokaprol Plus Surgical Suture, Poly[glycolide(75%-co-caprolactone(25%)](PGCL) Absorbable, Synthetic, Monofilament	87-92
Polidiox Plus Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament	93-96

Contents

Absorbable Sutures 97

Plain Catgut Absorbable Surgical Suture 98-106

Chromic Catgut Absorbable Surgical Suture 107-116

Non - Absorbable Sutures 117

Silk Surgical Suture, Silk, Silicone Coated Nonabsorbable, Noncapillary 118-128

Monoprolen Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament 129-141

Monoflorid Surgical Suture, Poly (vinylidene difluoride) (PVDF) Nonabsorbable, Synthetic, Monofilament 142-147

Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic 148-154

Multamid Surgical Suture, Multifilament Polyamide (PA) Silicone Coated Nonabsorbable, Synthetic 155-157

Polisil Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic 158-165

Monotef Surgical Suture, Polytetrafluoroethylene (PTFE) Nonabsorbable, Synthetic, Monofilament 166-169

Polytef Surgical Suture, Braided Polyester, (PTFE) Coated Nonabsorbable, Synthetic 170-171

Monowire Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament 172-174

Tempowire Surgical Stainless Steel Wire Nonabsorbable, Multifilament 175-177

PTFE pledget 178-179

Absorbable Haemostats 180

Reoxcel Oxidized Regenerated Cellulose (ORC) 181

Reoxcel - Sigma Knit Oxidized Regenerated Cellulose (ORC) 182


Reoxcel - Fibril Oxidized Regenerated Cellulose (ORC) 183

Zeocel Oxidized Cellulose (OC) 184

Specialist Products	185
Microsurgery Products Polyethylene	186
Oval Retention Plate Polypropylene	187
Buttressing Button Polypropylene	188
Locking Button Polypropylene Bead	189
Aluminum Ring	190


Needle Information & Use of Surgical Needles

Taper Point
Square Body, Taper Point
Blunt
Reverse Cutting
Inside Cutting
Diamond
Spatula
Lancet
CC Needle
Tapercutting
Trocars
Square Body Tapercutting




Tip Shapes According to Lengths

Round Body, Taper Point




Doesn't have a cutting edge. Has a conic body becoming sharp towards its tip. Is used in tissues showing low resistance to the needle.

Square Body, Taper Point




Doesn't have a cutting edge. Has a conic body becoming sharp towards its tip. In order to strengthen the body of the needle, the body is given a square shape.

Blunt




Doesn't have a property of cutting and drilling. Is used in parenchyma tissues that can be easily damaged like liver and kidney.

Reverse Cutting




Reverse (Outside apex) cutting needle has cutting edges at the bottom and in its two sides, becomes sharp towards its tip and has triangular point.

Inside Cutting




Inside (Inside apex) cutting needle has cutting edges at the top and in its two sides, becomes sharp towards its tip and has triangular point.

Diamond




This needle has a keen tip. It angled conically towards the needle body. It's generally used in hard tissues.

Spatula




This needle is flat and become thinner towards its tip. It's especially used in ophthalmic surgery.

Lancet




Its structure looks like spatulated needle. However, its tip isn't parallel but oblique to its body. Edges of its top are sharp. It's generally used in ophthalmic surgery.

CC Needle




It combines the properties of round and cutting needles. Three cutting edges of the needle does not combine with the round body. Sharp edges have single cutting type.

Tapercutting




Combine the features of the reverse cutting edge tip and taper point needles. Three cutting edges extend from the point. These blend into a round taper body. All three edges are sharpened to provide uniform cutting action.

Trocar






Trocar needle has a strong cutting tip. The tip of the needle have been designed for penetration of the even dense tissues.

Square Body Tapercutting



In order to strengthen the body of the tapercutting needle, the body is given a square shape.

Needle Point Profiles


Taper Point		Doesn't have a cutting edge. Has a conic body becoming sharp towards its tip. Is used in tissues showing low resistance to the needle.
Square Body, Taper Point		Doesn't have a cutting edge. Has a conic body becoming sharp towards its tip. In order to strengthen the body of the needle, the body is given a square shape.
Blunt		Doesn't have a property of cutting and drilling. Is used in parenchyma tissues that can be easily damaged like liver and kidney.
Reverse Cutting		Reverse (Outside apex) cutting needle has cutting edges at the bottom and in its two sides, becomes sharp towards its tip and has triangular point.
	 Premium	This reverse cutting needle is handhoned and has a long tip. Generally it's used in plastic surgery.
	 Premium Thin Line	Compared to premium needle, it's more grinded and sharp. Generally it's used in plastic surgery.
Inside Cutting		Inside (Inside apex) cutting needle has cutting edges at the top and in its two sides, becomes sharp towards its tip and has triangular point.
	 Premium	This cutting edge needle is handhoned and has a long tip. Generally it's used in plastic surgery.
	 Premium Thin Line	Compared to premium needle, it's more grinded and sharp. Generally it's used in plastic surgery.
Diamond		This needle has a keen tip. It angled conically towards the needle body. It's generally used in hard tissues.
	 Premium	It's handhoned.
Spatula		This needle is flat and become thinner towards its tip. It's especially used in ophthalmic surgery.
Lancet		It's structure looks like spatulated needle. However, its tip isn't parallel but oblique to its body. Edges of its top are sharp. It's generally used in ophthalmic surgery.
CC Needle		It combines the properties of round and cutting needles. Three cutting edges of the needle does not combine with the round body. Sharp edges have single cutting type.
Tapercutting		Combine the features of the reverse cutting edge tip and taper point needles. Three cutting edges extend from the point. These blend into a round taper body. All three edges are sharpened to provide uniform cutting action.
Trocar		Trocar needle has a strong cutting tip. The tip of the needle have been designed for penetration of the even dense tissues.
Square Body Tapercutting		In order to strengthen the body of the tapercutting needle, the body is given a square shape.

Needle Shapes According to Lengths


Taper Point




1/2




3/8



5/8



Düz /
Straigth




Needle Shapes According to Lengths


Square Body, Taper Point




1/2




3/8



5/8





Düz /
Straigh




Needle Shapes According to Lengths



Reverse Cutting



Düz / Straighth



Inside Cutting




Needle Shapes According to Lengths


Diamond



1/2




3/8



Düz /
Straigh

25



Premium


3/8



Blunt





3/8



90

Needle Shapes According to Lengths


Tapercutting



Spatula




Angle of Needles & Needle Shape




The choice of needle shape is frequently governed by the accessibility of the tissue to be sutured, and normally the more confined the operative site the greater the curvature required. The basic shapes involved are:


1/4 Circle




3/8 Circle




1/2 Circle




5/8 Circle



J Shape



Straight



Comparison of Suture Size

A

USP values are stated for the non-absorbable mono and multiflament, absorbable multiflament sutures in the columns under A.


USP	Metric	Suture Diameter, mm	% reduction
3	6	 0.600-0.699	
2	5	 0.500-0.599	28
1	4	 0.400-0.499	33
0	3.5	 0.350-0.399	30
2/0	3	 0.300-0.349	25
3/0	2	 0.200-0.249	43
4/0	1.5	 0.150-0.199	40
5/0	1	 0.100-0.149	50
6/0	0.7	 0.070-0.099	53
7/0	0.5	 0.050-0.069	49
8/0	0.4	 0.040-0.049	42
9/0		 0.030-0.039	
10/0		 0.020-0.029	

B

USP values are stated for the monofilament absorbable sutures in the columns under B.

USP	Metric	Suture Diameter, mm	% reduction
3		 0.600-0.699	
2	5	 0.571-0.610	
1	4	 0.500-0.570	18
0	3.5	 0.400-0.499	30
2/0	3	 0.340-0.399	32
3/0	2	 0.250-0.339	37
4/0	1.5	 0.200-0.249	35
5/0	1	 0.150-0.199	40
6/0	0.7	 0.095-0.149	52
7/0	0.5	 0.050-0.094	66
8/0		 0.030-0.039	
9/0		 0.020-0.029	
10/0		 0.010-0.019	

Use of Needle Holders




1- The needle holder should be carefully selected to match the size and strength of the needle to be used. The use of a needle holder larger than required can result in damage to the needle and in particular, distortion of the curvature.


2- The needle holder should be in good condition as worn jaws can result in needle rotation and instability in the needle holder. Nicks or defects in the needle holder jaws can also cause damage to the needle and loss of strength.

3- Needles should be grasped securely on the tip of the needle holder jaws.

4- Needles should only be held in the flatted area provided and should not be grasped on the attachment area or near to the needle point. Needles which are not flatted should be grasped for placement at a point approximately one third of the needle length from the butt. Excessive force should not be applied when gripping the needle, particularly with Tungsten Carbide jaw inserts, as this may damage the needle body and cause loss of strength or breakage.



Use of Surgical Needles



1- The force required to achieve passage of the needle through tissue should be applied in a direction following the curvature of the needle.

2- Care should be taken to match the size of the needle to the size of tissue bite required. The use of too small a needle for a given tissue bite can lead to bending.


3- Should the placement of the needle in tissue require to be readjusted, the needle should be removed and re-inserted. No attempt should be made to twist the needle in tissue.

4- The needles should normally be inserted separately into each side of the tissue to be approximated and should not be used to bridge a wound.

Packaging & Labelling Information




Packaging & Labelling Information




Suture Labelling Guide

Box Label



Sachet Label



Packaging

Sutures are packed in one or two / three dozen boxes.
Boxes are designed for easy usage and environmental friendly have been.
Box labels are colour coded for easy suture type defining.
Pealable double blister is the individual package.
Blister consists of transparent polyethylene film and medical grade of paper.
Paper is coloured with Ethylene Oxide (EO) indicator; yellow colour proves the product processed in EO gas.

Labelling

The Following symbols(from European Standard En980) are used on our labelling

- On each product and box; type of the suture, USP and EP size, needle length and following signs are stated.

Certificates Labelling



CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC

Medical Device Directive Labelling

- The following symbols (from European Standard EN980) are used on our labelling.



Do not reuse.



Use until year and month.



Sterile Product-Method of Sterilization-Ethylene Oxide.



Sterile Product-Method of Sterilization-Irradiation.



Batch Number



Attention, See instruction for use



Marking of Conformity. Notified Body number.

Additional Labelling

- For reordering products are defined by code numbers.



Product Catalogue (Code) Number.

- To supply easy informing needle shape is clearly defined by drawings.



Needlepoint slide view, cutting edge, taper etc.




Three dimensional view of the needle point.



Needle shape: 1/2 c. etc. and; length. Double needle shows double armed.

Multi - Lingual Symbols

The European Association of Surgical Suture Industry has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive and enables companies to provide information to the customer without having to provide multilingual translations.

Symbol	Description	Example Product
	Dyed, Absorbable, Braided, Coated	GLIKOLAK, GLIKOSORB, <i>Antibacterial</i> GLIKOLAK Plus, <i>Antibacterial</i> GLIKOSORB Plus
	Dyed, Absorbable, Monofilament	MONOKAPROL, POLIDIOX, <i>Antibacterial</i> MONOKAPROL Plus, <i>Antibacterial</i> POLIDIOX Plus
	Dyed, Non-absorbable, Braided, Coated	SILK, MULTAMID, POLISIL, POLYTEF
	Dyed, Non-absorbable, Monofilament	MONOPROLEN, MONOFLORID, MONAMID
	Dyed, Non-absorbable, Twisted, Coated	SILK
	Undyed, Absorbable, Braided, Coated	FASTLAK, FASTSORB, GLIKOLAK, <i>Antibacterial</i> GLIKOLAK Plus
	Undyed, Absorbable, Monofilament	MONOKAPROL, <i>Antibacterial</i> MONOKAPROL Plus
	Undyed, Non-absorbable, Braided, Coated	SILK
	Undyed, Non-absorbable, Monofilament	MONAMID, MONOTEF, MONOWIRE
	Undyed, Non-absorbable, Twisted, Coated	SILK, TEMPOWIRE
	Undyed, Non-absorbable Twisted	

Dozen Box



Standard individual package is a dozen box.
Twelve of sachets are put in drawer box, than the box is labelled.



Double or triple dozen boxes are used for packaging of 24/36 sachets on request.

Long Boxes



To avoid bending, monofilament stainless steel wires are put in individual envelopes which are as long as suture length. Than they are boxed in long boxes at the same size.

Shipping Container



Depending on order quantity boxes are put in containers which can store 40/60 or 100 standard boxes.

All boxes are preshrunk before packaging.



Comparison of Suture Size

Box Label

ORDER NO : ①		###		
COMPANY NAME : ②		XXX		
CITY : ③		XXX		
USP	CATALOG NO	TYPE	LOT NO. ④	AMOUNT
2/0	GS3858	GLIKOSORB 30 mm	1705017	100 ⑤
—	—	—	—	—
				100 ⑥
⑦ PERFORMER				

Each parcel is labelled with a sticker. This label defines the integrity of the parcel including order references.

- ① Order number given by manufacturer, refer to order history record.
- ② Company who ordered.
- ③ Destination: city if inland, country if abroad.
- ④ USP/Catalogue/Lot numbers and Type of the enclosed product/products.
- ⑤ Amount of each item enclosed, in dozen.
- ⑥ Total amount of enclosed products, in dozen.
- ⑦ Person who prepared parcel.



Boz Medical Device Industry and Trade Inc.

Absorbable Sutures

GLIKOSORB

GLIKOLAK

FASTSORB

FASTLAK

MONOKAPROL

POLIDIOX

Surgical Sutures Absorbable

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Boz®



Definition

Synthetic absorbable, coated, mid term wound support and mid term mass absorption.

Material

Polyglycolic Acid (PGA)

Coating

Calcium Stearate and polycaprolactone

Structure

Braided - Multifilament

Colour

Violet (dyed) or Beige (undyed)

Size

EP : 0.4 to 8

USP : 8/0 to 6

Tensile Strength Retention

2. Week 75 % - 3. Week 50 %

Mass Absorption

60-90 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery

Gynaecology

Ophthalmic Surgery

Orthopedics

Urology

Skin closer - intracutaneous, subcutaneous.

Ligatures








Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
6 mm							
3/8 Spatula 220 Micron			GS5401	8/0 0,4	45	12	Violet
3/8 Spatula 230 Micron			GS1156	8/0 0,4	45	12	Violet
			GS1252	7/0 0,5	30	12	Violet
6,6 mm							
1/2 Cobra 230 Micron			GS5519	8/0 0,4	60	12	Violet
1/2 Cobra 280 Micron			GS5518	7/0 0,5	60	12	Violet
8 mm							
1/2 Spatula 230 Micron			GS1193	8/0 0,4	30	12	Violet
			GS1195	8/0 0,4	45	12	Violet
1/2 Spatula 330 Micron			GS1272	7/0 0,5	45	12	Violet
			GS1395	6/0 0,7	45	12	Violet
1/2 Spatula 340 Micron			GS5526	5/0 1,0	45	12	Violet
1/4 Spatula 230 Micron			GS1190	8/0 0,4	45	12	Violet
1/4 Spatula 250 Micron			GS5876	7/0 0,5	45	12	Violet
1/4 Spatula 340 Micron			GS1268	7/0 0,5	45	12	Violet
			GS1404	6/0 0,7	30	12	Violet
			GS1405	6/0 0,7	45	12	Violet
			GS1406	6/0 0,7	45	12	Undyed
			GS1764	5/0 1,0	45	12	Violet











Glikosorb


Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 430 Micron			GS2362	4/0 1,5	45	12	Violet
1/4 Spatula 430 Micron			GS2363	4/0 1,5	45	12	Violet
3/8 Round Body, Taper Point 230 Micron			GS1207	8/0 0,4	30	12	Violet
3/8 Round Body, Taper Point 280 Micron			GS1413	6/0 0,7	45	12	Violet
			GS1426	6/0 0,7	45	12	Violet
			GS1768	5/0 1,0	45	12	Violet
10 mm							
3/8 Round Body, Taper Point 280 Micron			GS1480	6/0 0,7	75	12	Violet
3/8 Reverse Cutting 330 Micron			GS1496	6/0 0,7	45	12	Violet
3/8 Reverse Cutting 330 Micron			GS1492	6/0 0,7	45	12	Violet
3/8 Inside Cutting 280 Micron			GS1508	6/0 0,7	45	12	Violet
12 mm							
3/8 Reverse Cutting 330 Micron			GS1541	6/0 0,7	45	12	Violet
3/8 Spatula 280 Micron			GS1523	6/0 0,7	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 330 Micron			GS1560	6/0 0,7	45	12	Violet
			GS1564	6/0 0,7	75	12	Violet

Glikosorb








Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 330 Micron			GS1561	6/0 0,7	45	12	Violet
1/2 Round Body, Taper Point 380 Micron			GS1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			GS1894	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 330 Micron			GS1580	6/0 0,7	45	12	Violet
			GS1590	6/0 0,7	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			GS1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GS2432	4/0 1,5	45	12	Violet
			GS2434	4/0 1,5	60	12	Violet
3/8 Reverse Cutting 380 Micron			GS1928	5/0 1,0	75	12	Violet
3/8 Reverse Cutting 560 Micron			GS2456	4/0 1,5	75	12	Violet
3/8 Spatula 280 Micron			GS1605	6/0 0,7	45	12	Violet
3/8 Spatula 330 Micron			GS1873	5/0 1,0	60	12	Violet
3/8 Trocarn 330 Micron			GS2416	4/0 1,5	45	12	Violet

16 mm							
1/2 Round Body, Taper Point			GS1990	5/0 1,0	75	12	Violet
			GS1992	5/0 1,0	75	12	Undyed
			GS2503	4/0 1,5	75	12	Undyed
			GS2504	4/0 1,5	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			GS3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GS1674	6/0 0,7	75	12	Violet
			GS2014	5/0 1,0	45	12	Undyed
			GS2540	4/0 1,5	75	12	Violet
			GS3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GS2140	5/0 1,0	45	12	Violet
			GS2144	5/0 1,0	75	12	Violet
			GS2557	4/0 1,5	75	12	Undyed
			GS3039	3/0 2,0	75	12	Violet
17 mm							
1/2 Round Body, Taper Point			GS2584	4/0 1,5	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GS5525	5/0 1,0	75	12	Violet
			GS2613	4/0 1,5	75	12	Violet
			GS3624	2/0 3,0	75	12	Violet
3/8 Round Body, Taper Point			GS3085	3/0 2,0	75	12	Violet
19 mm							
3/8 Reverse Cutting			GS1706	6/0 0,7	75	12	Violet
			GS2640	4/0 1,5	75	12	Violet
3/8 Reverse Cutting Premium			GS2211	5/0 1,0	45	12	Violet


Glikosorb




Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			GS2229	5/0 1,0	75	12	Violet
			GS2685	4/0 1,5	75	12	Violet
			GS3142	3/0 2,0	75	12	Violet
1/2 Tapercutting			GS2673	4/0 1,5	75	12	Violet
3/8 Reverse Cutting			GS2250	5/0 1,0	75	12	Undyed
			GS2716	4/0 1,5	45	12	Undyed
			GS2719	4/0 1,5	75	12	Violet
			GS3164	3/0 2,0	45	12	Violet
			GS3168	3/0 2,0	75	12	Violet
			GS3663	2/0 3,0	75	12	Violet
22 mm							
1/2 Round Body, Taper Point			GS2747	4/0 1,5	75	12	Violet
			GS3196	3/0 2,0	75	12	Violet
			GS3671	2/0 3,0	75	12	Violet
			GS4575	1 4,0	75	12	Violet
3/8 Reverse Cutting			GS2763	4/0 1,5	75	12	Violet
			GS3214	3/0 2,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided





Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
24 mm					
1/2 Reverse Cutting  	GS5524	4/0 1,5	75	12	Violet
	GS5523	3/0 2,0	75	12	Violet
	GS3692	2/0 3,0	75	12	Violet

25 mm					
1/2 Round Body, Taper Point  	GS2788	4/0 1,5	75	12	Violet
	GS3252	3/0 2,0	75	12	Violet
	GS3746	2/0 3,0	100	12	Violet
	GS3724	2/0 3,0	75	12	Violet
1/2 Tapercutting  	GS3718	2/0 3,0	75	12	Violet
3/8 Reverse Cutting  	GS2804	4/0 1,5	75	12	Violet
	GS3757	2/0 3,0	75	12	Violet
3/8 Inside Cutting  	GS2812	4/0 1,5	45	12	Violet
	GS3284	3/0 2,0	45	12	Violet
	GS3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point  	GS2292	5/0 1,0	75	12	Violet
	GS2828	4/0 1,5	75	12	Violet
	GS3313	3/0 2,0	75	12	Violet
	GS3318	3/0 2,0	90	12	Violet
	GS3786	2/0 3,0	75	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point  	GS4215	0 3,5	75	12	Violet
	GS4592	1 4,0	90	12	Violet
1/2 Reverse Cutting  	GS4230	0 3,5	75	12	Violet
3/8 Round Body, Taper Point  	GS2841	4/0 1,5	75	12	Violet
3/8 Reverse Cutting  	GS3345	3/0 2,0	75	12	Undyed
	GS3346	3/0 2,0	75	12	Violet
	GS3821	2/0 3,0	75	12	Undyed
	GS3822	2/0 3,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point  	GS3369	3/0 2,0	75	12	Violet
	GS3374	3/0 2,0	90	12	Violet
	GS3858	2/0 3,0	75	12	Violet
	GS4244	0 3,5	75	12	Violet
	GS4608	1 4,0	75	12	Violet
1/2 Reverse Cutting  	GS3882	2/0 3,0	75	12	Violet
	GS4266	0 3,5	75	12	Violet
1/2 Tapercutting  	GS3843	2/0 3,0	75	12	Violet


Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting 		GS2887	4/0 1,5	75	12	Violet
		GS3397	3/0 2,0	75	12	Violet
		GS3891	2/0 3,0	75	12	Violet
		GS4273	0 3,5	75	12	Violet
32 mm						
1/2 Round Body, Taper Point 		GS5508	3/0 2,0	75	12	Violet
35 mm						
1/2 Round Body, Taper Point 		GS3411	3/0 2,0	75	12	Violet
		GS3916	2/0 3,0	75	12	Violet
		GS3922	2/0 3,0	90	12	Violet
		GS4300	0 3,5	75	12	Violet
		GS4644	1 4,0	75	12	Violet
3/8 Reverse Cutting 		GS3423	3/0 2,0	75	12	Violet
		GS3933	2/0 3,0	75	12	Violet
36 mm						
1/2 Round Body, Taper Point 		GS5509	3/0 2,0	75	12	Violet
1/2 Reverse Cutting 		GS4675	1 4,0	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
37 mm					
1/2 Round Body, Taper Point 	GS3950	2/0 3,0	75	12	Violet
	GS4337	0 3,5	75	12	Violet
	GS4681	1 4,0	75	12	Violet
1/2 Reverse Cutting 	GS5507	2/0 3,0	75	12	Violet
	GS4354	0 3,5	90	12	Violet
	GS4994	2 5,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	GS3982	2/0 3,0	75	12	Violet
	GS4368	0 3,5	75	12	Violet
	GS4375	0 3,5	90	12	Violet
	GS4713	1 4,0	75	12	Violet
	GS4719	1 4,0	90	12	Violet
	GS5003	2 5,0	75	12	Violet
	GS5006	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GS4397	0 3,5	75	12	Violet
	GS4745	1 4,0	75	12	Violet
	GS4749	1 4,0	90	12	Violet
	GS5019	2 5,0	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm						
1/2 Reverse Cutting		GS5023	2 5,0	90	12	Violet
1/2 Tapercutting		GS4363	0 3,5	75	12	Violet
3/8 Reverse Cutting		GS4756	1 4,0	75	12	Violet
		GS4408	0 3,5	75	12	Violet
45 mm						
1/2 Round Body,Taper Point		GS4422	0 3,5	75	12	Violet
		GS4772	1 4,0	75	12	Violet
		GS5035	2 5,0	75	12	Violet
		GS5040	2 5,0	90	12	Violet
1/2 Tapercutting		GS4766	1 4,0	75	12	Violet
		GS5738	2 5,0	75	12	Violet
3/8 Reverse Cutting		GS4432	0 3,5	75	12	Violet
		GS4787	1 4,0	75	12	Violet
Straight Reverse Cutting		GS2909	4/0 1,5	45	12	Violet
48 mm						
1/2 Round Body,Taper Point		GS4796	1 4,0	75	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm						
1/2 Round Body, Taper Point 		GS4454	0 3,5	75	12	Violet
		GS4459	0 3,5	90	12	Violet
		GS4812	1 4,0	75	12	Violet
		GS4818	1 4,0	90	12	Violet
		GS5066	2 5,0	75	12	Violet
		GS5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 		GS4471	0 3,5	75	12	Violet
60 mm						
Straight Reverse Cutting 		GS3492	3/0 2,0	75	12	Violet
		GS4074	2/0 3,0	75	12	Violet
63 mm						
3/8 Blunt 		GS4858	1 4,0	75	12	Violet
64 mm						
3/8 Blunt 		GS4500	0 3,5	90	12	Violet
		GS4861	1 4,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
3x45 cm					
GLIKOSORB Non-Needle Suture	GS2335	4/0 1,5	3x45	12	Violet

6x45 cm					
GLIKOSORB Non-Needle Suture	GS2950	3/0 2,0	6x45	12	Violet
	GS3543	2/0 3,0	6x45	12	Violet
	GS4137	0 3,5	6x45	12	Violet
	GS4563	1 4,0	6x45	12	Violet

10x50 cm					
GLIKOSORB Non-Needle Suture	GS2954	3/0 2,0	10x50	12	Violet
	GS3548	2/0 3,0	10x75	12	Violet

12x45 cm					
GLIKOSORB Non-Needle Suture	GS2344	4/0 1,5	12x45	12	Violet
	GS2958	3/0 2,0	12x45	12	Violet
	GS3552	2/0 3,0	12x45	12	Violet
	GS4145	0 3,5	12x45	12	Violet
	GS4554	1 4,0	12x45	12	Violet

17x45 cm					
GLIKOSORB Non-Needle Suture	GS2966	3/0 2,0	17x45	12	Violet
	GS3558	2/0 3,0	17x45	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
-------------------	------------	--------	-------------	----------	-------

150 cm

GLIKOSORB Non-Needle Suture	GS1737	5/0 1,0	150	12	Violet
	GS2326	4/0 1,5	150	12	Violet
	GS2935	3/0 2,0	150	12	Violet
	GS3511	2/0 3,0	150	12	Violet
	GS4114	0 3,5	150	12	Violet
	GS4532	1 4,0	150	12	Violet

250 cm

GLIKOSORB Non-Needle Suture	GS2332	4/0 1,5	250	12	Violet
	GS2942	3/0 2,0	250	12	Violet
	GS3518	2/0 3,0	250	12	Violet
	GS4121	0 3,5	250	12	Violet
	GS4541	1 4,0	250	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, coated, mid term wound support and mid term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Coating

Poly(glycolide-co-lactide) (30/70) & calcium stearate

Structure

Braided-Multifilament

Colour

Beige (undyed) or Violet (dyed)

Size

EP : 0.4 to 5

USP : 8/0 to 2

Tensile Strenght Retention

2. Week 75%

Mass Absorption

55-70 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery

Gynaecology

Ophthalmic Surgery

Orthopedics

Urology

Skin closer - intracutaneous, subcutaneous.

Ligatures










Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
6 mm							
3/8 Reverse Cutting 150 Micron			GL1165	8/0 0,4	45	12	Violet
3/8 Reverse Cutting 150 Micron			GL1256	7/0 0,5	45	12	Violet
3/8 Spatula 230 Micron			GL1155	8/0 0,4	30	12	Violet
3/8 Spatula 230 Micron			GL1253	7/0 0,5	45	12	Violet
8 mm							
1/2 Spatula 330 Micron			GL1394	6/0 0,7	30	12	Violet
1/4 Spatula 340 Micron			GL1405	6/0 0,7	45	12	Violet
1/4 Spatula 340 Micron			GL1404	6/0 0,7	30	12	Violet
			GL1764	5/0 1,0	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1299	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1432	6/0 0,7	75	12	Violet
3/8 Spatula 230 Micron			GL5827	8/0 0,4	45	12	Violet
3/8 Spatula 280 Micron			GL1413	6/0 0,7	45	12	Violet
10 mm							
3/8 Round Body, Taper Point 340 Micron			GL1803	5/0 1,0	45	12	Violet
3/8 Reverse Cutting 330 Micron			GL1496	6/0 0,7	45	12	Violet
3/8 Spatula 280 Micron			GL1450	6/0 0,7	45	12	Violet
3/8 Spatula 330 Micron			GL1790	5/0 1,0	75	12	Violet



Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
1/2 Round Body, Taper Point 430 Micron			GL1836	5/0 1,0	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 330 Micron			GL1560	6/0 0,7	45	12	Violet
			GL1564	6/0 0,7	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			GL1894	5/0 1,0	75	12	Violet
			GL1900	5/0 1,0	90	12	Violet
3/8 Round Body, Taper Point 330 Micron			GL1580	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 380 Micron			GL1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GL2434	4/0 1,5	60	12	Violet
3/8 Reverse Cutting Premium 330 Micron			GL1626	6/0 0,7	45	12	Violet
16 mm							
1/2 Round Body, Taper Point			GL1990	5/0 1,0	75	12	Violet
			GL1994	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point			GL2504	4/0 1,5	75	12	Violet
			GL3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GL2021	5/0 1,0	75	12	Violet
			GL2540	4/0 1,5	75	12	Violet
			GL3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GL2144	5/0 1,0	75	12	Violet


Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm					
3/8 Inside Cutting  	GL2557	4/0 1,5	75	12	Undyed
	GL2553	4/0 1,5	45	12	Violet
	GL3039	3/0 2,0	75	12	Violet
17 mm					
1/2 Round Body, Taper Point  	GL2165	5/0 1,0	75	12	Violet
18 mm					
1/2 Round Body, Taper Point  	GL2613	4/0 1,5	75	12	Violet
19 mm					
3/8 Reverse Cutting  	GL1706	6/0 0,7	75	12	Violet
	GL2206	5/0 1,0	75	12	Violet
	GL2640	4/0 1,5	75	12	Violet
	GL3109	3/0 2,0	75	12	Violet
20 mm					
1/2 Round Body, Taper Point  	GL2229	5/0 1,0	75	12	Violet
	GL2685	4/0 1,5	75	12	Violet
	GL2690	4/0 1,5	90	12	Violet
	GL3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting  	GL2716	4/0 1,5	45	12	Undyed
	GL3164	3/0 2,0	45	12	Violet
	GL3663	2/0 3,0	75	12	Violet


Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
25 mm							
1/2 Round Body, Taper Point			GL2788	4/0 1,5	75	12	Violet
			GL3252	3/0 2,0	75	12	Violet
			GL3746	2/0 3,0	100	12	Violet
			GL3724	2/0 3,0	75	12	Violet
25 mm							
3/8 Reverse Cutting			GL2804	4/0 1,5	75	12	Violet
3/8 Inside Cutting			GL2812	4/0 1,5	45	12	Violet
			GL3284	3/0 2,0	45	12	Violet
			GL3763	2/0 3,0	75	12	Violet
26 mm							
1/2 Round Body, Taper Point			GL3313	3/0 2,0	75	12	Violet
			GL4220	0 3,5	90	12	Violet
3/8 Reverse Cutting			GL2852	4/0 1,5	75	12	Violet
			GL3346	3/0 2,0	75	12	Violet
			GL3822	2/0 3,0	75	12	Violet
30 mm							
1/2 Round Body, Taper Point			GL3369	3/0 2,0	75	12	Violet
			GL3858	2/0 3,0	75	12	Violet
			GL3865	2/0 3,0	90	12	Violet




Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	GL4244	0 3,5	75	12	Violet
	GL4250	0 3,5	90	12	Violet
	GL4608	1 4,0	75	12	Violet
	GL4614	1 4,0	90	12	Violet
1/2 Round Body, Taper Point 	GL3854	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GL3397	3/0 2,0	75	12	Violet
	GL3891	2/0 3,0	75	12	Violet
	GL3895	2/0 3,0	90	12	Violet
	GL4273	0 3,5	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GL3916	2/0 3,0	75	12	Violet
	GL4300	0 3,5	75	12	Violet
	GL4644	1 4,0	75	12	Violet
	GL4649	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GL3423	3/0 2,0	75	12	Violet
	GL3933	2/0 3,0	75	12	Violet
	GL5881	1 4,0	75	12	Violet




Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	GL3982	2/0 3,0	75	12	Violet
	GL4368	0 3,5	75	12	Violet
	GL4375	0 3,5	90	12	Violet
	GL4713	1 4,0	75	12	Violet
	GL4719	1 4,0	90	12	Violet
	GL5003	2 5,0	75	12	Violet
1/2 Reverse Cutting 	GL4399	0 3,5	90	12	Violet
	GL4749	1 4,0	90	12	Violet
	GL5019	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4408	0 3,5	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GL4772	1 4,0	75	12	Violet
	GL4775	1 4,0	90	12	Violet
	GL5035	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4022	2/0 3,0	75	12	Violet
	GL4432	0 3,5	75	12	Violet
	GL4787	1 4,0	75	12	Violet
Straight Reverse Cutting 	GL2909	4/0 1,5	45	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GL4454	0 3,5	75	12	Violet
	GL4459	0 3,5	90	12	Violet
	GL4812	1 4,0	75	12	Violet
	GL4818	1 4,0	90	12	Violet
	GL5066	2 5,0	75	12	Violet
	GL5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GL5078	2 5,0	75	12	Violet
60 mm					
Straight Reverse Cutting 	GL2923	4/0 1,5	75	12	Violet
	GL3492	3/0 2,0	75	12	Violet
	GL4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GL4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GL4861	1 4,0	90	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOLAK Non-Needle Suture	GL2950	3/0 2,0	6x45	12	Violet
	GL4137	0 3,5	6x45	12	Violet
	GL4563	1 4,0	6x45	12	Violet

12x45 cm					
GLIKOLAK Non-Needle Suture	GL2344	4/0 1,5	12x45	12	Violet
	GL2958	3/0 2,0	12x45	12	Violet
	GL3552	2/0 3,0	12x45	12	Violet
	GL4145	0 3,5	12x45	12	Violet
	GL4554	1 4,0	12x45	12	Violet

17x45 cm					
GLIKOLAK Non-Needle Suture	GL2966	3/0 2,0	17x45	12	Violet
	GL4559	1 4,0	17x45	12	Violet

10x50 cm					
GLIKOLAK Non-Needle Suture	GL4140	0 3,5	10x50	12	Violet
	GL2955	3/0 2,0	10x75	12	Violet

10x75 cm					
GLIKOLAK Non-Needle Suture	GL4550	1 4,0	10x75	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
GLIKOLAK Non-Needle Suture	GL1737	5/0 1,0	150	12	Violet
	GL2326	4/0 1,5	150	12	Violet
	GL2935	3/0 2,0	150	12	Violet
	GL3511	2/0 3,0	150	12	Violet
	GL4114	0 3,5	150	12	Violet
	GL4532	1 4,0	150	12	Violet

Fastsorb

Surgical Suture, Polyglycolicacid (PGA)
Fast Absorbable, Synthetic, Braided, Undyed

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

Material

Polyglycolic Acid (PGA)

Coating

Calcium Stearate and polycaprolactone

Structure

Braided

Colour

Beige (undyed)

Size

EP : 0.4 to 8

USP : 8/0 to 6

Tensile Strenght Retention

7 Days 75 %

Mass Absorption

42 days

Sterilization Method

Gamma Irradiation

Main Indication

Gynaecology (episiotomies)

Urology

Ophtalmic surgery (conjunctive suturing)

Odontology (oral mucosa)












Pediatric surgery

Skin closer - intracutaneous, subcutaneous.

Ligatures








Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
3/8 Reverse Cutting 330 Micron			FS1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 330 Micron			FS1563	6/0 0,7	75	12	Undyed
1/2 Round Body, Taper Point 380 Micron			FS1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 330 Micron			FS5912	6/0 0,7	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FS1929	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting Premium 380 Micron			FS1941	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FS1989	5/0 1,0	75	12	Undyed
			FS2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FS1673	6/0 0,7	75	12	Undyed
			FS2019	5/0 1,0	75	12	Undyed
			FS2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			FS2557	4/0 1,5	75	12	Undyed
			FS3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FS2639	4/0 1,5	75	12	Undyed
			FS3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			FS2228	5/0 1,0	75	12	Undyed

Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			FS2683	4/0 1,5	75	12	Undyed
			FS3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS2250	5/0 1,0	75	12	Undyed
22 mm							
1/2 Round Body, Taper Point			FS2745	4/0 1,5	75	12	Undyed
25 mm							
1/2 Round Body, Taper Point			FS2787	4/0 1,5	75	12	Undyed
			FS3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS3273	3/0 2,0	75	12	Undyed
			FS3756	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body, Taper Point			FS3785	2/0 3,0	75	12	Undyed
			FS4214	0 3,5	75	12	Undyed
3/8 Reverse Cutting			FS3821	2/0 3,0	75	12	Undyed
			FS4236	0 3,5	75	12	Undyed
			FS3345	3/0 2,0	75	12	Undyed
30 mm							
1/2 Round Body, Taper Point			FS3368	3/0 2,0	75	12	Undyed
			FS3857	2/0 3,0	75	12	Undyed
			FS4243	0 3,5	75	12	Undyed
			FS4607	1 4,0	75	12	Undyed






Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting		FS2885	4/0 1,5	75	12	Undyed
		FS3396	3/0 2,0	75	12	Undyed
		FS4628	1 4,0	75	12	Undyed
35 mm						
1/2 Round Body, Taper Point		FS3915	2/0 3,0	75	12	Undyed
		FS4299	0 3,5	75	12	Undyed
		FS4305	0 3,5	90	12	Undyed
		FS4643	1 4,0	75	12	Undyed
		FS4648	1 4,0	90	12	Undyed
3/8 Reverse Cutting		FS3422	3/0 2,0	75	12	Undyed
		FS3932	2/0 3,0	75	12	Undyed
37 mm						
1/2 Round Body, Taper Point		FS3954	2/0 3,0	90	12	Undyed
40 mm						
1/2 Round Body, Taper Point		FS3981	2/0 3,0	75	12	Undyed
		FS3988	2/0 3,0	90	12	Undyed
		FS4367	0 3,5	75	12	Undyed
		FS4374	0 3,5	90	12	Undyed

Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	FS4712	1 4,0	75	12	Undyed
	FS4718	1 4,0	90	12	Undyed
	FS5001	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	FS4005	2/0 3,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	FS4421	0 3,5	75	12	Undyed
	FS4771	1 4,0	75	12	Undyed
	FS5034	2 5,0	75	12	Undyed
Straight Reverse Cutting 	FS2911	4/0 1,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	FS4037	2/0 3,0	75	12	Undyed
	FS4453	0 3,5	75	12	Undyed
	FS4811	1 4,0	75	12	Undyed
	FS5065	2 5,0	75	12	Undyed
65 mm					
1/2 Round Body, Taper Point 	FS4864	1 4,0	100	12	Undyed

Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
FASTSORB Non-Needle Suture	FS1373	6/0 0,7	150	12	Undyed
	FS2325	4/0 1,5	150	12	Undyed
	FS2934	3/0 2,0	150	12	Undyed
	FS3510	2/0 3,0	150	12	Undyed
	FS4113	0 3,5	150	12	Undyed
	FS4531	1 4,0	150	12	Undyed
	FS4907	2 5,0	150	12	Undyed

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)

Fast Absorbable, Synthetic, Braided, Undyed

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Coating

Poly(glycolide-co-lactide) (30/70) & calcium stearate

Structure

Braided

Colour

Beige (undyed)

Size

EP : 0.4 to 5

USP : 8/0 to 2

Tensile Strength Retention

7 Days 50 % - 10-14 Days 0 %

Mass Absorption

42 days

Sterilization Method

Gamma Irradiation






Main Indication

Gynaecology (episiotomies)
Urology
Ophtalmic surgery (conjunctive suturing)
Odontology (oral mucosa)
Pediatric surgery
Skin closer - intracutaneous, subcutaneous.
Ligatures

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed



Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
11 mm							
3/8 Reverse Cutting 330 Micron			FL1545	6/0 0,7	75	12	Undyed
12 mm							
3/8 Reverse Cutting 330 Micron			FL1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 380 Micron			FL1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FL1926	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FL1989	5/0 1,0	75	12	Undyed
			FL2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FL1673	6/0 0,7	75	12	Undyed
			FL2019	5/0 1,0	75	12	Undyed
			FL2539	4/0 1,5	75	12	Undyed
			FL3024	3/0 2,0	75	12	Undyed
3/8 Inside Cutting			FL2557	4/0 1,5	75	12	Undyed
			FL3038	3/0 2,0	75	12	Undyed
18 mm							
1/2 Round Body, Taper Point			FL3080	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FL2639	4/0 1,5	75	12	Undyed
			FL3108	3/0 2,0	75	12	Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm						
1/2 Round Body, Taper Point  		FL2228	5/0 1,0	75	12	Undyed
		FL2683	4/0 1,5	75	12	Undyed
		FL3141	3/0 2,0	75	12	Undyed
		FL3150	3/0 2,0	90	12	Undyed
3/8 Reverse Cutting  		FL2250	5/0 1,0	75	12	Undyed
22 mm						
1/2 Round Body, Taper Point  		FL2745	4/0 1,5	75	12	Undyed
		FL3195	3/0 2,0	75	12	Undyed
1/2 Inside Cutting  		FL3216	3/0 2,0	45	12	Undyed
3/8 Reverse Cutting  		FL3210	3/0 2,0	75	12	Undyed
25 mm						
1/2 Round Body, Taper Point  		FL2787	4/0 1,5	75	12	Undyed
		FL3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting  		FL2805	4/0 1,5	75	12	Undyed
26 mm						
1/2 Round Body, Taper Point  		FL3312	3/0 2,0	75	12	Undyed
		FL3785	2/0 3,0	75	12	Undyed
3/8 Reverse Cutting  		FL3345	3/0 2,0	75	12	Undyed
		FL3821	2/0 3,0	75	12	Undyed

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
30 mm							
1/2 Round Body, Taper Point			FL3368	3/0 2,0	75	12	Undyed
			FL3857	2/0 3,0	75	12	Undyed
			FL4243	0 3,5	75	12	Undyed
			FL4607	1 4,0	75	12	Undyed
3/8 Reverse Cutting			FL2885	4/0 1,5	75	12	Undyed
			FL3396	3/0 2,0	75	12	Undyed
			FL3890	2/0 3,0	75	12	Undyed
35 mm							
1/2 Round Body, Taper Point			FL3915	2/0 3,0	75	12	Undyed
			FL4299	0 3,5	75	12	Undyed
			FL4305	0 3,5	90	12	Undyed
			FL4643	1 4,0	75	12	Undyed
3/8 Reverse Cutting			FL3422	3/0 2,0	75	12	Undyed
			FL3932	2/0 3,0	75	12	Undyed
36 mm							
1/2 Round Body, Taper Point			FL4667	1 4,0	75	12	Undyed
40 mm							
1/2 Round Body, Taper Point			FL3981	2/0 3,0	75	12	Undyed
			FL3988	2/0 3,0	90	12	Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	FL4367	0 3,5	75	12	Undyed
	FL4374	0 3,5	90	12	Undyed
	FL4712	1 4,0	75	12	Undyed
	FL4718	1 4,0	90	12	Undyed
	FL5001	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	FL4005	2/0 3,0	75	12	Undyed
	FL4755	1 4,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	FL4424	0 3,5	90	12	Undyed
	FL4771	1 4,0	75	12	Undyed
	FL5034	2 5,0	75	12	Undyed
Straight Reverse Cutting 	FL2911	4/0 1,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	FL4037	2/0 3,0	75	12	Undyed
	FL4453	0 3,5	75	12	Undyed
	FL4811	1 4,0	75	12	Undyed
	FL5065	2 5,0	75	12	Undyed

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
FASTLAK Non-Needle Suture	FL1736	5/0 1,0	150	12	Undyed
	FL2325	4/0 1,5	150	12	Undyed
	FL2934	3/0 2,0	150	12	Undyed
	FL3510	2/0 3,0	150	12	Undyed
	FL4113	0 3,5	150	12	Undyed
	FL4531	1 4,0	150	12	Undyed
	FL4907	2 5,0	150	12	Undyed

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Surgical Sutures
Absorbable



Boz®

Definition

Short-term synthetic absorbable, closure of superficial wounds

Material

Poly(glycolide (75%)- co-caprolactone (25%)) (PGCL)

Coating

No coating

Structure

Monofilament

Colour

Undyed Beige (undyed) or Violet (dyed)

Size

EP : 0.7 to 5
USP : 6/0 to 2

Tensile Strength Retention

1. Week 45 % - 70 %
2. Week 20 %

Mass Absorption

90 - 120 Days

Sterilization Method





Ethylene Oxide Gas

Main Indication

Gynaecology
Obstetrics
Urology
Plastic Surgery
Odontology (oral mucosa)
Skin closer - intracutaneous, subcutaneous.
Reconstructive Surgery
Gastrointestinal Surgery


Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 330 Micron			MK1398	6/0 0,7	45	12	Undyed
			MK1759	5/0 1,0	45	12	Undyed
10 mm							
3/8 Reverse Cutting 330 Micron			MK1494	6/0 0,7	45	12	Undyed
13 mm							
3/8 Round Body, Taper Point 330 Micron			MK1578	6/0 0,7	45	12	Undyed
3/8 Round Body, Taper Point 330 Micron			MK1905	5/0 1,0	50	12	Undyed
16 mm							
1/2 Round Body, Taper Point			MK1989	5/0 1,0	75	12	Undyed
			MK2503	4/0 1,5	75	12	Undyed
			MK3004	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			MK2139	5/0 1,0	45	12	Undyed
			MK2557	4/0 1,5	75	12	Undyed
			MK3038	3/0 2,0	75	12	Undyed
18 mm							
3/8 Reverse Cutting			MK1704	6/0 0,7	45	12	Undyed
19 mm							
3/8 Reverse Cutting			MK2639	4/0 1,5	75	12	Undyed
			MK3108	3/0 2,0	75	12	Undyed

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body,Taper Point			MK2683	4/0 1,5	75	12	Undyed
			MK3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2718	4/0 1,5	75	12	Undyed
			MK3167	3/0 2,0	75	12	Undyed
24 mm							
3/8 Inside Cutting			MK3236	3/0 2,0	75	12	Undyed
25 mm							
1/2 Round Body,Taper Point			MK2787	4/0 1,5	75	12	Undyed
			MK3251	3/0 2,0	75	12	Undyed
			MK3252	3/0 2,0	75	12	Violet
			MK3723	2/0 3,0	75	12	Undyed
3/8 Reverse Cutting			MK2805	4/0 1,5	75	12	Undyed
			MK3273	3/0 2,0	75	12	Undyed
3/8 Inside Cutting			MK3762	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body,Taper Point			MK3312	3/0 2,0	75	12	Undyed
			MK3313	3/0 2,0	75	12	Violet
			MK3785	2/0 3,0	75	12	Undyed
1/2 Round Body,Taper Point			MK5846	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting			MK3335	3/0 2,0	75	12	Undyed

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament


Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
26 mm							
1/2 Reverse Cutting			MK4229	0 3,5	75	12	Undyed
			MK2851	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			MK3345	3/0 2,0	75	12	Undyed
			MK3821	2/0 3,0	75	12	Undyed
30 mm							
1/2 Round Body,Taper Point			MK3368	3/0 2,0	75	12	Undyed
			MK3369	3/0 2,0	75	12	Violet
			MK3857	2/0 3,0	75	12	Undyed
			MK3858	2/0 3,0	75	12	Violet
			MK4243	0 3,5	75	12	Undyed
			MK4607	1 4,0	75	12	Undyed
3/8 Round Body,Taper Point			MK4263	0 3,5	75	12	Undyed
3/8 Reverse Cutting			MK3890	2/0 3,0	75	12	Undyed
			MK4272	0 3,5	75	12	Undyed
			MK4628	1 4,0	75	12	Undyed
35 mm							
1/2 Round Body,Taper Point			MK3915	2/0 3,0	75	12	Undyed
			MK4299	0 3,5	75	12	Undyed
			MK4643	1 4,0	75	12	Undyed


Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
35 mm						
1/2 Reverse Cutting		MK4658	1 4,0	75	12	Undyed
3/8 Reverse Cutting		MK3932	2/0 3,0	75	12	Undyed




40 mm						
1/2 Round Body, Taper Point 	MK3981	2/0 3,0	75	12	Undyed	
	MK4367	0 3,5	75	12	Undyed	
	MK4374	0 3,5	90	12	Undyed	
	MK4712	1 4,0	75	12	Undyed	
	MK4713	1 4,0	75	12	Violet	
	MK4718	1 4,0	90	12	Undyed	
	MK5001	2 5,0	75	12	Undyed	
1/2 Reverse Cutting 	MK4744	1 4,0	75	12	Undyed	
	MK5018	2 5,0	75	12	Undyed	
3/8 Reverse Cutting 	MK4407	0 3,5	75	12	Undyed	

45 mm						
1/2 Round Body, Taper Point 	MK4771	1 4,0	75	12	Undyed	
	MK5034	2 5,0	75	12	Undyed	

50 mm					
1/2 Round Body,Taper Point 	MK4037	2/0 3,0	75	12	Undyed
	MK4811	1 4,0	75	12	Undyed

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament


Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	MK4817	1 4,0	90	12	Undyed
	MK5065	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	MK4474	0 3,5	75	12	Undyed
	MK4840	1 4,0	75	12	Undyed
63 mm					
3/8 Blunt 	MK4859	1 4,0	100	12	Undyed
	MK4857	1 4,0	75	12	Undyed
Non-Needle Suture	Product No	USP EP	Length (cm)	Unit	Color
150 cm					
MONOKAPROL Non-Needle Suture	MK2934	3/0 2,0	150	12	Undyed
	MK3510	2/0 3,0	150	12	Undyed
	MK4113	0 3,5	150	12	Undyed
	MK4531	1 4,0	150	12	Undyed
	MK4907	2 5,0	150	12	Undyed

Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable,
Synthetic, Monofilament, Violet

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, wound support through an extended healing period

Material

Polydioxanone (PDO)

Coating

No coating

Structure

Monofilament

Colour

Violet (dyed)

Size

EP : 0.5 to 5

USP : 7/0 to 2

Tensile Strength Retention

2. Week 70 %

Mass Absorption

180-240 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Orthopaedics

General Surgery









Abdominal Wall Closure

Paediatric cardiovascular tissues

Ophthalmic Surgery












Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
3/8 Round Body, Taper Point 220 Micron			PX1289	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX5444	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX1426	6/0 0,7	45	12	Violet
10 mm							
3/8 Reverse Cutting 330 Micron			PX5850	6/0 0,7	45	12	Violet
11 mm							
3/8 Round Body, Taper Point 220 Micron			PX1509	6/0 0,7	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 330 Micron			PX1564	6/0 0,7	75	12	Violet
1/2 Round Body, Taper Point 380 Micron			PX1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			PX1897	5/0 1,0	75	12	Violet
1/2 Spatula 280 Micron			PX1599	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX1362	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1575	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1576	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1128	6/0 0,7	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			PX1914	5/0 0,7	75	12	Violet
3/8 Tapercutting 280 Micron			PX1609	6/0 0,7	45	12	Violet









Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			PX1990	5/0 1,0	75	12	Violet
			PX2504	4/0 1,5	75	12	Violet
			PX3005	3/0 2,0	75	12	Violet
1/2 Round Body, Taper Point			PX1999	5/0 1,0	90	12	Violet
			PX2513	4/0 1,5	90	12	Violet
3/8 Round Body, Taper Point			PX2005	5/0 1,0	75	12	Violet
			PX2520	4/0 1,5	75	12	Violet
3/8 Round Body, Taper Point			PX2523	4/0 1,5	90	12	Violet
3/8 Reverse Cutting			PX2021	5/0 1,0	75	12	Violet
			PX2540	4/0 1,5	75	12	Violet
17 mm							
1/2 Round Body, Taper Point			PX5517	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point			PX5516	4/0 1,5	75	12	Violet
19 mm							
3/8 Reverse Cutting			PX2635	4/0 1,5	45	12	Violet
			PX3109	3/0 2,0	75	12	Violet
20 mm							
1/2 Round Body, Taper Point			PX2229	5/0 1,0	75	12	Violet
			PX2685	4/0 1,5	75	12	Violet



Polidiox



Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet


Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			PX3142	3/0 2,0	75	12	Violet
1/2 Inside Cutting			PX3175	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			PX2719	4/0 1,5	75	12	Violet
22 mm							
1/2 Inside Cutting			PX3217	3/0 2,0	45	12	Violet
25 mm							
1/2 Round Body, Taper Point			PX3252	3/0 2,0	75	12	Violet
			PX3724	2/0 3,0	75	12	Violet
			PX4190	0 3,5	75	12	Violet
26 mm							
1/2 Round Body, Taper Point			PX2828	4/0 1,5	75	12	Violet
			PX3313	3/0 2,0	75	12	Violet
			PX3786	2/0 3,0	75	12	Violet
1/2 Reverse Cutting			PX3336	3/0 2,0	75	12	Violet
			PX3813	2/0 3,0	75	12	Violet
3/8 Reverse Cutting			PX3341	3/0 2,0	45	12	Violet
			PX3346	3/0 2,0	75	12	Violet
3/8 Inside Cutting			PX3360	3/0 2,0	75	12	Violet


Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point  	PX3369	3/0 2,0	75	12	Violet
	PX3858	2/0 3,0	75	12	Violet
	PX4244	0 3,5	75	12	Violet
3/8 Reverse Cutting  	PX3891	2/0 3,0	75	12	Violet



40 mm					
1/2 Round Body, Taper Point  	PX3982	2/0 3,0	75	12	Violet
	PX4368	0 3,5	75	12	Violet
	PX4713	1 4,0	75	12	Violet
	PX4719	1 4,0	90	12	Violet
	PX5006	2 5,0	90	12	Violet
1/2 Round Body, Taper Point LOOP  	PX4388	0 3,5	150	12	Violet
	PX4736	1 4,0	150	12	Violet

45 mm					
1/2 Round Body, Taper Point Heavy  	PX4779	1 4,0	90	12	Violet

50 mm					
1/2 Round Body Taper Point  	PX4459	0 3,5	90	12	Violet
	PX4818	1 4,0	90	12	Violet
	PX5070	2 5,0	90	12	Violet

Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point LOOP 	PX4827	1 4,0	150	12	Violet
1/2 Round Body, Taper Point Heavy, LOOP 	PX4466	0 3,5	150	12	Violet



Boz Medical Device Industry and Trade Inc.

Antibacterial **Absorbable Sutures**

GLIKOSORB Plus

GLIKOLAK Plus

MONOKAPROL Plus

POLIDIOX Plus

web : www.boztibbi.com
email: boz@boztibbi.com

Surgical Sutures Absorbable

Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Boz®



Definition

Synthetic absorbable, coated, antibacterial mid term wound support and mid term mass absorption.

Material

Polyglycolic Acid (PGA)

Antibacterial Agent

CHA (Chlorhexidine diacetate)

Coating

Calcium Stearate and poly(glycolide-co-lactide)

Structure

Braided - Multifilament

Colour

Violet (dyed)

Size

EP : 1 to 5

USP : 5/0 to 2

Tensile Strenght Retention

2. Week 90 % - 3. Week 50 %

Mass Absorption

70-90 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery

Gynaecology

Orthopedics


Urology

Skin closer - intracutaneous, subcutaneous.

Ligatures

Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 430 Micron			GSP2363	4/0 1,5	45	12	Violet
3/8 Spatula 280 Micron			GSP1768	5/0 1,0	45	12	Violet
3/8 Spatula 340 Micron			GSP1764	5/0 1,0	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 430 Micron			GSP1897	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			GSP1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GSP1894	5/0 1,0	75	12	Violet
			GSP2434	4/0 1,5	60	12	Violet
			GSP2436	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 560 Micron			GSP2456	4/0 1,5	75	12	Violet
3/8 Spatula 330 Micron			GSP1873	5/0 1,0	60	12	Violet
16 mm							
1/2 Round Body, Taper Point			GSP1992	5/0 1,0	75	12	Undyed
			GSP2503	4/0 1,5	75	12	Undyed
			GSP3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GSP2021	5/0 1,0	75	12	Violet
			GSP2540	4/0 1,5	75	12	Violet
			GSP3025	3/0 2,0	75	12	Violet






Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Inside Cutting			GSP2140	5/0 1,0	45	12	Violet
			GSP2557	4/0 1,5	75	12	Undyed
			GSP3039	3/0 2,0	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GSP2613	4/0 1,5	75	12	Violet
3/8 Round Body, Taper Point			GSP3085	3/0 2,0	75	12	Violet
19 mm							
1/2 Round Body, Taper Point			GSP2197	5/0 1,0	75	12	Violet
3/8 Reverse Cutting			GSP2206	5/0 1,0	75	12	Violet
			GSP2640	4/0 1,5	75	12	Violet
3/8 Reverse Cutting Premium			GSP2211	5/0 1,0	45	12	Violet
20 mm							
1/2 Round Body, Taper Point			GSP2229	5/0 1,0	75	12	Violet
			GSP2685	4/0 1,5	75	12	Violet
			GSP3142	3/0 2,0	75	12	Violet
1/2 Tapercutting			GSP2673	4/0 1,5	75	12	Violet
3/8 Reverse Cutting			GSP2250	5/0 1,0	75	12	Undyed
			GSP2716	4/0 1,5	45	12	Undyed
			GSP2719	4/0 1,5	75	12	Violet

Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm					
3/8 Reverse Cutting  	GSP3168	3/0 2,0	75	12	Violet
	GSP3663	2/0 3,0	75	12	Violet
22 mm					
1/2 Round Body, Taper Point  	GSP2747	4/0 1,5	75	12	Violet
3/8 Reverse Cutting  	GSP3214	3/0 2,0	90	12	Violet
25 mm					
1/2 Round Body, Taper Point  	GSP2788	4/0 1,5	75	12	Violet
	GSP3252	3/0 2,0	75	12	Violet
	GSP3746	2/0 3,0	100	12	Violet
	GSP3724	2/0 3,0	75	12	Violet
1/2 Tapercutting  	GSP3718	2/0 3,0	75	12	Violet
3/8 Reverse Cutting  	GSP2804	4/0 1,5	75	12	Violet
	GSP3274	3/0 2,0	75	12	Violet
	GSP3757	2/0 3,0	75	12	Violet
3/8 Inside Cutting  	GSP2812	4/0 1,5	45	12	Violet
	GSP3284	3/0 2,0	45	12	Violet
26 mm					
1/2 Round Body, Taper Point  	GSP2292	5/0 1,0	75	12	Violet
	GSP2828	4/0 1,5	75	12	Violet




Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point 	GSP3313	3/0 2,0	75	12	Violet
	GSP3318	3/0 2,0	90	12	Violet
	GSP3786	2/0 3,0	75	12	Violet
3/8 Round Body, Taper Point 	GSP2841	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 	GSP3345	3/0 2,0	75	12	Undyed
	GSP3346	3/0 2,0	75	12	Violet
	GSP3821	2/0 3,0	75	12	Undyed
	GSP3822	2/0 3,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point 	GSP3858	2/0 3,0	75	12	Violet
	GSP4244	0 3,5	75	12	Violet
	GSP4608	1 4,0	75	12	Violet
1/2 Reverse Cutting 	GSP3882	2/0 3,0	75	12	Violet
	GSP4266	0 3,5	75	12	Violet
1/2 Tapercutting 	GSP3843	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GSP2887	4/0 1,5	75	12	Violet
	GSP3397	3/0 2,0	75	12	Violet
	GSP3891	2/0 3,0	75	12	Violet
	GSP4273	0 3,5	75	12	Violet




Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
1/2 Round Body, Taper Point 	GSP3411	3/0 2,0	75	12	Violet
	GSP3916	2/0 3,0	75	12	Violet
	GSP3922	2/0 3,0	90	12	Violet
	GSP4300	0 3,5	75	12	Violet
	GSP4644	1 4,0	75	12	Violet
3/8 Reverse Cutting 	GSP3423	3/0 2,0	75	12	Violet
	GSP3933	2/0 3,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	GSP3982	2/0 3,0	75	12	Violet
	GSP4368	0 3,5	75	12	Violet
	GSP4375	0 3,5	90	12	Violet
	GSP4713	1 4,0	75	12	Violet
	GSP4719	1 4,0	90	12	Violet
	GSP5003	2 5,0	75	12	Violet
	GSP5006	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GSP4397	0 3,5	75	12	Violet
	GSP4745	1 4,0	75	12	Violet
	GSP4749	1 4,0	90	12	Violet
	GSP5023	2 5,0	90	12	Violet





Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm						
1/2 Tapercutting		GSP4363	0 3,5	75	12	Violet
		GSP4408	0 3,5	75	12	Violet
3/8 Reverse Cutting		GSP4756	1 4,0	75	12	Violet
45 mm						
1/2 Round Body, Taper Point		GSP4422	0 3,5	75	12	Violet
		GSP4772	1 4,0	75	12	Violet
		GSP5035	2 5,0	75	12	Violet
		GSP5040	2 5,0	90	12	Violet
1/2 Tapercutting		GSP4766	1 4,0	75	12	Violet
3/8 Reverse Cutting		GSP4432	0 3,5	75	12	Violet
		GSP4787	1 4,0	75	12	Violet
Straight Reverse Cutting		GSP2909	4/0 1,5	45	12	Violet
50 mm						
1/2 Round Body, Taper Point		GSP4454	0 3,5	75	12	Violet
		GSP4459	0 3,5	90	12	Violet
		GSP4812	1 4,0	75	12	Violet
		GSP4818	1 4,0	90	12	Violet

Antibacterial Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GSP5066	2 5,0	75	12	Violet
	GSP5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GSP4471	0 3,5	75	12	Violet
60 mm					
Straight Reverse Cutting 	GSP3492	3/0 2,0	75	12	Violet
	GSP4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GSP4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GSP4500	0 3,5	90	12	Violet
	GSP4861	1 4,0	90	12	Violet

Antibacterial

Glikosorb Plus

Surgical Suture, Polyglycolic acid (PGA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOSORB Plus Non-Needle Suture	GSP2950	3/0 2,0	6x45	12	Violet
	GSP3543	2/0 3,0	6x45	12	Violet
	GSP4137	0 3,5	6x45	12	Violet
	GSP4563	1 4,0	6x45	12	Violet
12x45 cm					
GLIKOSORB Plus Non-Needle Suture	GSP2344	4/0 1,5	12x45	12	Violet
	GSP2958	3/0 2,0	12x45	12	Violet
	GSP3552	2/0 3,0	12x45	12	Violet
	GSP4145	0 3,5	12x45	12	Violet
7x75 cm					
GLIKOSORB Plus Non-Needle Suture	GSP2952	3/0 2,0	7x75	12	Violet
10x75 cm					
GLIKOSORB Plus Non-Needle Suture	GSP4550	1 4,0	10x75	12	Violet
150 cm					
GLIKOSORB Plus Non-Needle Suture	GSP1737	5/0 1,0	150	12	Violet
	GSP2326	4/0 1,5	150	12	Violet
	GSP2935	3/0 2,0	150	12	Violet
	GSP3511	2/0 3,0	150	12	Violet
	GSP4114	0 3,5	150	12	Violet
	GSP4532	1 4,0	150	12	Violet

Antibacterial Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, coated, antibacterial mid term wound support and mid term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

Poly(glycolide-co-lactide) & calcium stearate

Structure

Braided-Multifilament

Colour

Beige (undyed) or Violet (dyed)

Size

EP : 1 to 5
USP : 5/0 to 2

Tensile Strenght Retention

2. Week 75%

Mass Absorption

55-70 days

Sterilization Method

Ethylene Oxide Gas

Main Indication





Gastrointestinal Surgery
Gynaecology
Orthopedics
Urology
Skin closer - intracutaneous, subcutaneous.
Ligatures


Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
8 mm							
1/4 Spatula 340 Micron			GLP1764	5/0 1,0	45	12	Violet
3/8 Spatula 230 Micron			GLP1768	5/0 1,0	45	12	Violet
3/8 Round Body, Taper Point 380 Micron			GLP1778	5/0 1,0	75	12	Violet

13 mm							
1/2 Round Body, Taper Point 380 Micron			GLP1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			GLP1897	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			GLP1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GLP2434	4/0 1,5	60	12	Violet
3/8 Round Body, Taper Point 430 Micron			GLP2436	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 380 Micron			GLP1928	5/0 1,0	75	12	Violet

16 mm							
1/2 Round Body,Taper Point			GLP1990	5/0 1,0	75	12	Violet
			GLP2504	4/0 1,5	75	12	Violet
			GLP3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GLP2021	5/0 1,0	75	12	Violet
			GLP2540	4/0 1,5	75	12	Violet
			GLP3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GLP2144	5/0 1,0	75	12	Violet
			GLP2553	4/0 1,5	45	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Inside Cutting			GLP2557	4/0 1,5	75	12	Undyed
			GLP3039	3/0 2,0	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GLP2613	4/0 1,5	75	12	Violet
19 mm							
3/8 Reverse Cutting			GLP2206	5/0 1,0	75	12	Violet
			GLP2640	4/0 1,5	75	12	Violet
20 mm							
1/2 Round Body, Taper Point			GLP2229	5/0 1,0	75	12	Violet
			GLP2685	4/0 1,5	75	12	Violet
			GLP2690	4/0 1,5	90	12	Violet
			GLP3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GLP2716	4/0 1,5	45	12	Undyed
			GLP3164	3/0 2,0	45	12	Violet
			GLP3168	3/0 2,0	75	12	Violet
22 mm							
3/8 Reverse Cutting			GLP2763	4/0 1,5	75	12	Violet
25 mm							
1/2 Round Body, Taper Point			GLP2788	4/0 1,5	75	12	Violet
			GLP3252	3/0 2,0	75	12	Violet

Antibacterial

Glikolak Plus


Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	GLP3746	2/0 3,0	100	12	Violet
	GLP3724	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GLP2804	4/0 1,5	75	12	Violet
3/8 Inside Cutting 	GLP2812	4/0 1,5	45	12	Violet
	GLP3284	3/0 2,0	45	12	Violet
	GLP3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point 	GLP2828	4/0 1,5	75	12	Violet
	GLP3313	3/0 2,0	75	12	Violet
	GLP3786	2/0 3,0	75	12	Violet
	GLP4220	0 3,5	90	12	Violet
3/8 Reverse Cutting 	GLP2852	4/0 1,5	75	12	Violet
	GLP3346	3/0 2,0	75	12	Violet
	GLP3822	2/0 3,0	75	12	Violet

30 mm					
1/2 Round Body, Taper Point 	GLP3369	3/0 2,0	75	12	Violet
	GLP3858	2/0 3,0	75	12	Violet
	GLP3865	2/0 3,0	90	12	Violet
	GLP4244	0 3,5	75	12	Violet





Antibacterial**Glikolak Plus**Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	GLP4250	0 3,5	90	12	Violet
	GLP4608	1 4,0	75	12	Violet
	GLP4614	1 4,0	90	12	Violet
1/2 Round Body, Taper Point 	GLP3854	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GLP3397	3/0 2,0	75	12	Violet
	GLP3891	2/0 3,0	75	12	Violet
	GLP3895	2/0 3,0	90	12	Violet
	GLP4273	0 3,5	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GLP3411	3/0 2,0	75	12	Violet
	GLP3916	2/0 3,0	75	12	Violet
	GLP4300	0 3,5	75	12	Violet
	GLP4644	1 4,0	75	12	Violet
	GLP4649	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GLP3423	3/0 2,0	75	12	Violet
	GLP3933	2/0 3,0	75	12	Violet
3/8 Inside Cutting 	GLP3938	2/0 3,0	75	12	Violet

Antibacterial

Glikolak Plus




Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	GLP3982	2/0 3,0	75	12	Violet
	GLP4368	0 3,5	75	12	Violet
	GLP4375	0 3,5	90	12	Violet
	GLP4713	1 4,0	75	12	Violet
	GLP4719	1 4,0	90	12	Violet
	GLP5003	2 5,0	75	12	Violet
1/2 Reverse Cutting 	GLP4399	0 3,5	90	12	Violet
	GLP4749	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GLP4756	1 4,0	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GLP4422	0 3,5	75	12	Violet
	GLP4772	1 4,0	75	12	Violet
	GLP4775	1 4,0	90	12	Violet
	GLP5035	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GLP4787	1 4,0	75	12	Violet
50 mm					
1/2 Round Body, Taper Point 	GLP4454	0 3,5	75	12	Violet
	GLP4459	0 3,5	90	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GLP4812	1 4,0	75	12	Violet
	GLP4818	1 4,0	90	12	Violet
	GLP5066	2 5,0	75	12	Violet
	GLP5070	2 5,0	90	12	Violet
60 mm					
Straight Reverse Cutting 	GLP2923	4/0 1,5	75	12	Violet
	GLP3492	3/0 2,0	75	12	Violet
	GLP3498	3/0 2,0	90	12	Violet
	GLP4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GLP4085	2/0 3,0	75	12	Violet
	GLP4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GLP4861	1 4,0	90	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2950	3/0 2,0	6x45	12	Violet
	GLP4137	0 3,5	6x45	12	Violet

12x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2344	4/0 1,5	12x45	12	Violet
	GLP2958	3/0 2,0	12x45	12	Violet
	GLP3552	2/0 3,0	12x45	12	Violet
	GLP4145	0 3,5	12x45	12	Violet

17x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2349	4/0 1,5	17x45	12	Violet
	GLP2966	3/0 2,0	17x45	12	Violet
	GLP3558	2/0 3,0	17x45	12	Violet
	GLP4149	0 3,5	17x45	12	Violet

10x50 cm					
GLIKOLAK Plus Non-Needle Suture	GLP3547	2/0 3,0	10x50	12	Violet
	GLP4140	0 3,5	10x50	12	Violet
	GLP4549	1 4,0	10x50	12	Violet

7x75 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2952	3/0 2,0	7x75	12	Violet
	GLP3545	2/0 3,0	7x75	12	Violet
	GLP4139	0 3,5	7x75	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
10x75 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2341	4/0 1,5	10x75	12	Violet
	GLP4141	0 3,5	10x75	12	Violet

150 cm					
GLIKOLAK Plus Non-Needle Suture	GLP1737	5/0 1,0	150	12	Violet
	GLP2326	4/0 1,5	150	12	Violet
	GLP2934	3/0 2,0	150	12	Undyed
	GLP2935	3/0 2,0	150	12	Violet
	GLP3511	2/0 3,0	150	12	Violet
	GLP4114	0 3,5	150	12	Violet

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, antibacterial, closure of superficial wounds

Material

Poly(glycolide (75%)- co-caprolactone (25%) (PGCL)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

No coating

Structure

Monofilament

Colour

Undyed Beige (undyed) or Violet (dyed)

Size

EP : 1 to 5

USP : 5/0 to 2

Tensile Strenght Retention

1. Week 70 % - 75 % - 2. Week 40 % - 50%

Mass Absorption

90 - 120 Days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gynaecology

Obstetrics

Urology


Plastic Surgery

Odontology (oral mucosa)

Skin closer - intracutaneous, subcutaneous.





Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 330 Micron			MKP1759	5/0 1,0	45	12	Undyed
13 mm							
3/8 Round Body, Taper Point 330 Micron			MKP1905	5/0 1,0	50	12	Undyed
16 mm							
1/2 Round Body, Taper Point			MKP1985	5/0 1,0	45	12	Undyed
			MKP1989	5/0 1,0	75	12	Undyed
			MKP2503	4/0 1,5	75	12	Undyed
			MKP3004	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MKP2019	5/0 1,0	75	12	Undyed
			MKP2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			MKP2139	5/0 1,0	45	12	Undyed
			MKP2557	4/0 1,5	75	12	Undyed
			MKP3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			MKP2639	4/0 1,5	75	12	Undyed
			MKP3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			MKP2683	4/0 1,5	75	12	Undyed
			MKP3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MKP2718	4/0 1,5	75	12	Undyed
			MKP3167	3/0 2,0	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament


Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point  	MKP2787	4/0 1,5	75	12	Undyed
	MKP3251	3/0 2,0	75	12	Undyed
	MKP3723	2/0 3,0	75	12	Undyed
3/8 Inside Cutting  	MKP3762	2/0 3,0	75	12	Undyed
26 mm					
1/2 Round Body, Taper Point  	MKP3312	3/0 2,0	75	12	Undyed
	MKP3785	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting  	MKP3335	3/0 2,0	75	12	Undyed
30 mm					
1/2 Round Body, Taper Point  	MKP3368	3/0 2,0	75	12	Undyed
	MKP3857	2/0 3,0	75	12	Undyed
	MKP4243	0 3,5	75	12	Undyed
	MKP4249	0 3,5	90	12	Undyed
	MKP4607	1 4,0	75	12	Undyed
1/2 Round Body, Taper Point  	MKP3853	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting  	MKP4265	0 3,5	75	12	Undyed
3/8 Round Body, Taper Point  	MKP4263	0 3,5	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
3/8 Reverse Cutting 	MKP3890	2/0 3,0	75	12	Undyed
	MKP4272	0 3,5	75	12	Undyed

35 mm					
1/2 Round Body, Taper Point 	MKP3915	2/0 3,0	75	12	Undyed
	MKP4299	0 3,5	75	12	Undyed
	MKP4643	1 4,0	75	12	Undyed
3/8 Reverse Cutting 	MKP3932	2/0 3,0	75	12	Undyed
	MKP4319	0 3,5	75	12	Undyed


40 mm					
1/2 Round Body, Taper Point 	MKP3981	2/0 3,0	75	12	Undyed
	MKP4367	0 3,5	75	12	Undyed
	MKP4374	0 3,5	90	12	Undyed
	MKP4723	1 4,0	100	12	Undyed
	MKP4712	1 4,0	75	12	Undyed
	MKP4718	1 4,0	90	12	Undyed
	MKP5001	2 5,0	75	12	Undyed
1/2 Reverse Cutting 	MKP4744	1 4,0	75	12	Undyed
	MKP5018	2 5,0	75	12	Undyed


Antibacterial


Monokaprol Plus


Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
3/8 Reverse Cutting 	MKP4407	0 3,5	75	12	Undyed
	MKP4755	1 4,0	75	12	Undyed

45 mm					
1/2 Round Body, Taper Point 	MKP4421	0 3,5	75	12	Undyed
	MKP4771	1 4,0	75	12	Undyed
3/8 Reverse Cutting 	MKP4431	0 3,5	75	12	Undyed

50 mm					
1/2 Round Body, Taper Point 	MKP4037	2/0 3,0	75	12	Undyed
	MKP4453	0 3,5	75	12	Undyed
	MKP4811	1 4,0	75	12	Undyed

50 mm					
1/2 Round Body, Taper Point 	MKP4817	1 4,0	90	12	Undyed
	MKP5065	2 5,0	75	12	Undyed

63 mm					
3/8 Blunt 	MKP4857	1 4,0	75	12	Undyed

Antibacterial

Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament


Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
MONOKAPROL Plus Non-Needle Suture	MKP2934	3/0 2,0	150	12	Undyed
	MKP3510	2/0 3,0	150	12	Undyed
	MKP4113	0 3,5	150	12	Undyed
	MKP4531	1 4,0	150	12	Undyed
	MKP4907	2 5,0	150	12	Undyed

Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO) Absorbable,
Synthetic, Monofilament, Violet

Surgical Sutures
Absorbable

Boz®



Definition

Short-term synthetic absorbable, antibacterial, wound support through an extended healing period.

Material

Polydioxanone (PDO)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

No coating

Structure

Monofilament

Colour

Violet (dyed)

Size

EP : 1 to 5

USP : 5/0 to 2

Tensile Strenght Retention

2. Week 70%

4. Week 50%

Mass Absorption

180 - 240 Days

Sterilization Method

Ethylene Oxide Gas

Main Indication










Orthopaedics

General Surgery

Abdominal Wall Closure


Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 430 Micron			PXP1897	5/0 1,0	75	12	Violet
16 mm							
1/2 Round Body, Taper Point			PXP1990	5/0 1,0	75	12	Violet
			PXP2504	4/0 1,5	75	12	Violet
			PXP3005	3/0 2,0	75	12	Violet
1/2 Round Body, Taper Point			PXP1999	5/0 1,0	90	12	Violet
			PXP2513	4/0 1,5	90	12	Violet
3/8 Round Body, Taper Point			PXP2005	5/0 1,0	75	12	Violet
			PXP2520	4/0 1,5	75	12	Violet
3/8 Round Body, Taper Point			PXP2523	4/0 1,5	90	12	Violet
3/8 Reverse Cutting			PXP2021	5/0 1,0	75	12	Violet
19 mm							
3/8 Reverse Cutting			PXP2635	4/0 1,5	45	12	Violet
20 mm							
1/2 Round Body, Taper Point			PXP2685	4/0 1,5	75	12	Violet
			PXP3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			PXP2719	4/0 1,5	75	12	Violet
22 mm							
1/2 Inside Cutting			PXP3217	3/0 2,0	45	12	Violet







Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	PXP3252	3/0 2,0	75	12	Violet
	PXP3724	2/0 3,0	75	12	Violet
	PXP4190	0 3,5	75	12	Violet
26 mm					
1/2 Round Body, Taper Point 	PXP2828	4/0 1,5	75	12	Violet
	PXP3313	3/0 2,0	75	12	Violet
	PXP3786	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	PXP3346	3/0 2,0	75	12	Violet
3/8 Inside Cutting 	PXP3360	3/0 2,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point 	PXP3369	3/0 2,0	75	12	Violet
	PXP3858	2/0 3,0	75	12	Violet
	PXP4244	0 3,5	75	12	Violet
3/8 Reverse Cutting 	PXP3891	2/0 3,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	PXP3982	2/0 3,0	75	12	Violet
	PXP4368	0 3,5	75	12	Violet
	PXP4375	0 3,5	90	12	Violet
	PXP4713	1 4,0	75	12	Violet

Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	PXP4719	1 4,0	90	12	Violet
	PXP5006	2 5,0	90	12	Violet
1/2 Round Body, Taper Point, LOOP 	PXP4388	0 3,5	150	12	Violet
	PXP4736	1 4,0	150	12	Violet
45 mm					
1/2 Round Body, Taper Point Heavy 	PXP4779	1 4,0	90	12	Violet
50 mm					
1/2 Round Body, Taper Point 	PXP4459	0 3,5	90	12	Violet
	PXP4818	1 4,0	90	12	Violet
	PXP5070	2 5,0	90	12	Violet
1/2 Round Body, Taper Point LOOP 	PXP4827	1 4,0	150	12	Violet
1/2 Round Body, Taper Point Heavy LOOP 	PXP4466	0 3,5	150	12	Violet



Boz Medical Device Industry and Trade Inc.

Absorbable Sutures

PLAIN CATGUT

CHROMIC CATGUT

Plain Catgut

Absorbable Surgical Suture



Boz®

Definition

Collagen based biomaterial, absorbable.

Material

Collagen tissue derived from purified sheep intestine

Coating

No coating

Structure

Monofilament

Colour

Yellowish

Size

EP : 1 to 8

USP : 6/0 to 4

Sterilization Method

Ethylene Oxide Gas

Gamma Irradiation


Main Indication

Ligatures

Ophthalmic Surgery

Plain Catgut

Absorbable Surgical Suture




Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 330 Micron			CP1563	6/0 1,0	75	12	Yellowish
1/2 Round Body, Taper Point 430 Micron			CP1892	5/0 1,5	75	12	Yellowish
1/2 Round Body, Taper Point 330 Micron With Alcohol			CP1563A	6/0 1,0	75	12	Yellowish
3/8 Reverse Cutting 430 Micron			CP1933	5/0 1,5	75	12	Yellowish
3/8 Reverse Cutting 430 Micron With Alcohol			CP1933A	5/0 1,5	75	12	Yellowish
16 mm							
1/2 Round Body, Taper Point			CP1989	5/0 1,5	75	12	Yellowish
			CP2503	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol			CP1989A	5/0 1,5	75	12	Yellowish
			CP2503A	4/0 2,0	75	12	Yellowish
			CP3004A	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting			CP2019	5/0 1,5	75	12	Yellowish
			CP2539	4/0 2,0	75	12	Yellowish
			CP3024	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol			CP2019A	5/0 1,5	75	12	Yellowish
			CP2539A	4/0 2,0	75	12	Yellowish
			CP3024A	3/0 3,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point	 		CP2227	5/0 1,5	50	12	Yellowish
			CP2228	5/0 1,5	75	12	Yellowish
			CP2683	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol	 		CP2683A	4/0 2,0	75	12	Yellowish
3/8 Reverse Cutting	 		CP2718	4/0 2,0	75	12	Yellowish
			CP3167	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol	 		CP2718A	4/0 2,0	75	12	Yellowish
25 mm							
1/2 Round Body, Taper Point With Alcohol	 		CP2787A	4/0 2,0	75	12	Yellowish
3/8 Reverse Cutting	 		CP2276	5/0 1,5	75	12	Yellowish
26 mm							
1/2 Round Body, Taper Point	 		CP2827	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol	 		CP2827A	4/0 2,0	75	12	Yellowish
			CP3312A	3/0 3,0	75	12	Yellowish
			CP3785A	2/0 3,5	75	12	Yellowish
3/8 Reverse Cutting	 		CP2851	4/0 2,0	75	12	Yellowish
			CP3345	3/0 3,0	75	12	Yellowish
			CP3821	2/0 3,5	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
3/8 Reverse Cutting With Alcohol 	CP3345A	3/0 3,0	75	12	Yellowish
	CP3821A	2/0 3,5	75	12	Yellowish
	CP4236A	0 4,0	75	12	Yellowish
30 mm					
1/2 Round Body, Taper Point 	CP3368	3/0 3,0	75	12	Yellowish
	CP4243	0 4,0	75	12	Yellowish
	CP4607	1 5,0	75	12	Yellowish
	CP4954	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP2874A	4/0 2,0	75	12	Yellowish
	CP3368A	3/0 3,0	75	12	Yellowish
	CP3857A	2/0 3,5	75	12	Yellowish
	CP4243A	0 4,0	75	12	Yellowish
	CP4607A	1 5,0	75	12	Yellowish
3/8 Reverse Cutting 	CP3396	3/0 3,0	75	12	Yellowish
	CP3890	2/0 3,5	75	12	Yellowish
	CP4272	0 4,0	75	12	Yellowish
	CP4628	1 5,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture








Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
3/8 Reverse Cutting With Alcohol 	CP3396A	3/0 3,0	75	12	Yellowish
	CP3890A	2/0 3,5	75	12	Yellowish
	CP4272A	0 4,0	75	12	Yellowish
	CP4628A	1 5,0	75	12	Yellowish
35 mm					
1/2 Round Body, Taper Point 	CP3410	3/0 3,0	75	12	Yellowish
	CP3915	2/0 3,5	75	12	Yellowish
	CP4299	0 4,0	75	12	Yellowish
	CP4969	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP3410A	3/0 3,0	75	12	Yellowish
	CP4299A	0 4,0	75	12	Yellowish
	CP4643A	1 5,0	75	12	Yellowish
	CP4969A	2 6,0	75	12	Yellowish
1/2 Reverse Cutting 	CP4976	2 6,0	75	12	Yellowish
3/8 Reverse Cutting 	CP3422	3/0 3,0	75	12	Yellowish
	CP3932	2/0 3,5	75	12	Yellowish
	CP4319	0 4,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP3422A	3/0 3,0	75	12	Yellowish
	CP4319A	0 4,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture



Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	CP4367	0 4,0	75	12	Yellowish
	CP4374	0 4,0	90	12	Yellowish
	CP4712	1 5,0	75	12	Yellowish
	CP5001	2 6,0	75	12	Yellowish
	CP5118	3 7,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4367A	0 4,0	75	12	Yellowish
	CP4712A	1 5,0	75	12	Yellowish
	CP5001A	2 6,0	75	12	Yellowish
	CP5118A	3 7,0	75	12	Yellowish
1/2 Reverse Cutting With Alcohol 	CP5018A	2 6,0	75	12	Yellowish
3/8 Reverse Cutting 	CP4407	0 4,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP4407A	0 4,0	75	12	Yellowish
45 mm					
1/2 Round Body, Taper Point 	CP4421	0 4,0	75	12	Yellowish
	CP4771	1 5,0	75	12	Yellowish
	CP5034	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4771A	1 5,0	75	12	Yellowish
	CP5034A	2 6,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
45 mm					
3/8 Reverse Cutting 	CP4787	1 5,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP4787A	1 5,0	75	12	Yellowish
50 mm					
1/2 Round Body, Taper Point 	CP4453	0 4,0	75	12	Yellowish
	CP4811	1 5,0	75	12	Yellowish
	CP5065	2 6,0	75	12	Yellowish
	CP5127	3 7,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4453A	0 4,0	75	12	Yellowish
	CP4811A	1 5,0	75	12	Yellowish
	CP5065A	2 6,0	75	12	Yellowish
	CP5127A	3 7,0	75	12	Yellowish
Straight Reverse Cutting 	CP3475	3/0 3,0	75	12	Yellowish
Straight Reverse Cutting With Alcohol 	CP3475A	3/0 3,0	75	12	Yellowish
60 mm					
Straight Reverse Cutting 	CP4073	2/0 3,5	75	12	Yellowish
Straight Reverse Cutting With Alcohol 	CP4073A	2/0 3,5	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
75 cm					
PLAIN CATGUT Non-Needle Suture	CP1748	5/0 1,5	75	12	Yellowish
	CP2316	4/0 2,0	75	12	Yellowish
	CP2925	3/0 3,0	75	12	Yellowish
	CP4104	0 4,0	75	12	Yellowish
	CP4521	1 5,0	75	12	Yellowish
	CP4898	2 6,0	75	12	Yellowish
75 cm With Alcohol					
PLAIN CATGUT Non-Needle Suture With Alcohol	CP1748A	5/0 1,5	75	12	Yellowish
	CP2316A	4/0 2,0	75	12	Yellowish
	CP2925A	3/0 3,0	75	12	Yellowish
	CP3530A	2/0 3,5	75	12	Yellowish
	CP4104A	0 4,0	75	12	Yellowish
	CP4898A	2 6,0	75	12	Yellowish
150 cm					
PLAIN CATGUT Non-Needle Suture	CP1736	5/0 1,5	150	12	Yellowish
	CP2325	4/0 2,0	150	12	Yellowish
	CP2934	3/0 3,0	150	12	Yellowish
	CP3510	2/0 3,5	150	12	Yellowish
	CP4113	0 4,0	150	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
PLAIN CATGUT Non-Needle Suture	CP4531	1 5,0	150	12	Yellowish
	CP4907	2 6,0	150	12	Yellowish
	CP5106	3 7,0	150	12	Yellowish
	CP5134	4 8,0	150	12	Yellowish
150 cm With Alcohol					
PLAIN CATGUT Non-Needle Suture With Alcohol	CP1736A	5/0 4,5	150	12	Yellowish
	CP2325A	4/0 2,0	150	12	Yellowish
	CP2934A	3/0 3,0	150	12	Yellowish
	CP3510A	2/0 3,5	150	12	Yellowish
	CP4113A	0 3,5	150	12	Yellowish
	CP4531A	1 4,0	150	12	Yellowish
	CP4907A	2 6,0	150	12	Yellowish
	CP5106A	3 7,0	150	12	Yellowish
	CP5134A	4 8,0	150	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotos the product is "wet" e.g.

*After code number CP####A

Chromic Catgut

Absorbable Surgical Suture



Definition

Collagen based biomaterial, absorbable.

Material

Collagen tissue derived from purified sheep intestine

Coating

No coating

Structure

Monofilament

Colour

Brown

Size

EP : 1 to 8

USP : 6/0 to 4

Sterilization Method

Ethylene Oxide Gas

Gamma Irradiation




Main Indication

Ligatures

Ophthalmic Surgery

Chromic Catgut

Absorbable Surgical Suture








Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 430 Micron With Alchol			CK1892A	5/0 1,5	75	12	Brown
3/8 Reverse Cutting 430 Micron			CK1933	5/0 1,5	75	12	Brown
16 mm							
1/2 Round Body, Taper Point			CK1989	5/0 1,5	75	12	Brown
			CK2503	4/0 2,0	75	12	Brown
			CK3004	3/0 3,0	75	12	Brown
1/2 Round Body, Taper Point With Alchol			CK1989A	5/0 1,5	75	12	Brown
			CK2503A	4/0 2,0	75	12	Brown
			CK3004A	3/0 3,0	75	12	Brown
3/8 Reverse Cutting			CK2019	5/0 1,5	75	12	Brown
			CK2539	4/0 2,0	75	12	Brown
			CK3024	3/0 3,0	75	12	Brown
3/8 Reverse Cutting With Alchol			CK1670A	6/0 1,0	75	12	Brown
			CK2539A	4/0 2,0	75	12	Brown
			CK3024A	3/0 3,0	75	12	Brown
19 mm							
3/8 Reverse Cutting			CK2639	4/0 2,0	75	12	Brown
3/8 Reverse Cutting With Alchol			CK2205A	5/0 1,5	75	12	Brown
			CK2639A	4/0 2,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm					
1/2 Round Body, Taper Point  	CK2228	5/0 1,5	75	12	Brown
	CK2683	4/0 2,0	75	12	Brown
	CK3141	3/0 3,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol  	CK2683A	4/0 2,0	75	12	Brown
	CK3141A	3/0 3,0	75	12	Brown
3/8 Reverse Cutting  	CK2250	5/0 1,5	75	12	Brown
	CK2718	4/0 2,0	75	12	Brown
	CK3167	3/0 3,0	75	12	Brown
3/8 Reverse Cutting With Alcohol  	CK2250A	5/0 1,5	75	12	Brown
	CK2718A	4/0 2,0	75	12	Brown
	CK3167A	3/0 3,0	75	12	Brown
25 mm					
3/8 Reverse Cutting  	CK2800	4/0 2,0	75	12	Brown
3/8 Reverse Cutting With Alcohol  	CK2800A	4/0 2,0	75	12	Brown
26 mm					
1/2 Round Body, Taper Point  	CK2827	4/0 2,0	75	12	Brown
	CK3312	3/0 3,0	75	12	Brown
	CK3785	2/0 3,5	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point 	CK4214	0 4,0	75	12	Brown
	CK4589	1 5,0	75	12	Brown
	CK4946	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK2827A	4/0 2,0	75	12	Brown
	CK3312A	3/0 3,0	75	12	Brown
	CK3785A	2/0 3,5	75	12	Brown
	CK4214A	0 4,0	75	12	Brown
	CK4589A	1 5,0	75	12	Brown
3/8 Reverse Cutting 	CK2851	4/0 2,0	75	12	Brown
	CK3345	3/0 3,0	75	12	Brown
	CK3821	2/0 3,5	75	12	Brown
	CK4236	0 4,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK2851A	4/0 2,0	75	12	Brown
	CK3345A	3/0 3,0	75	12	Brown
	CK3821A	2/0 3,5	75	12	Brown
	CK4236A	0 4,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	CK3368	3/0 3,0	75	12	Brown
	CK3857	2/0 3,5	75	12	Brown
	CK4243	0 4,0	75	12	Brown
	CK4954	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK3857A	2/0 3,5	75	12	Brown
	CK4243A	0 4,0	75	12	Brown
	CK4607A	1 5,0	75	12	Brown
	CK4954A	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK3396	3/0 3,0	75	12	Brown
	CK3890	2/0 3,5	75	12	Brown
	CK4628	1 5,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK3396A	3/0 3,0	75	12	Brown
	CK3890A	2/0 3,5	75	12	Brown
	CK4272A	0 4,0	75	12	Brown
35 mm					
1/2 Round Body, Taper Point 	CK3410	3/0 3,0	75	12	Brown
	CK3915	2/0 3,5	75	12	Brown
	CK4299	0 4,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture


Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
1/2 Round Body, Taper Point 	CK4643	1 5,0	75	12	Brown
	CK4969	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4299A	0 4,0	75	12	Brown
	CK4643A	1 5,0	75	12	Brown
1/2 Reverse Cutting 	CK4976	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK3422	3/0 3,0	75	12	Brown
	CK4319	0 4,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4319A	0 4,0	75	12	Brown
40 mm					
1/2 Round Body, Taper Point 	CK4712	1 5,0	75	12	Brown
	CK5001	2 6,0	75	12	Brown
	CK5118	3 7,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4367A	0 4,0	75	12	Brown
	CK5001A	2 6,0	75	12	Brown
	CK5118A	3 7,0	75	12	Brown
1/2 Reverse Cutting 	CK5018	2 6,0	75	12	Brown
	CK5120	3 7,0	75	12	Brown
1/2 Reverse Cutting With Alcohol 	CK5120A	3 7,0	75	12	Brown


*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
3/8 Reverse Cutting 	CK4407	0 4,0	75	12	Brown
	CK4755	1 5,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4407A	0 4,0	75	12	Brown
	CK4755A	1 5,0	75	12	Brown





45 mm					
1/2 Round Body, Taper Point 	CK4421	0 4,0	75	12	Brown
	CK5034	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4421A	0 4,0	75	12	Brown
	CK4771A	1 5,0	75	12	Brown
	CK5034A	2 6,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4431A	0 4,0	75	12	Brown
Straight Reverse Cutting 	CK2911	4/0 2,0	75	12	Brown
Straight Reverse Cutting With Alcohol 	CK2911A	4/0 2,0	75	12	Brown

50 mm					
1/2 Round Body, Taper Point 	CK4811	1 5,0	75	12	Brown
	CK5065	2 6,0	75	12	Brown
	CK5127	3 7,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point With Alcohol 	CK4811A	1 5,0	75	12	Brown
	CK5065A	2 6,0	75	12	Brown
	CK5127A	3 7,0	75	12	Brown
1/2 Reverse Cutting 	CK5077	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK4474	0 4,0	75	12	Brown
	CK4840	1 5,0	75	12	Brown
	CK5082	2 6,0	75	12	Brown
	CK5130	3 7,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4840A	1 5,0	75	12	Brown
	CK5082A	2 6,0	75	12	Brown
	CK5130A	3 7,0	75	12	Brown
60 mm					
Straight Reverse Cutting 	CK3491	3/0 3,0	75	12	Brown
	CK4073	2/0 3,5	75	12	Brown
Straight Reverse Cutting With Alcohol 	CK4073A	2/0 3,5	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denetos the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
CHROMIC CATGUT Non-Needle Suture	CK1373	6/0 1,0	150	12	Brown
	CK1736	5/0 1,5	150	12	Brown
	CK2325	4/0 2,0	150	12	Brown
	CK2934	3/0 3,0	150	12	Brown
	CK3510	2/0 3,5	150	12	Brown
	CK4113	0 4,0	150	12	Brown
	CK4531	1 5,0	150	12	Brown
	CK4907	2 6,0	150	12	Brown
	CK5106	3 7,0	150	12	Brown
	CK5134	4 8,0	150	12	Brown

150 cm With Alcohol					
CHROMIC CATGUT Non-Needle Suture With Alcohol	CK1373A	6/0 1,0	150	12	Brown
	CK1736A	5/0 1,5	150	12	Brown
	CK2325A	4/0 2,0	150	12	Brown
	CK2934A	3/0 3,0	150	12	Brown
	CK3510A	2/0 3,5	150	12	Brown
	CK4113A	0 4,0	150	12	Brown
	CK4531A	1 5,0	150	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denetos the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm With Alcohol					
CHROMIC CATGUT Non-Needle Suture With Alcohol	CK4907A	2 6,0	150	12	Brown
	CK5106A	3 7,0	150	12	Brown
	CK5134A	4 8,0	150	12	Brown

75 cm					
CHROMIC CATGUT Non-Needle Suture	CK1748	5/0 1,5	75	12	Brown
	CK2316	4/0 2,0	75	12	Brown
	CK2925	3/0 3,0	75	12	Brown
	CK3530	2/0 3,5	75	12	Brown
	CK4104	0 4,0	75	12	Brown
	CK4521	1 5,0	75	12	Brown
	CK4898	2 6,0	75	12	Brown

75 cm With Alcohol					
CHROMIC CATGUT Non-Needle Suture With Alcohol	CK1748A	5/0 1,5	75	12	Brown
	CK2316A	4/0 2,0	75	12	Brown
	CK2925A	3/0 3,0	75	12	Brown
	CK3530A	2/0 3,5	75	12	Brown
	CK4104A	0 4,0	75	12	Brown
	CK4521A	1 5,0	75	12	Brown
	CK4898A	2 6,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A



Boz Medical Device Industry and Trade Inc.

Non - Absorbable Sutures

SILK

MONOPROLEN

MONOFLORID

MONAMID

MULTAMID

POLISIL

MONOTEF

POLYTEF

MONOWIRE

TEMPOWIRE

web : www.boztibbi.com
email: boz@boztibbi.com

Silk

Surgical Suture, Silk, Silicone Coated
Nonabsorbable, Noncapillary

Surgical Sutures Non-Absorbable

Boz®



Definition

Natural Non-absorbable

Material

Natural Silk

Coating

Silicon Coated, Paraffin Wax, Uncoated

Structure

Braided, Twisted

Colour

Black, Ivory, methylene blue

Size

EP : 0.2 to 7

USP : 10/0 to 5

Tensile Strenght Retention

Moderate Loss

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

Ligatures

Skin Closure

Oral Surgery




Neurosurgery




Cardiovascular procedures







Ophthalmic procedures



Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
6 mm							
3/8 Spatula 200 Micron			SB1120	9/0 0,3	30	12	Black
3/8 Spatula 220 Micron			SB1153	8/0 0,4	45	12	Black
6,55 mm							
3/8 Spatula 220 Micron			SB1174	8/0 0,4	45	12	Black
7,6 mm							
1/2 Reverse Cutting 330 Micron			SB5522	6/0 0,7	75	12	Black
8 mm							
1/2 Spatula 230 Micron			SB1197	8/0 0,4	45	12	Black
1/4 Spatula 330 Micron			SB1399	6/0 0,7	45	12	Black
			SB1761	5/0 1,0	45	12	Black
3/8 Round Body, Taper Point 230 Micron			SB1297	7/0 0,5	45	12	Black
3/8 Round Body, Taper Point 230 Micron			SB1209	8/0 0,4	45	12	Black
3/8 Reverse Cutting 280 Micron			SB1308	7/0 0,5	45	12	Black
3/8 Spatula 230 Micron			SB1201	8/0 0,4	45	12	Black
			SB1282	7/0 0,5	45	12	Black
3/8 Spatula 280 Micron			SB1414	6/0 0,7	45	12	Black
10 mm							
3/8 Round Body, Taper Point 330 Micron			SB1810	5/0 1,0	75	12	Black
3/8 Reverse Cutting 330 Micron			SB1495	6/0 0,7	45	12	Black


Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm					
3/8 Round Body, Taper Point 330 Micron  	SB1581	6/0 0,7	45	12	Black
	SB1591	6/0 0,7	75	12	Black
3/8 Round Body, Taper Point 380 Micron  	SB1913	5/0 1,0	75	12	Black
3/8 Round Body, Taper Point 430 Micron  	SB2437	4/0 1,5	75	12	Black
3/8 Reverse Cutting 330 Micron  	SB1613	6/0 0,7	45	12	Black
3/8 Reverse Cutting 330 Micron  	SB1617	6/0 0,7	45	12	Black
3/8 Reverse Cutting 430 Micron  	SB1931	5/0 1,0	45	12	Black
	SB2446	4/0 1,5	45	12	Black
16 mm					
1/2 Round Body, Taper Point  	SB1653	6/0 0,7	45	12	Black
	SB1993	5/0 1,0	75	12	Black
	SB2505	4/0 1,5	75	12	Black
	SB3006	3/0 2,0	75	12	Black
1/2 Reverse Cutting  	SB3019	3/0 2,0	75	12	Black
3/8 Reverse Cutting  	SB1667	6/0 0,7	45	12	Black
	SB1671	6/0 0,7	75	12	Black
	SB2015	5/0 1,0	45	12	Black

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Reverse Cutting			SB2023	5/0 1,0	75	12	Black
			SB2542	4/0 1,5	75	12	Black
3/8 Inside Cutting			SB1683	6/0 0,7	75	12	Black
			SB2145	5/0 1,0	75	12	Ivory
			SB2146	5/0 1,0	75	12	Black
			SB2558	4/0 1,5	75	12	Black
			SB3041	3/0 2,0	75	12	Black
18 mm							
1/2 Round Body, Taper Point			SB2176	5/0 1,0	75	12	Black
			SB2614	4/0 1,5	75	12	Black
			SB3082	3/0 2,0	75	12	Black
3/8 Round Body, Taper Point			SB4164	0 3,5	75	12	Black
3/8 Inside Cutting			SB2192	5/0 1,0	75	12	Black
			SB3094	3/0 2,0	75	12	Black
19 mm							
3/8 Reverse Cutting			SB5418	5/0 1,0	75	12	Black
			SB2641	4/0 1,5	75	12	Black
3/8 Inside Cutting			SB2662	4/0 1,5	75	12	Black

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			SB2231	5/0 1,0	75	12	Black
			SB2687	4/0 1,5	75	12	Black
			SB3145	3/0 2,0	75	12	Black
			SB3649	2/0 3,0	75	12	Black
1/2 Reverse Cutting			SB2711	4/0 1,5	75	12	Black
			SB3163	3/0 2,0	75	12	Black
1/2 Inside Cutting			SB2725	4/0 1,5	75	12	Black
3/8 Reverse Cutting			SB2251	5/0 1,0	75	12	Black
			SB2720	4/0 1,5	75	12	Black
			SB3171	3/0 2,0	75	12	Black
			SB3664	2/0 3,0	75	12	Black
3/8 Inside Cutting			SB2257	5/0 1,0	75	12	Black
			SB2732	4/0 1,5	75	12	Black
			SB3183	3/0 2,0	75	12	Black
22 mm							
1/2 Round Body, Taper Point			SB2749	4/0 1,5	75	12	Black
1/2 Inside Cutting			SB3221	3/0 2,0	75	12	Black
			SB3687	2/0 3,0	75	12	Black



Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm						
3/8 Reverse Cutting  	SB2807	4/0 1,5	75	12	Black	
	SB2279	5/0 1,0	45	12	Black	
	SB2282	5/0 1,0	75	12	Black	
	SB2813	4/0 1,5	45	12	Black	
	SB3286	3/0 2,0	75	12	Black	
3/8 Inside Cutting  	SB3764	2/0 3,0	75	12	Black	
26 mm						
1/2 Round Body, Taper Point  	SB2831	4/0 1,5	75	12	Black	
	SB3315	3/0 2,0	75	12	Black	
	SB3789	2/0 3,0	75	12	Black	
	SB4218	0 3,5	75	12	Black	
1/2 Reverse Cutting  	SB3337	3/0 2,0	75	12	Black	
	SB4231	0 3,5	75	12	Black	
3/8 Reverse Cutting  	SB2854	4/0 1,5	75	12	Black	
	SB3347	3/0 2,0	75	12	Black	
	SB3825	2/0 3,0	75	12	Black	
	SB4238	0 3,5	75	12	Black	

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	SB3371	3/0 2,0	75	12	Black
	SB3859	2/0 3,0	75	12	Black
	SB4245	0 3,5	75	12	Black
	SB4609	1 4,0	75	12	Black
1/2 Reverse Cutting 	SB3883	2/0 3,0	75	12	Black
	SB4267	0 3,5	75	12	Black
	SB4624	1 4,0	75	12	Black
3/8 Reverse Cutting 	SB2886	4/0 1,5	75	12	Black
	SB3399	3/0 2,0	75	12	Black
30 mm					
3/8 Reverse Cutting 	SB3893	2/0 3,0	75	12	Black
	SB4274	0 3,5	75	12	Black
	SB4630	1 4,0	75	12	Black
35 mm					
1/2 Round Body, Taper Point 	SB3917	2/0 3,0	75	12	Black
	SB4301	0 3,5	75	12	Black
	SB4646	1 4,0	75	12	Black

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
3/8 Reverse Cutting 	SB3425	3/0 2,0	75	12	Black
	SB3934	2/0 3,0	75	12	Black
	SB4663	1 4,0	75	12	Black
37 mm					
1/2 Round Body, Taper Point 	SB4339	0 3,5	75	12	Black
	SB4683	1 4,0	75	12	Black
1/2 Reverse Cutting 	SB4353	0 3,5	75	12	Black
	SB4695	1 4,0	75	12	Black
40 mm					
1/2 Round Body, Taper Point 	SB3983	2/0 3,0	75	12	Black
	SB4370	0 3,5	75	12	Black
	SB4716	1 4,0	75	12	Black
	SB5004	2 5,0	75	12	Black
1/2 Reverse Cutting 	SB4746	1 4,0	75	12	Black
	SB5020	2 5,0	75	12	Black
3/8 Reverse Cutting 	SB4409	0 3,5	75	12	Black
	SB4757	1 4,0	75	12	Black
	SB4759	1 4,0	90	12	Black

Silk

Surgical Suture, Silk, Silicone Coated Nonabsorbable, Noncapillary

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
45 mm					
1/2 Round Body, Taper Point 	SB5037	2 5,0	75	12	Black
3/8 Reverse Cutting 	SB4788	1 4,0	75	12	Black
Straight Reverse Cutting 	SB3463	3/0 2,0	75	12	Black

50 mm					
Straight Reverse Cutting 	SB4048	2/0 3,0	75	12	Black

60 mm					
Straight Reverse Cutting 	SB3493	3/0 2,0	75	12	Black
	SB4076	2/0 3,0	75	12	Black

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
10x75 cm					
SILK Non-Needle Suture	SB2340	4/0 1,5	10x75	12	Black
	SB2956	3/0 2,0	10x75	12	Black
	SB3549	2/0 3,0	10x75	12	Black
	SB4142	0 3,5	10x75	12	Black
	SB4552	1 4,0	10x75	12	Black
	SB4929	2 5,0	10x75	12	Black

12x45 cm					
SILK Non-Needle Suture	SB5840	0 3,5	12x45	12	Black
	SB5839	1 4,0	12x45	12	Black

13x60 cm					
SILK Non-Needle Suture	SB2348	4/0 1,5	13x60	12	Black
	SB2963	3/0 2,0	13x60	12	Black
	SB3555	2/0 3,0	13x60	12	Black
	SB4148	0 3,5	13x60	12	Black
	SB4558	1 4,0	13x60	12	Black

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
17x45 cm					
SILK Non-Needle Suture	SB1740	5/0 1,0	17x45	12	Black
	SB2351	4/0 1,5	17x45	12	Black
	SB2964	3/0 2,0	17x45	12	Black
	SB3560	2/0 3,0	17x45	12	Black
	SB4151	0 3,5	17x45	12	Black
	SB4938	2 5,0	17x45	12	Black
	SB4561	1 4,0	17x45	12	Black

150 cm					
SILK Non-Needle Suture	SB2936	3/0 2,0	150	12	Black

180 cm					
SILK Non-Needle Suture	SB2940	3/0 2,0	180	12	Black
	SB3516	2/0 3,0	180	12	Black
	SB4119	0 3,5	180	12	Black
	SB4539	1 4,0	180	12	Black
	SB4914	2 5,0	180	12	Black
	SB5107	3 6,0	180	12	Black
	SB5133	4 6,0	180	12	Black

Monoprolen

Surgical Suture, Polypropylene (PP)
Nonabsorbable, Synthetic, Monofilament

Boz®



Definition

Synthetic Non-Absorbable

Material

Polypropylene (PP)

Coating

No Coating

Structure

Monofilament

Colour

Blue, Black, Undyed

Size

EP : 0.2 to 5

USP : 10/0 to 2

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

Microsugery

Skin closure

Cosmetic Surgery

Cardiac Surgery




Vascular Surgery

Ophthalmic procedures

Neurosurgical procedures

Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
5 mm							
3/8 Round Body, Taper Point 130 Micron			ML1047	10/0 0,2	13	12	Black
5.1 mm							
3/8 Round Body, Taper Point 100 Micron			ML1110	9/0 0,3	13	12	Blue
6 mm							
3/8 Round Body, Taper Point 130 Micron			ML1161	8/0 0,4	60	12	Blue
3/8 Round Body, Taper Point 220 Micron			ML1255	7/0 0,5	75	12	Blue
3/8 Round Body, Taper Point 150 Micron			ML1387	6/0 0,7	45	12	Blue
6.5 mm							
3/8 Spatula 130 Micron			ML1172	8/0 0,4	45	12	Blue
8 mm							
3/8 Round Body, Taper Point 220 Micron			ML1292	7/0 0,5	60	12	Blue
3/8 Round Body, Taper Point 220 Micron Black Needle			ML1422	6/0 0,7	60	12	Blue
3/8 Round Body, Taper Point 220 Micron Black Needle			ML1293	7/0 0,5	60	12	Blue
3/8 Round Body, Taper Point 230 Micron			ML1212	8/0 0,4	60	12	Blue
			ML1213	8/0 0,4	75	12	Blue
3/8 Round Body, Taper Point 230 Micron Black Needle			ML1211	8/0 0,4	60	12	Blue
3/8 Round Body, Taper Point 280 Micron			ML1300	7/0 0,5	45	12	Blue
			ML1305	7/0 0,5	60	12	Blue

Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm					
3/8 Round Body, Taper Point 280 Micron 	ML1425	6/0 0,7	45	12	Blue
	ML1430	6/0 0,7	60	12	Blue
3/8 Spatula 220 Micron 	ML1284	7/0 0,5	75	12	Blue













9 mm					
3/8 Round Body, Taper Point 230 Micron 	ML1313	7/0 0,5	75	12	Blue

10 mm					
1/2 Round Body, Taper Point 280 Micron 	ML1463	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 220 Micron 	ML1326	7/0 0,5	60	12	Blue
3/8 Round Body, Taper Point 230 Micron Black Needle 	ML1324	7/0 0,5	75	12	Blue
3/8 Round Body, Taper Point 280 Micron 	ML1477	6/0 0,7	60	12	Blue
	ML1483	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 280 Micron Black Needle 	ML1476	6/0 0,7	60	12	Blue
3/8 Round Body, Taper Point 330 Micron 	ML1801	5/0 1,0	60	12	Blue
	ML1802	5/0 1,0	75	12	Blue
3/8 Reverse Cutting 340 Micron 	ML1497	6/0 0,7	45	12	Blue

12 mm					
3/8 Round Body, Taper Point 280 Micron 	ML1346	7/0 0,5	75	12	Blue





Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 320 Micron			ML1882	5/0 1,0	75	12	Blue
1/2 Round Body, Taper Point 330 Micron			ML1565	6/0 0,7	75	12	Blue
1/2 Round Body, Taper Point 330 Micron			ML1566	6/0 0,7	75	12	Blue
			ML1569	6/0 0,7	90	12	Blue
1/2 Round Body, Taper Point 380 Micron			ML1884	5/0 1,0	75	12	Blue
1/2 Round Body, Taper Point 430 Micron			ML2427	4/0 1,5	75	12	Blue
1/2 Round Body, Taper Point 430 Micron			ML1898	5/0 1,0	75	12	Blue
			ML2428	4/0 1,5	75	12	Blue
3/8 Round Body, Taper Point 230 Micron			ML1363	7/0 0,5	60	12	Blue
3/8 Round Body, Taper Point 280 Micron			ML1571	6/0 0,7	60	12	Blue
3/8 Round Body, Taper Point 280 Micron Black Needle			ML1365	7/0 0,5	75	12	Blue
			ML1572	6/0 0,7	60	12	Blue
3/8 Round Body, Taper Point 330 Micron			ML1582	6/0 0,7	45	12	Blue
			ML1592	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 330 Micron			ML1584	6/0 0,7	60	12	Blue
			ML1588	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 330 Micron Black Needle			ML1593	6/0 0,7	90	12	Blue






Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
3/8 Round Body, Taper Point 380 Micron			ML1910	5/0 1,0	60	12	Blue
3/8 Round Body, Taper Point 380 Micron			ML1917	5/0 1,0	75	12	Blue
3/8 Reverse Cutting 330 Micron			ML1612	6/0 0,7	45	12	Blue
3/8 Reverse Cutting 330 Micron			ML1621	6/0 0,7	75	12	Blue
3/8 Reverse Cutting 380 Micron			ML1927	5/0 1,0	45	12	Blue
3/8 Reverse Cutting 560 Micron			ML2454	4/0 1,5	45	12	Blue
3/8 Spatula 280 Micron			ML1606	6/0 0,7	75	12	Blue
16 mm							
1/2 Round Body, Taper Point			ML1987	5/0 1,0	45	12	Blue
			ML2506	4/0 1,5	75	12	Blue
			ML3007	3/0 2,0	75	12	Blue
1/2 Round Body, Taper Point			ML1658	6/0 0,7	90	12	Blue
			ML1998	5/0 1,0	90	12	Blue
			ML2514	4/0 1,5	90	12	Blue
			ML3010	3/0 2,0	90	12	Blue
1/2 Round Body, Taper Point Black Needle			ML1997	5/0 1,0	90	12	Blue
			ML2499	4/0 1,5	45	12	Blue








Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Round Body, Taper Point			ML2007	5/0 1,0	75	12	Blue
3/8 Round Body, Taper Point			ML2009	5/0 1,0	90	12	Blue
			ML2524	4/0 1,5	90	12	Blue
3/8 Reverse Cutting			ML1669	6/0 0,7	45	12	Blue
			ML1672	6/0 0,7	75	12	Blue
			ML2017	5/0 1,0	45	12	Blue
			ML2126	5/0 1,0	90	12	Blue
			ML2536	4/0 1,5	45	12	Blue
			ML2541	4/0 1,5	75	12	Blue
			ML3027	3/0 2,0	75	12	Blue
3/8 Inside Cutting			ML1681	6/0 0,7	45	12	Blue
			ML2141	5/0 1,0	45	12	Blue
			ML2556	4/0 1,5	45	12	Blue
			ML3041	3/0 2,0	75	12	Blue
17 mm							
1/2 Round Body, Taper Point			ML1690	6/0 0,7	75	12	Blue
			ML2164	5/0 1,0	75	12	Blue
			ML3060	3/0 2,0	75	12	Blue









Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
17 mm							
3/8 Round Body, Taper Point			ML2169	5/0 1,0	90	12	Blue
			ML2595	4/0 1,5	90	12	Blue
18 mm							
1/2 Round Body, Taper Point			ML2615	4/0 1,5	75	12	Blue
19 mm							
3/8 Reverse Cutting			ML2207	5/0 1,0	75	12	Blue
			ML2642	4/0 1,5	75	12	Blue
20 mm							
1/2 Round Body, Taper Point			ML2230	5/0 1,0	75	12	Blue
			ML2684	4/0 1,5	75	12	Blue
			ML3143	3/0 2,0	75	12	Blue
			ML3647	2/0 3,0	75	12	Blue
1/2 Round Body, Taper Point			ML2234	5/0 1,0	90	12	Blue
			ML2696	4/0 1,5	90	12	Blue
			ML3155	3/0 2,0	90	12	Blue
1/2 Round Body, Taper Point Black Needle			ML2697	4/0 1,5	90	12	Blue
3/8 Reverse Cutting			ML2249	5/0 1,0	45	12	Blue
			ML2715	4/0 1,5	45	12	Blue



Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
3/8 Reverse Cutting			ML2721	4/0 1,5	75	12	Blue
			ML3166	3/0 2,0	45	12	Blue
			ML3169	3/0 2,0	75	12	Blue
3/8 Inside Cutting			ML2253	5/0 1,0	45	12	Blue
			ML2726	4/0 1,5	45	12	Blue
			ML3178	3/0 2,0	45	12	Blue
			ML3182	3/0 2,0	75	12	Blue
22 mm							
1/2 Round Body, Taper Point			ML2267	5/0 1,0	75	12	Blue
1/2 Round Body, Taper Point			ML2753	4/0 1,5	90	12	Blue
1/2 Inside Cutting			ML3218	3/0 2,0	45	12	Blue
3/8 Round Body, Taper Point			ML2270	5/0 1,0	90	12	Blue
			ML2758	4/0 1,5	90	12	Blue
25 mm							
1/2 Round Body, Taper Point			ML2785	4/0 1,5	45	12	Blue
			ML3253	3/0 2,0	75	12	Blue
			ML3725	2/0 3,0	75	12	Blue
			ML4194	0 3,5	100	12	Blue


Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	ML2796	4/0 1,5	90	12	Blue
	ML3261	3/0 2,0	90	12	Blue
	ML3735	2/0 3,0	75	12	Blue
	ML3744	2/0 3,0	90	12	Blue
1/2 Round Body, Taper Point Black Needle 	ML2795	4/0 1,5	90	12	Blue
	ML3260	3/0 2,0	90	12	Blue
3/8 Reverse Cutting 	ML2806	4/0 1,5	75	12	Blue
	ML3272	3/0 2,0	45	12	Blue
3/8 Inside Cutting 	ML2811	4/0 1,5	45	12	Blue
	ML3283	3/0 2,0	45	12	Blue
	ML3761	2/0 3,0	45	12	Blue
26 mm					
1/2 Round Body, Taper Point 	ML2830	4/0 1,5	75	12	Blue
	ML3787	2/0 3,0	75	12	Blue
	ML4216	0 3,5	75	12	Blue
	ML4591	1 4,0	75	12	Blue
3/8 Reverse Cutting 	ML2853	4/0 1,5	75	12	Blue
	ML3348	3/0 2,0	75	12	Blue
	ML3823	2/0 3,0	75	12	Blue



Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	ML2875	4/0 1,5	75	12	Blue
	ML3372	3/0 2,0	75	12	Blue
	ML3861	2/0 3,0	75	12	Blue
	ML4246	0 3,5	75	12	Blue
	ML4615	1 4,0	100	12	Blue
	ML4610	1 4,0	75	12	Blue
1/2 Round Body, Taper Point 	ML3378	3/0 2,0	90	12	Blue
	ML3867	2/0 3,0	90	12	Blue
1/2 Reverse Cutting 	ML3390	3/0 2,0	75	12	Blue
	ML3884	2/0 3,0	75	12	Blue
3/8 Reverse Cutting 	ML3398	3/0 2,0	75	12	Blue
	ML3894	2/0 3,0	75	12	Blue
	ML4276	0 3,5	75	12	Blue
	ML4631	1 4,0	75	12	Blue
35 mm					
1/2 Round Body, Taper Point 	ML3413	3/0 2,0	75	12	Blue
	ML3919	2/0 3,0	75	12	Blue






Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
1/2 Round Body, Taper Point 	ML4302	0 3,5	75	12	Blue
	ML4652	1 4,0	100	12	Blue
3/8 Reverse Cutting 	ML3426	3/0 2,0	75	12	Blue
	ML3936	2/0 3,0	75	12	Blue
	ML4321	0 3,5	75	12	Blue
	ML4662	1 4,0	75	12	Blue
40 mm					
1/2 Round Body, Taper Point 	ML3990	2/0 3,0	100	12	Blue
	ML4377	0 3,5	100	12	Blue
	ML4371	0 3,5	75	12	Blue
	ML4726	1 4,0	100	12	Blue
	ML4715	1 4,0	75	12	Blue
	ML5011	2 5,0	100	12	Blue
	ML5002	2 5,0	75	12	Blue
1/2 Round Body, Taper Point LOOP 	ML4387	0 3,5	150	12	Blue
	ML4737	1 4,0	150	12	Blue
1/2 Round Body, Taper Point Heavy 	ML4733	1 4,0	100	12	Blue


Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Reverse Cutting 	ML4402	0 3,5	100	12	Blue
	ML4751	1 4,0	100	12	Blue
3/8 Reverse Cutting 	ML4009	2/0 3,0	100	12	Blue
	ML4413	0 3,5	100	12	Blue
	ML4760	1 4,0	100	12	Blue
45 mm					
1/2 Round Body, Taper Point LOOP 	ML4781	1 4,0	150	12	Blue
Straight Reverse Cutting 	ML4028	2/0 3,0	75	12	Blue
50 mm					
1/2 Round Body, Taper Point 	ML4820	1 4,0	100	12	Blue
	ML5071	2 5,0	100	12	Blue
1/2 Round Body, Taper Point LOOP 	ML4465	0 3,5	150	12	Blue
	ML4829	1 4,0	150	12	Blue
1/2 Reverse Cutting 	ML4837	1 4,0	100	12	Blue
60 mm					
Straight Reverse Cutting 	ML3495	3/0 2,0	75	12	Blue
	ML4077	2/0 3,0	75	12	Blue

Monoprolen

Surgical Suture, Polypropylene (PP) Nonabsorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
90 mm					
 3/8 Inside Cutting Vinyl Tube ØOutside: 4 mm ØInside: 2.5 mm Length: 40 mm	ML4887	1 4,0	2 X 50	12	Blue

Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Boz®



Definition

Synthetic Non-Absorbable

Material

Polyvinylidene difluoride (PVDF)

Coating

No Coating

Structure

Monofilament

Colour

Blue, Black

Size

EP : 0.2 to 5

USP : 10/0 to 2

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

Cosmetic Surgery

Skin closure

Cardiovascular procedures

Ophthalmic procedures

Neurosurgical procedures








Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
9 mm							
3/8 Round Body, Taper Point 230 Micron Black Needle			MD1311	7/0 0,5	60	12	Blue
10 mm							
3/8 Round Body, Taper Point 220 Micron			MD1325	7/0 0,5	45	12	Blue
13 mm							
1/2 Round Body, Taper Point 330 Micron			MD1565	6/0 0,7	75	12	Blue
1/2 Round Body, Taper Point 330 Micron			MD1566	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 280 Micron			MD1363	7/0 0,5	60	12	Blue
3/8 Round Body, Taper Point 330 Micron Black Needle			MD1583	6/0 0,7	45	12	Blue
3/8 Round Body, Taper Point 330 Micron Black Needle			MD1589	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 330 Micron			MD1584	6/0 0,7	60	12	Blue
			MD1588	6/0 0,7	75	12	Blue
3/8 Round Body, Taper Point 380 Micron Black Needle			MD1907	5/0 1,0	45	12	Blue
3/8 Round Body, Taper Point 380 Micron			MD1912	5/0 1,0	60	12	Blue
			MD1917	5/0 1,0	75	12	Blue
16 mm							
1/2 Round Body, Taper Point			MD1654	6/0 0,7	60	12	Blue
3/8 Reverse Cutting			MD1669	6/0 0,7	45	12	Blue
			MD2017	5/0 1,0	45	12	Blue
			MD2536	4/0 1,5	45	12	Blue



Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Reverse Cutting			MD2541	4/0 1,5	75	12	Blue
			MD3023	3/0 2,0	45	12	Blue
3/8 Inside Cutting			MD2141	5/0 1,0	45	12	Blue
			MD2556	4/0 1,5	45	12	Blue
17 mm							
1/2 Round Body, Taper Point			MD2164	5/0 1,0	75	12	Blue
			MD2585	4/0 1,5	75	12	Blue
3/8 Round Body, Taper Point			MD2169	5/0 1,0	90	12	Blue
18 mm							
1/2 Round Body, Taper Point			MD2180	5/0 1,0	90	12	Blue
			MD2616	4/0 1,5	90	12	Blue
20 mm							
1/2 Round Body, Taper Point			MD2230	5/0 1,0	75	12	Blue
			MD2684	4/0 1,5	75	12	Blue
			MD3143	3/0 2,0	75	12	Blue
			MD3647	2/0 3,0	75	12	Blue
1/2 Round Body, Taper Point			MD2234	5/0 1,0	90	12	Blue
			MD2696	4/0 1,5	90	12	Blue
			MD3149	3/0 2,0	90	12	Blue

Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
3/8 Reverse Cutting			MD2715	4/0 1,5	45	12	Blue
			MD3169	3/0 2,0	75	12	Blue
3/8 Inside Cutting			MD2726	4/0 1,5	45	12	Blue
			MD3182	3/0 2,0	75	12	Blue
25 mm							
1/2 Round Body, Taper Point			MD3253	3/0 2,0	75	12	Blue
			MD3725	2/0 3,0	75	12	Blue
3/8 Reverse Cutting			MD3275	3/0 2,0	75	12	Blue
			MD3758	2/0 3,0	75	12	Blue
26 mm							
1/2 Round Body, Taper Point			MD2830	4/0 1,5	75	12	Blue
			MD3320	3/0 2,0	90	12	Blue
1/2 Round Body, Taper Point			MD2836	4/0 1,5	90	12	Blue
			MD3323	3/0 2,0	90	12	Blue
			MD3796	2/0 3,0	90	12	Blue
1/2 Reverse Cutting			MD2300	5/0 1,0	75	12	Blue
30 mm							
1/2 Round Body, Taper Point			MD3372	3/0 2,0	75	12	Blue
			MD3861	2/0 3,0	75	12	Blue



Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
1/2 Round Body, Taper Point		MD4246	0 3,5	75	12	Blue
		MD4610	1 4,0	75	12	Blue
1/2 Reverse Cutting		MD3884	2/0 3,0	75	12	Blue
3/8 Reverse Cutting		MD3398	3/0 2,0	75	12	Blue
		MD3894	2/0 3,0	75	12	Blue
		MD4276	0 3,5	75	12	Blue
35 mm						
1/2 Round Body, Taper Point		MD4302	0 3,5	75	12	Blue
3/8 Reverse Cutting		MD4321	0 3,5	75	12	Blue
40 mm						
1/2 Round Body, Taper Point		MD4377	0 3,5	100	12	Blue
		MD4371	0 3,5	75	12	Blue
		MD4726	1 4,0	100	12	Blue
		MD4715	1 4,0	75	12	Blue
1/2 Round Body, Taper Point, LOOP		MD4387	0 3,5	150	12	Blue
		MD4737	1 4,0	150	12	Blue
1/2 Reverse Cutting		MD4402	0 3,5	100	12	Blue
		MD4751	1 4,0	100	12	Blue

Monoflorid

Surgical Suture, Poly(vinylidene difluoride) (PVDF)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm						
3/8 Reverse Cutting 	MD4006	2/0 3,0	75	12	Blue	
45 mm						
1/2 Round Body, Taper Point, LOOP 	MD4427	0 3,5	150	12	Blue	
50 mm						
1/2 Round Body, Taper Point 	MD4813	1 4,0	75	12	Blue	
1/2 Round Body, Taper Point, LOOP 	MD4465	0 3,5	150	12	Blue	
	MD4829	1 4,0	150	12	Blue	
60 mm						
Straight Reverse Cutting 	MD3495	3/0 2,0	75	12	Blue	


Surgical Sutures

Non-Absorbable

Monamid

Surgical Suture, Monofilament Polyamide (PA)
Nonabsorbable, Synthetic

Boz®



Definition

Synthetic Non-Absorbable

Material

Nylon 6 and Nylon 6.6 (PA)

Coating

Uncoated

Structure

Monofilament

Colour

Blue, Black, Undyed

Size

EP : 0.2 to 5

USP : 10/0 to 2

Tensile Strenght Retention

Moderate Loss

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

General Soft tissue approximation or ligation

Cardiovascular

Ophthalmic Surgery

Neurosurgical Surgery

Plastic Surgery

Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
4 mm							
3/8 Round Body,Taper Point 70 Micron			MA1001	11/0 0,1	10	12	Black
3/8 Round Body,Taper Point 100 Micron			MA1023	10/0 0,2	15	12	Black
3/8 Round Body,Taper Point 150 Micron			MA1099	9/0 0,3	30	12	Black
			MA1129	8/0 0,4	13	12	Black
5 mm							
3/8 Round Body,Taper Point 100 Micron			MA1045	10/0 0,2	13	12	Black
			MA1107	9/0 0,3	13	12	Black
3/8 Round Body,Taper Point 130 Micron			MA1133	8/0 0,4	13	12	Black
6 mm							
3/8 Round Body,Taper Point 140 Micron			MA1162	8/0 0,4	15	12	Black
3/8 Spatula 150 Micron Micropoint			MA1065	10/0 0,2	30	12	Black
3/8 Spatula 150 Micron			MA1066	10/0 0,2	30	12	Black
			MA1068	10/0 0,2	45	12	Black
			MA1118	9/0 0,3	30	12	Black
6,6 mm							
1/2 Cobra 150 Micron			MA5520	10/0 0,2	60	12	Blue
			MA5521	9/0 0,3	60	12	Blue
10 mm							
1/2 Reverse Cutting 330 Micron			MA1488	6/0 0,7	45	12	Black
3/8 Reverse Cutting 330 Micron			MA1495	6/0 0,7	45	12	Black


Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
13 mm							
3/8 Reverse Cutting 330 Micron			MA1612	6/0 0,7	45	12	Blue
			MA1619	6/0 0,7	75	12	Blue
16 mm							
1/2 Reverse Cutting			MA2012	5/0 1,0	75	12	Blue
3/8 Reverse Cutting			MA1672	6/0 0,7	75	12	Blue
			MA1667	6/0 0,7	45	12	Black
			MA2017	5/0 1,0	45	12	Blue
			MA2536	4/0 1,5	45	12	Blue
			MA2542	4/0 1,5	75	12	Black
			MA3021	3/0 2,0	45	12	Black
			MA3027	3/0 2,0	75	12	Blue
19 mm							
3/8 Reverse Cutting			MA2201	5/0 1,0	45	12	Blue
			MA2207	5/0 1,0	75	12	Blue
			MA2638	4/0 1,5	45	12	Blue
			MA2642	4/0 1,5	75	12	Blue
3/8 Inside Cutting			MA3120	3/0 2,0	45	12	Blue
			MA3121	3/0 2,0	45	12	Black

Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
19 mm					
3/8 Inside Cutting  	MA3125	3/0 2,0	75	12	Blue
	MA2659	4/0 1,5	45	12	Blue
	MA2663	4/0 1,5	75	12	Blue
20 mm					
1/2 Reverse Cutting  	MA3162	3/0 2,0	75	12	Blue
3/8 Reverse Cutting  	MA3165	3/0 2,0	45	12	Black
24 mm					
3/8 Inside Cutting  	MA3235	3/0 2,0	45	12	Blue
25 mm					
1/2 Tapercutting  	MA5513	3/0 2,0	45	12	Blue
3/8 Reverse Cutting  	MA2802	4/0 1,5	45	12	Blue
	MA3758	2/0 3,0	75	12	Blue
3/8 Inside Cutting  	MA2811	4/0 1,5	45	12	Blue
26 mm					
1/2 Reverse Cutting  	MA3812	2/0 3,0	75	12	Blue
3/8 Reverse Cutting  	MA3344	3/0 2,0	45	12	Blue
	MA3347	3/0 2,0	75	12	Black
	MA3348	3/0 2,0	75	12	Blue


Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
3/8 Reverse Cutting 	MA3820	2/0 3,0	45	12	Blue
	MA3823	2/0 3,0	75	12	Blue
3/8 Inside Cutting 	MA3358	3/0 2,0	45	12	Blue
30 mm					
1/2 Round Body, Taper Point 	MA4254	0 3,5	100	12	Blue
	MA4255	0 3,5	100	12	Black
	MA4615	1 4,0	100	12	Blue
	MA4610	1 4,0	75	12	Blue
1/2 Reverse Cutting 	MA3392	3/0 2,0	90	12	Blue
	MA3884	2/0 3,0	75	12	Blue
3/8 Reverse Cutting 	MA3398	3/0 2,0	75	12	Blue
	MA3399	3/0 2,0	75	12	Black
	MA3893	2/0 3,0	75	12	Black
	MA3894	2/0 3,0	75	12	Blue
	MA4276	0 3,5	75	12	Blue
	MA4631	1 4,0	75	12	Blue
3/8 Inside Cutting 	MA3902	2/0 3,0	75	12	Blue






Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
1/2 Round Body, Taper Point 	MA4309	0 3,5	100	12	Black
	MA4298	0 3,5	45	12	Blue
	MA4652	1 4,0	100	12	Black
1/2 Reverse Cutting 	MA4317	0 3,5	75	12	Black
3/8 Reverse Cutting 	MA3936	2/0 3,0	75	12	Blue
	MA4664	1 4,0	100	12	Blue
38 mm					
1/2 Reverse Cutting 	MA5512	2/0 3,0	75	12	Black
	MA5510	0 3,5	75	12	Black
40 mm					
1/2 Round Body, Taper Point 	MA4377	0 3,5	100	12	Blue
	MA4370	0 3,5	75	12	Black
	MA4725	1 4,0	100	12	Black
	MA4726	1 4,0	100	12	Blue
1/2 Round Body, Taper Point LOOP 	MA4387	0 3,5	150	12	Blue
1/2 Reverse Cutting 	MA4401	0 3,5	100	12	Black
	MA4402	0 3,5	100	12	Blue

Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Reverse Cutting Heavy 	MA4405	0 3,5	100	12	Black
3/8 Reverse Cutting 	MA4758	1 4,0	75	12	Blue
45 mm					
1/2 Round Body, Taper Point, LOOP 	MA4781	1 4,0	150	12	Blue
3/8 Reverse Cutting 	MA4789	1 4,0	75	12	Blue
50 mm					
1/2 Round Body, Taper Point LOOP 	MA4464	0 3,5	150	12	Black
	MA4465	0 3,5	150	12	Blue
	MA4828	1 4,0	150	12	Black
	MA4829	1 4,0	150	12	Blue
60 mm					
Straight Reverse Cutting 	MA4077	2/0 3,0	75	12	Blue
61 mm					
Straight Reverse Cutting 	MA4083	2/0 3,0	100	12	Blue

Surgical Sutures Non-Absorbable

Multamid

Surgical Suture, Multifilament Polyamide (PA) Silicone Coated
Nonabsorbable, Synthetic, Black

Boz®



Definition

Synthetic Non-Absorbable

Material

Nylon 6 and Nylon 6.6 (PA)

Coating

Silicone, Paraffin Wax Uncoated

Structure

Braided

Colour

Black, Blue, Undyed

Size

EP : 0.01 to 7

USP : 12/0 to 5

Tensile Strenght Retention

Moderate Loss

Mass Absorption

Permanent

Sterilization Method




Ethylene Oxide Gas

Main Indication

General Soft tissue
approximation or ligation
Cardiovascular
Ophthalmic Surgery
Neurosurgical Surgery
Plastic Surgery


Multamid

Surgical Suture, Multifilament Polyamide (PA)
Silicone Coated Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body,Taper Point			MU3006	3/0 2,0	75	12	Black
26 mm							
1/2 Round Body,Taper Point			MU2824	4/0 1,5	8 X 45	6	Black
			MU3315	3/0 2,0	75	12	Black
			MU3789	2/0 3,0	75	12	Black
			MU3792	2/0 3,0	8 X 75	6	Black
			MU4218	0 3,5	75	12	Black
			MU4590	1 4,0	75	12	Black
3/8 Reverse Cutting			MU3343	3/0 2,0	45	12	Black
			MU4238	0 3,5	75	12	Black
30 mm							
1/2 Round Body,Taper Point			MU3371	3/0 2,0	75	12	Black
			MU3859	2/0 3,0	75	12	Black
			MU4609	1 4,0	75	12	Black
35 mm							
1/2 Round Body,Taper Point			MU4301	0 3,5	75	12	Black
3/8 Reverse Cutting			MU3425	3/0 2,0	75	12	Black
			MU3934	2/0 3,0	75	12	Black
			MU4663	1 4,0	75	12	Black

Multamid

Surgical Suture, Multifilament Polyamide (PA)
Silicone Coated Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
90 mm					
 ØOutside: 4 mm ØInside: 2.5 mm Length: 40 mm 3/8 Inside Cutting Vinyl Tube	MU4886	1 4,0	2 X 50	12	Black

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit	Color
17x45 cm					
MULTAMID Non-Needle Suture	MU1740	5/0 1,0	17x45	12	Black
	MU2351	4/0 1,5	17x45	12	Black
	MU2964	3/0 2,0	17x45	12	Black
	MU3560	2/0 3,0	17x45	12	Black
	MU4561	1 4,0	17x45	12	Black


10x75 cm					
MULTAMID Non-Needle Suture	MU2340	4/0 1,5	10x75	12	Black
	MU2956	3/0 2,0	10x75	12	Black
	MU3549	2/0 3,0	10x75	12	Black
	MU4552	1 4,0	10x75	12	Black

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET)
Nonabsorbable, Synthetic

Surgical Sutures
Non-Absorbable

Boz®



Definition

Synthetic Non- Absorbable

Material

Polyester Polyethylene Terephthalate (PET)

Coating

Silicone, Uncoated

Structure

Braided

Colour

Green, White

Size

EP : 0.5 to 7
USP : 7/0 to 5

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method




Ethylene Oxide Gas

Main Indication

Cardiac Valve Repair
General Surgery
Cardiovascular Procedures
Neurosurgical Procedures









Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/2 Spatula 340 Micron			PS1756	5/0 1,0	45	12	Green
1/4 Spatula 340 Micron			PS1408	6/0 0,7	45	12	Green
			PS1765	5/0 1,0	45	12	White
			PS1766	5/0 1,0	45	12	Green
3/8 Round Body, Taper Point 280 Micron			PS1302	7/0 0,5	45	12	Green
8,1 mm							
1/4 Sabre 380 Micron			PS5531	4/0 1,5	60	12	Green
10 mm							
3/8 Round Body, Taper Point 280 Micron			PS1482	6/0 0,7	75	12	Green
3/8 Spatula 280 Micron			PS1452	6/0 0,7	45	12	Green
13 mm							
3/8 Round Body, Taper Point 280 Micron			PS1573	6/0 0,7	75	12	Green
3/8 Round Body, Taper Point 320 Micron			PS1903	5/0 1,0	75	12	Green
3/8 Round Body, Taper Point 330 Micron			PS1587	6/0 0,7	75	12	Green
3/8 Round Body, Taper Point 330 Micron 6mm x 3mm x 3mm Teflon Pledget			PS1916	5/0 1,0	75	12	Green
3/8 Reverse Cutting 560 Micron			PS2455	4/0 1,5	45	12	White
3/8 Spatula 280 Micron			PS1604	6/0 0,7	45	12	Green








Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body,Taper Point			PS3008	3/0 2,0	75	12	Green
1/2 Round Body,Taper Point			PS2515	4/0 1,5	90	12	Green
			PS3011	3/0 2,0	90	12	Green
1/2 Round Body,Taper Point			PS3012	3/0 2,0	90	6	5 Green 5 White
			PS5923	3/0 2,0	90	12	Green
6mm x 3mm x 1.5mm Teflon Pledget 							
3/8 Inside Cutting			PS2560	4/0 1,5	75	12	Green
			PS3040	3/0 2,0	75	12	Green
17 mm							
1/2 Round Body,Taper Point			PS2582	4/0 1,5	75	12	Green
			PS3609	2/0 3,0	75	12	Green
1/2 Round Body,Taper Point			PS2587	4/0 1,5	90	12	Green
			PS3066	3/0 2,0	90	12	Green
			PS3613	2/0 3,0	90	12	Green
			PS3615	2/0 3,0	90	6	5 Green 5 White
1/2 Round Body,Taper Point			PS2579	4/0 1,5	75	12	Green
			PS3606	2/0 3,0	75	6	5 Green 5 White
			3mm x 3mm x 1.5mm Teflon Pledget 			PS3617	2/0 3,0

Polisil





Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
17 mm							
1/2 Round Body, Taper Point 6mm x 3mm x 1.5mm Teflon Pledget	 		PS3067	3/0 2,0	90	6	5 Green 5 White
			PS3616	2/0 3,0	90	6	5 Green 5 White
			PS3618	2/0 3,0	90	12	Green
1/2 Tapercutting			PS2570	4/0 1,5	90	12	Green
1/2 Tapercutting 7mm x 3mm x 1.5mm Teflon Pledget			PS3592	2/0 3,0	75	6	5 Green 5 White
3/8 Round Body, Taper Point			PS5762	3/0 2,0	75	12	Green
20 mm							
1/2 Round Body, Taper Point			PS3648	2/0 3,0	75	12	Green
1/2 Round Body, Taper Point			PS3650	2/0 3,0	75	12	Green
1/2 Round Body, Taper Point			PS2688	4/0 1,5	75	12	Green
			PS2695	4/0 1,5	90	12	Green
			PS5843	3/0 2,0	75	12	Green
			PS3148	3/0 2,0	90	12	Green
1/2 Round Body, Taper Point 6mm x 3mm x 1.5mm Teflon Pledget			PS3651	2/0 3,0	75	6	5 Green 5 White
			PS3653	2/0 3,0	90	12	Green
1/2 Round Body, Taper Point 6mm x 3mm x 3mm Teflon Pledget			PS2706	4/0 1,5	75	12	Green

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic


Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm						
1/2 Round Body, Taper Point		PS3727	2/0 3,0	75	12	Green
		PS4191	0 3,5	75	12	Green
1/2 Round Body, Taper Point		PS3258	3/0 2,0	75	12	Green
		PS3262	3/0 2,0	90	12	Green
		PS3730	2/0 3,0	75	6	5 Green 5 White
		PS3733	2/0 3,0	75	12	Green
		PS3737	2/0 3,0	90	6	5 Green 5 White
		PS3741	2/0 3,0	90	12	Green

25 mm							
1/2 Round Body, Taper Point			PS3738	2/0 3,0	90	12	Green
3mm x 3mm x 1.5mm Teflon Pledget							
1/2 Round Body, Taper Point			PS3263	3/0 2,0	90	12	Green
			PS3729	2/0 3,0	75	6	5 Green 5 White
			PS3732	2/0 3,0	75	12	Green
6mm x 3mm x 1.5mm Teflon Pledget							
1/2 Tapercutting			PS5826	2/0 3,0	90	6	5 Green 5 White
3mm x 3mm x 1.5mm Teflon Pledget							

26 mm						
1/2 Round Body, Taper Point		PS3788	2/0 3,0	75	12	Green
		PS4217	0 3,5	75	12	Green
1/2 Round Body, Taper Point		PS3322	3/0 2,0	90	12	Green
		PS4221	0 3,5	90	12	Green




Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm							
1/2 Round Body, Taper Point			PS3370	3/0 2,0	75	12	Green
			PS3860	2/0 3,0	75	12	Green
			PS4247	0 3,5	75	12	Green
			PS4611	1 4,0	75	12	Green
1/2 Reverse Cutting			PS4268	0 3,5	75	12	Green
			PS4623	1 4,0	75	12	Green
3/8 Reverse Cutting			PS3400	3/0 2,0	75	12	Green
			PS3892	2/0 3,0	75	12	Green
35 mm							
3/8 Reverse Cutting			PS3935	2/0 3,0	75	12	Green
36 mm							
1/2 Round Body, Taper Point			PS5515	1 4,0	75	12	Green
1/2 Reverse Cutting			PS5514	0 3,5	75	12	Green
40 mm							
1/2 Round Body, Taper Point			PS4369	0 3,5	75	12	Green
			PS4714	1 4,0	75	12	Green
			PS5007	2 5,0	90	12	Green
45 mm							
1/2 Round Body, Taper Point			PS5036	2 5,0	75	12	Green

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	PS4814	1 4,0	75	12	Green
	PS5067	2 5,0	75	12	Green
Straight Reverse Cutting 	PS4049	2/0 3,0	75	12	Green
55 mm					
1/2 Round Body, Taper Point 	PS5181	5 7,0	75	12	Green
1/2 Taper cutting 	PS5173	5 7,0	75	12	Green

Polisil


Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
17x45 cm					
POLISIL Non-Needle Suture	PS2965	3/0 2,0	17x45	12	Green
	PS3559	2/0 3,0	17x45	12	Green
	PS4150	0 3,5	17x45	12	Green
	PS4560	1 4,0	17x45	12	Green
13x60 cm					
POLISIL Non-Needle Suture	PS2962	3/0 2,0	13x60	12	Green
	PS3554	2/0 3,0	13x60	12	Green
250 cm					
POLISIL Non-Needle Suture	PS4543	1 4,0	250	12	Green

Surgical Sutures Non-Absorbable

Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament



Boz®

Definition

Synthetic Non-Absorbable

Material

(PTFE) Polytetrafluoroethylene

Coating

No Coating

Structure

Monofilament

Colour

White

Size

EP : 0.5 to 2

USP : 7/0 to 3/0

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

General Soft Tissue Approximation and/or Ligation





Cardiovascular Surgery

Dental and General Surgeries as

Well as Repair of the Dura Mater




Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/2 Round Body, Taper Point 280 Micron			MT5770	6/0 0,7	45	12	White
8.5 mm							
3/8 Reverse Cutting Premium 280 Micron			MT5768	7/0 0,5	55	12	White
10 mm							
3/8 Reverse Cutting Premium 330 Micron			MT5767	6/0 0,7	55	12	White
12 mm							
3/8 Reverse Cutting 330 Micron			MT1540	6/0 0,7	45	12	White
12.3 mm							
3/8 Reverse Cutting Premium 380 Micron			MT5763	5/0 1,0	55	12	White
3/8 Reverse Cutting Premium 430 Micron			MT5765	4/0 1,5	75	12	White
13 mm							
3/8 Round Body, Taper Point 380 Micron			MT5852	5/0 1,0	75	12	White
3/8 Round Body, Taper Point 430 Micron			MT5854	4/0 1,5	75	12	White
3/8 CC Needle 380 Micron			MT5853	5/0 1,0	75	12	White
3/8 CC Needle 430 Micron			MT5855	4/0 1,5	75	12	White
15 mm							
3/8 Round Body, Taper Point			MT5856	3/0 2,0	75	12	White








Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
15 mm							
3/8 Reverse Cutting			MT5769	6/0 0,7	45	12	White
3/8 CC Needle			MT5857	3/0 2,0	75	12	White
16 mm							
3/8 Reverse Cutting			MT5771	5/0 1,0	45	12	White
			MT5772	4/0 1,5	45	12	White
3/8 Inside Cutting			MT5755	5/0 1,0	45	12	White
16.2 mm							
3/8 Reverse Cutting Premium 430 Micron			MT5766	4/0 1,5	75	12	White
17 mm							
1/2 Round Body, Taper Point			MT5860	5/0 1,0	75	12	White
1/2 CC Needle			MT5861	5/0 1,0	75	12	White
18 mm							
3/8 Reverse Cutting			MT5829	3/0 2,0	45	12	White
			MT5558	3/0 2,0	75	12	White
18.7 mm							
3/8 Reverse Cutting 560 Micron			MT5764	3/0 2,0	75	12	White
19 mm							
3/8 Reverse Cutting			MT5773	4/0 1,5	45	12	White

Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament


Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm						
3/8 Reverse Cutting	 	MT5566	2/0 3,0	75	12	White
25 mm						
3/8 Round Body, Taper Point	 	MT5858	2/0 3,0	75	12	White
3/8 CC Needle	 	MT5859	2/0 3,0	75	12	White
30 mm						
1/2 Round Body, Taper Point	 	MT5862	0 3,5	75	12	White
1/2 CC Needle	 	MT5863	0 3,5	75	12	White

Polytef

Surgical Suture, Braided Polyester, PTFE Coated Nonabsorbable, Synthetic

Surgical Sutures
Non-Absorbable

Boz®



Definition

Synthetic Non-Absorbable

Material

Polyester Polyethylene Terephthalate (PET)

Coating

PTFE (polytetrafluoroethylene)

Structure

Braided

Colour

Green

Size

EP : 0.5 to 3

USP : 7/0 to 2/0

Tensile Strength Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

General Soft tissue approximation or ligation

Cardiovascular, Vascular Surgery


Ophthalmic Surgery

Neurosurgical Surgery



Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
10 mm							
3/8 Reverse Cutting Micropoint 280 Micron			PT5557	7/0 0,5	55	12	Green
3/8 Reverse Cutting Micropoint 330 Micron			PT5556	6/0 0,7	55	12	Green

12 mm							
3/8 Reverse Cutting Micropoint 380 Micron			PT5555	5/0 1,0	55	12	Green

16 mm							
3/8 Reverse Cutting Micropoint			PT5554	5/0 1,0	75	12	Green
3/8 Reverse Cutting Micropoint			PT5553	4/0 1,5	75	12	Green

17 mm							
1/2 Trocar			PT5552	4/0 1,5	75	12	Green

18 mm							
3/8 Reverse Cutting			PT5551	3/0 2,0	75	12	Green

20 mm							
1/2 Trocar Point			PT5550	3/0 2,0	75	12	Green
3/8 Reverse Cutting			PT5565	2/0 3,0	75	12	Green

Monowire

Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament



Boz®

Definition

Non-Absorbable

Material

316L Stainless Steel conforming to
ISO 5832 Part-1

Coating

FEP or PE

Structure

Monofilament or Multifilament

Colour

Metallic, Blue, Transperant, Orange

Size

EP : 0.7 to 10

USP : 6/0 to 8

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

Sternum Closure




Orthopaedics

Abdominal wound closure





Hernia repair

Monowire

Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament


Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm							
1/2 Reverse Cutting			MW3881	2/0 3,0	75	12	Mono
40 mm							
1/2 Reverse Cutting			MW4743	1 4,0	4 X 45	12	Mono
1/2 Tapercutting			MW4707	1 4,0	75	12	Mono
48 mm							
1/2 Reverse Cutting			MW5056	2 5,0	4 X 45	12	Mono
			MW5057	2 5,0	45	12	Mono
1/2 Reverse Cutting Heavy			MW5126	3 6,0	75	12	Mono
			MW5153	4 6,0	75	12	Mono
1/2 Tapercutting Heavy			MW5163	5 7,0	45	12	Mono
			MW5198	6 8,0	45	12	Mono
1/2 Tapercutting Heavy, Rotating			MW5146	4 6,0	45	12	Mono
55 mm							
1/2 Reverse Cutting Heavy			MW5131	3 6,0	75	12	Mono
1/2 Reverse Cutting Rotating			MW5087	2 5,0	75	12	Mono
1/2 Reverse Cutting Heavy, Rotating			MW5156	4 6,0	2 X 75	12	Mono
			MW5187	5 7,0	45	12	Mono

Monowire Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
55 mm					
1/2 Reverse Cutting Heavy, Rotating 	MW5188	5 7,0	2 X 45	12	Mono
	MW5189	5 7,0	4 X 45	12	Mono
	MW5190	5 7,0	75	12	Mono
	MW5191	5 7,0	2 X 75	12	Mono
	MW5192	5 7,0	4 X 75	12	Mono
	MW5206	6 8,0	4 X 75	12	Mono
	MW5205	6 8,0	75	12	Mono
1/2 Reverse Cutting Heavy, Rotating 	MW5193	5 7,0	75	12	Mono
1/2 Trocac Point Heavy, Rotating 	MW5777	5 7,0	4 X 45	12	Mono
	MW5179	5 7,0	75	12	Mono
	MW5180	5 7,0	2 X 75	12	Mono
	MW5203	6 8,0	75	12	Mono
	MW5204	6 8,0	2 X 75	12	Mono
	MW5208	7 9,0	75	12	Mono
	MW5209	7 9,0	2 X 75	12	Mono
57 mm					
1/2 Reverse Cutting Heavy, Rotating 	MW5195	5 7,0	75	12	Mono
	MW5196	5 7,0	2 X 75	12	Mono

Tempowire

Surgical Stainless Steel Wire Nonabsorbable, Multifilament



Boz[®]

Definition

Non-Absorbable

Material

316L Stainless Steel conforming to
(ISO 5832 Part-1)

Coating

FEP, Polyethylene

Structure

Braided, Twisted

Colour

Blue, Transparent, Orange

Size

EP : 2 to 3.5

USP : 3/0 to 0

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method










Ethylene Oxide Gas












Main Indication

Cardiac Surgery

Tempowire

Surgical Stainless Steel Wire Nonabsorbable, Multifilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
17 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle	TW3068	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5918	3/0 2,0	60	12	PE-Blue
22 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Zig - Zag type	TW3194	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5874	3/0 2,0	60	12	PE-Blue
25 mm					
1/2 Round Body, Taper Point  60 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5869	2/0 3,0	60	12	PE-Blue
26 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle	TW3327	3/0 2,0	60	12	PE-Blue

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point  85 mm Straight Reverse Cutting   	TW5877	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW3326	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW4208	0 3,5	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW5576	3/0 2,0	60	12	PE-Blue
1/2 Tapercutting  63 mm Straight Reverse Cutting  	TW3296	3/0 2,0	60	12	PE-Blue
63 mm Straight Reverse Cutting  	TW4204	0 3,5	60	12	FEP
1/2 Tapercutting  63 mm Straight Reverse Cutting   	TW3297	3/0 2,0	60	12	PE-Blue
3/8 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW5907	3/0 2,0	60	12	PE-Blue

PTFE pledget

(Accessory for sutures)



Surgical Pledget PTFE

Description

Boz® Pledgets are precut pieces from a sheet of polytetra-fluoroethylene (PTFE) non-woven fabric that is used to cushioning the load of the suture against host tissue. Pledgets are indicated as sterile, single use, implantable devices.

Indications

To buttress sutures in cardiovascular surgery.

Absorbability-Tissue Reaction

Non-Absorbable. It's inert to body fluids and enzymes.



web : www.boztibbi.com
email: boz@boztibbi.com

Teflon Pledgets Firm, White, Non-Woven, Non-Absorbable

Product No	Dimensions	Unit / Pcs
TPFE01	 6x3x1.5mm	6 pledgets per 1 pack
TPFE02	 6x3x1.5mm	5 pledgets per 1 pack
TPFE03	 6x3x1.5mm	10 pledgets per 1 pack
TPFE04	 7x3x1.5mm	6 pledgets per 1 pack
TPFE10	 3x3x1.5mm	10 pledgets per 1 pack
TPFE11	 3x3x1.5mm	5 pledgets per 1 pack
TPFE20	 6x4x1.5mm	10 pledgets per 1 pack
TPFE28	 6x8x1.5mm	4 pledgets per 1 pack
TPFE30	 6x10x1.5mm	5 pledgets per 1 pack
TPFE50	 6x6x1.5mm	10 pledgets per 1 pack
TPFE60	 7x3x1.5mm	10 pledgets per 1 pack

****Soft and firm consistencies in several geometries and sizes.
and in any required size shape**

ABSORBABLE HAEMOSTATS



Boz Tibbi Malz.San. Tic. A.Ş.

Driven by Surgical Need,
Committed to Innovative
Solution

REOXCEL

Oxidized Regenerated Cellulose (ORC)

REOXCEL - Sigma Knit

Oxidized Regenerated Cellulose (ORC)

REOXCEL - Fibril

Oxidized Regenerated Cellulose (ORC)

ZEOCCEL

Oxidized Cellulose (OC)

www.boztibbi.com
boz@boztibbi.com

REOXCEL

Absorbable Haemostat,
Oxidized Regenerated Cellulose (ORC)

Reliable Performance

Safe Haemostasis
Effective Against Bacteria



- ▶ Light knit formation covers uneven tissue
- ▶ Applicable on wide range of operations

DESCRIPTION

REOXCEL Absorbable Haemostat is a sterile haemostatic preparation available in the flat form with textile structure and made of Oxidized Regenerated Cellulose (ORC) (polyanhydro glucuronic acid). REOXCEL Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Regenerated Cellulose". REOXCEL Absorbable Haemostat does not contain any animal or collagen additives. REOXCEL Absorbable Haemostat can be sutured and cut without fraying. Its structure is stable and should be stored at controlled room temperature. The fabric is pale yellow and has a faint, caramel-like aroma. Although a slight discoloration may occur with time, this does not affect its performance.

ACTIONS

REOXCEL Absorbable Haemostat can be laid on, held against, closed exactly or sutured to a bleeding surface. REOXCEL Absorbable Haemostat is applied as dry. Action mechanism of REOXCEL Absorbable Haemostat is independent from blood coagulation mechanism of the body. REOXCEL Absorbable Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes when it contacts with blood. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. REOXCEL Absorbable Haemostat is absorbed from the sites of implantation without tissue reaction when it is used properly in minimal amounts. Absorption depends on; tissue type, haemostat amount used and saturation degree of blood. It is completely absorbed in 7-14 days. REOXCEL Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr (-) (gram negative) microorganisms including those of.

Product Code	Description of Product	Unit of Sale
RX0502	1,25cm x 5cm	12 Sachets per box
RX1001	2,5cm x 2,5cm	12 Sachets per box
RX1002	2,5cm x 5cm	12 Sachets per box
RX1003	2,5cm x 7,5cm	12 Sachets per box
RX1063	2,5cm x 9cm	12 Sachets per box
RX2003	5cm x 7,5cm	12 Sachets per box
RX2004	5cm x 10cm	12 Sachets per box
RX2014	5cm x 35cm	12 Sachets per box
RX3004	7,5cm x 10cm	12 Sachets per box
RX4004	10cm x 10cm	12 Sachets per box
RX4008	10cm x 20cm	12 Sachets per box
RX6009	15cm x 23cm	12 Sachets per box

Product Code	Description of Product	Unit of Sale
RX0502X	1,25cm x 5cm	10 Sachets per box
RX1001X	2,5cm x 2,5cm	10 Sachets per box
RX1002X	2,5cm x 5cm	10 Sachets per box
RX2003X	5cm x 7,5cm	10 Sachets per box
RX2004X	5cm x 10cm	10 Sachets per box
RX2014X	5cm x 35cm	10 Sachets per box
RX3004X	7,5cm x 10cm	10 Sachets per box
RX4004X	10cm x 10cm	10 Sachets per box
RX4008X	10cm x 20cm	10 Sachets per box
RX6009X	15cm x 23cm	10 Sachets per box

REOXCEL - Sigma Knit

Absorbable Haemostat,
Oxidized Regenerated Cellulose (ORC)

**Unmatched
Stability Against
Heavier Bleeding**
Safe Haemostasis
Effective Against Bacteria



- Dense knit web, sustain surgical suture use
- Ideal web structure where wrapping is required

DESCRIPTION

REOXCEL-Sigma Knit Absorbable Haemostat is a sterile haemostatic preparation available in the flat form with textile structure and made of Oxidized Regenerated Cellulose (ORC) (polyanhydro glucuronic acid). REOXCEL-Sigma Knit Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Regenerated Cellulose". REOXCEL-Sigma Knit Absorbable Haemostat does not contain any animal or collagen additives. REOXCEL-Sigma Knit Absorbable Haemostat can be sutured and cut without fraying. Its structure is stable and should be stored at controlled room temperature. The fabric is pale yellow and has a faint, caramel-like aroma. Although a slight discoloration may occur with time, this does not affect its performance.

ACTIONS

REOXCEL-Sigma Knit Absorbable Haemostat can be laid on, held against, closed exactly or sutured to a bleeding surface. REOXCEL-Sigma Knit Absorbable Haemostat is applied as dry. Action mechanism of REOXCEL-Sigma Knit Absorbable Haemostat is independent from blood coagulation mechanism of the body. REOXCEL-Sigma Knit Absorbable Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes when it contacts with blood. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. REOXCEL-Sigma Knit Absorbable Haemostat is absorbed from the sites of implantation without tissue reaction when it is used properly in minimal amounts. Absorption depends on; tissue type, haemostat amount used and saturation degree of blood. It is completely absorbed in 7-14 days. REOXCEL-Sigma Knit Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr (-) (gram negative) microorganisms including those of.

Product code	Description of Product	Unit of Sale
RK1001	2,6cm x 2,6cm	12 Sachets per box
RK1003	2,6cm x 10cm	12 Sachets per box
RK2003	5cm x 7,5cm	12 Sachets per box
RK3004	7,6cm x 10,2cm	12 Sachets per box
RK6009	15,2cm x 23cm	10 Sachets per box

Product code	Description of Product	Unit of Sale
RK1001X	2,6cm x 2,6cm	10 Sachets per box
RK1003X	2,6cm x 10cm	10 Sachets per box
RK2003X	5cm x 7,5cm	10 Sachets per box
RK3004X	7,6cm x 10,2cm	10 Sachets per box

REOXCEL - Fibril

Absorbable Haemostat,
Oxidized Regenerated Cellulose (ORC)

All-Round Performance

Safe Haemostasis
Effective Against Bacteria



- ▶ Thin layers of web adapt easily to uneven and hard to
- ▶ Moistened light layer web allows bipolar cauter use
- ▶ Range of diverse forms; Light single layers, Roll, Pad, Tufts

DESCRIPTION

REOXCEL FIBRIL Absorbable Haemostat is a sterile haemostatic preparation available in the fibrous form and made of Oxidized Regenerated Cellulose (ORC) (polyanhydro glucuronic acid). REOXCEL FIBRIL Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Regenerated Cellulose". REOXCEL FIBRIL Absorbable Haemostat does not contain any animal or collagen additives. REOXCEL FIBRIL Absorbable Haemostat can be cut without fraying. Its structure is stable and should be stored at controlled room temperature. It is pale yellow and has a faint, caramel-like aroma. Although a slight discoloration may occur with time, this does not affect its performance. Under favour of its fibrous form surgeon can grasp with forceps any amount of REOXCEL FIBRIL Absorbable Haemostat needed to achieve haemostasis at a particular bleeding site. This provides surgeon to reach irregularly shaped bleeding sites. It is easy to pull the desired amount of REOXCEL FIBRIL Absorbable Haemostat and the group of selected fibres continue to cohere to one another. It does not disperse over the operative site unwantedly.

ACTIONS


REOXCEL FIBRIL Absorbable Haemostat can be laid on, held against, closed exactly to a bleeding surface. REOXCEL FIBRIL Absorbable Haemostat is applied as dry. Action mechanism of REOXCEL FIBRIL Absorbable Haemostat is independent from blood coagulation mechanism of the body. REOXCEL FIBRIL Absorbable Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes when it contacts with blood. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. REOXCEL FIBRIL Absorbable Haemostat is absorbed from the sites of implantation without tissue reaction when it is used properly in minimal amounts. Absorption depends on; tissue type, haemostat amount used and saturation degree of blood. It is completely absorbed in 7-14 days. REOXCEL FIBRIL Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr(-) (gram negative) microorganisms including those of.

Product Code	Description of Product	Unit of Sale
RF1002	2.6cm x 5.1cm	10 Sachets per box
RF2003	5cm x 7.5cm	10 Sachets per box
RF2004	5.1cm x 10.2cm	10 Sachets per box
RF4004	10.2cm x 10.2cm	10 Sachets per box

ZEOCEL

Absorbable Haemostat,
Oxidized Cellulose (OC)

**Safe Haemostasis
Effective Against
Bacteria**



DESCRIPTION

ZEOCEL Absorbable Haemostat is a woven sterile hemostatic agent formed of oxidized cellulose (OC) (polyanhydro glucuronic acid). It has been manufactured from alpha grade cotton. Agent is a pale yellow material and has a caramel like smell. As it is stored in controlled room temperature, there might be a slight change in color in time but this does not effect the performance of the material. ZEOCEL Absorbable Haemostat complies with the requirements of the United States Pharmacopoeia for "Oxidized Cellulose". ZEOCEL Absorbable Haemostat does not contain any animal or collagen additives.

ACTIONS

ZEOCEL Absorbable Haemostat can be applied to bleeding surface in different forms. The haemostatic agent material can be closed exactly, pressured against or sutured to the desired usage area. ZEOCEL Absorbable Haemostat should be used as dry. After the application of the material to the bleeding surface, haemostatic agent will form a black gelatin like mass. The gelatin like mass will aid haemostasis process by allowing the platelet aggregation and formation of platelet-fibrin plug. The coagulation will be achieved in 3-4 minutes after ZEOCEL's contact with the bleeding area. When ZEOCEL Absorbable Haemostat is used only as needed, it will be absorbed from the site of usage without tissue reaction. Absorption of the haemostatic agent from the used site depends on the amount of haemostatic agent used, tissue type and blood saturation degree. The material will be fully absorbed in 7-14 days.

ZEOCEL Absorbable Haemostat has been shown to be bactericidal in vitro against 40 types of Gr(+) (gram positive) and Gr(-) (gram negative) microorganisms including those of.

Product code	Description of Product	Unit of Sale
ZC0502	1.25cm x 5cm	12 Sachets per box
ZC1001	2.5cm x 2.5cm	12 Sachets per box
ZC1002	2.5cm x 5cm	12 Sachets per box
ZC2003	5cm x 7.5cm	12 Sachets per box
ZC2004	5cm x 7cm	12 Sachets per box
ZC2014	5cm x 35cm	12 Sachets per box
ZC3004	7.5cm x 10cm	12 Sachets per box
ZC4004	10.2cm x 10.2cm	12 Sachets per box
ZC4008	10cm x 20cm	12 Sachets per box
ZC6009	15cm x 23cm	12 Sachets per box

Product code	Description of Product	Unit of Sale
ZC0502X	1.25cm x 5cm	10 Sachets per box
ZC1001X	2.5cm x 2.5cm	10 Sachets per box
ZC1002X	2.5cm x 5cm	10 Sachets per box
ZC2003X	5cm x 7.5cm	10 Sachets per box
ZC2004X	5cm x 7cm	10 Sachets per box
ZC2014X	5cm x 35cm	10 Sachets per box
ZC3004X	7.5cm x 10cm	10 Sachets per box
ZC4004X	10.2cm x 10.2cm	10 Sachets per box
ZC4008X	10cm x 20cm	10 Sachets per box
ZC6009X	15cm x 23cm	10 Sachets per box

Creative

ideas

suggestions

Boz®

Boz Tıbbi Malz.San. Tic. A.Ş.

Specialist Products



Microsurgery
Products



Retention
Padded Plates



Polypropylene
Button




Polypropylene
Locking Button



Polypropylene Bead
Aluminium Ring

www.boztibbi.com
boz@boztibbi.com

Surgical Indications



Eye, Hand, Plastic, Brain and Nerve operations.

Needle



Hardened stainless steel for custom purposes, no bending, no breaking, microscopically inspected available in taperpoint, spatulated, lancet.

Thread

Excellent resistance, easy tying, safe knot, non-capillary, extremely smooth surface, just the right elasticity. Virgin Silk, Nylon, Polypropylene, Polyvinylidene fluoride, Polyglycolic acid.


Packaging



Keeps suture material free of tangles. Protects the sharpness of needle. Provides double armed suture for single arming by locating the mid-point. Stores the second suture until needed. Provides a disposable pad for needle count.

Description

For the treatment and prevention of abdominal ruptures.




Description

Polypropylene button, White,
Ø = 12 mm, thickness.....mm, Ø hole =.....mm.



Indications




Polypropylene button in the Boz tendon repair kit is intended to be used in Bunnell Technique in tendon surgery to act as knot-fastening support over the skin to avoid suture cut and to offer pull-out suture opportunity after healing.

Absorbability-Tissue Reaction

Nonabsorbable, skin contact.

Description

Polypropylene locking button.



Indications



The Boz polypropylene locking button is intended to be used as anchor device instead of knot in subcuticular suturing .

Absorbability-Tissue Reaction

Nonabsorbable, skin surface application.

Description




Polypropylene Bead, White.
 $\varnothing = 6 \text{ mm}$, $\varnothing \text{ hole} = 1.7 \text{ mm}$,



Aluminium Ring
 $\varnothing = 3 \text{ mm}$, $h = 3 \text{ mm}$, $\varnothing \text{ hole} = 1.7 \text{ mm}$,

Indications



Absorbability-Tissue Reaction

Nonabsorbable, skin contact.



Boz Medical Device Industry and Trade Inc.

Sağlık Mahallesi Sağlık Sokak No:33/5 Sıhhiye/Ankara/TÜRKİYE

Tel: + 90 (312) 254 03 40 Fax: +90 (312) 254 03 49

web: www.boztibbi.com - email: boz@boztibbi.com

GenelKatalog-rev-00-26-10-2017



Product
Catalog Link



www.boztibbi.com/bozkatalog.pdf