

Surgical Sutures Absorbable

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

BOZ[®]



Definition

Synthetic absorbable, coated, mid term wound support and mid term mass absorption.

Material

Polyglycolic Acid (PGA)

Coating

Calcium Stearate and polycaprolactone

Structure

Braided - Multifilament

Colour

Violet (dyed) or Beige (undyed)

Size

EP : 0.4 to 8

USP : 8/0 to 6

Tensile Strength Retention

2. Week 75 % - 3. Week 50 %

Mass Absorption

60-90 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery

Gynaecology

Ophthalmic Surgery

Orthopedics

Urology

Skin closer - intracutaneous, subcutaneous.

Ligatures

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
6 mm						
3/8 Spatula 220 Micron	 	GS5401	8/0 0,4	45	12	Violet
3/8 Spatula 230 Micron	 	GS1156	8/0 0,4	45	12	Violet
		GS1252	7/0 0,5	30	12	Violet
6,6 mm						
1/2 Cobra 230 Micron	 	GS5519	8/0 0,4	60	12	Violet
1/2 Cobra 280 Micron	 	GS5518	7/0 0,5	60	12	Violet
8 mm						
1/2 Spatula 230 Micron	 	GS1193	8/0 0,4	30	12	Violet
		GS1195	8/0 0,4	45	12	Violet
1/2 Spatula 330 Micron	 	GS1272	7/0 0,5	45	12	Violet
		GS1395	6/0 0,7	45	12	Violet
1/2 Spatula 340 Micron	 	GS5526	5/0 1,0	45	12	Violet
1/4 Spatula 230 Micron	 	GS1190	8/0 0,4	45	12	Violet
1/4 Spatula 250 Micron	 	GS5876	7/0 0,5	45	12	Violet
1/4 Spatula 340 Micron	 	GS1268	7/0 0,5	45	12	Violet
		GS1404	6/0 0,7	30	12	Violet
		GS1405	6/0 0,7	45	12	Violet
		GS1406	6/0 0,7	45	12	Undyed
		GS1764	5/0 1,0	45	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
8 mm						
1/4 Spatula 430 Micron	 	GS2362	4/0 1,5	45	12	Violet
1/4 Spatula 430 Micron	 	GS2363	4/0 1,5	45	12	Violet
3/8 Round Body, Taper Point 230 Micron	 	GS1207	8/0 0,4	30	12	Violet
3/8 Round Body, Taper Point 280 Micron	 	GS1413	6/0 0,7	45	12	Violet
		GS1426	6/0 0,7	45	12	Violet
		GS1768	5/0 1,0	45	12	Violet
10 mm						
3/8 Round Body, Taper Point 280 Micron	 	GS1480	6/0 0,7	75	12	Violet
3/8 Reverse Cutting 330 Micron	 	GS1496	6/0 0,7	45	12	Violet
3/8 Reverse Cutting 330 Micron	 	GS1492	6/0 0,7	45	12	Violet
3/8 Inside Cutting 280 Micron	 	GS1508	6/0 0,7	45	12	Violet
12 mm						
3/8 Reverse Cutting 330 Micron	 	GS1541	6/0 0,7	45	12	Violet
3/8 Spatula 280 Micron	 	GS1523	6/0 0,7	45	12	Violet
13 mm						
1/2 Round Body, Taper Point 330 Micron	 	GS1560	6/0 0,7	45	12	Violet
		GS1564	6/0 0,7	75	12	Violet

















Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
13 mm						
1/2 Round Body, Taper Point 330 Micron	 	GS1561	6/0 0,7	45	12	Violet
1/2 Round Body, Taper Point 380 Micron	 	GS1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron	 	GS1894	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 330 Micron	 	GS1580	6/0 0,7	45	12	Violet
		GS1590	6/0 0,7	75	12	Violet
3/8 Round Body, Taper Point 380 Micron	 	GS1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron	 	GS2432	4/0 1,5	45	12	Violet
		GS2434	4/0 1,5	60	12	Violet
3/8 Reverse Cutting 380 Micron	 	GS1928	5/0 1,0	75	12	Violet
3/8 Reverse Cutting 560 Micron	 	GS2456	4/0 1,5	75	12	Violet
3/8 Spatula 280 Micron	 	GS1605	6/0 0,7	45	12	Violet
3/8 Spatula 330 Micron	 	GS1873	5/0 1,0	60	12	Violet
3/8 Trocar 330 Micron	 	GS2416	4/0 1,5	45	12	Violet
16 mm						
1/2 Round Body, Taper Point	 	GS1990	5/0 1,0	75	12	Violet
		GS1992	5/0 1,0	75	12	Undyed
		GS2503	4/0 1,5	75	12	Undyed
		GS2504	4/0 1,5	75	12	Violet






Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			GS3005	3/0 2,0	75	12	Violet
			GS1674	6/0 0,7	75	12	Violet
3/8 Reverse Cutting			GS2014	5/0 1,0	45	12	Undyed
			GS2540	4/0 1,5	75	12	Violet
			GS3025	3/0 2,0	75	12	Violet
			GS2140	5/0 1,0	45	12	Violet
3/8 Inside Cutting			GS2144	5/0 1,0	75	12	Violet
			GS2557	4/0 1,5	75	12	Undyed
			GS3039	3/0 2,0	75	12	Violet
			GS2140	5/0 1,0	45	12	Violet
17 mm							
1/2 Round Body, Taper Point			GS2584	4/0 1,5	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GS5525	5/0 1,0	75	12	Violet
			GS2613	4/0 1,5	75	12	Violet
			GS3624	2/0 3,0	75	12	Violet
3/8 Round Body, Taper Point			GS3085	3/0 2,0	75	12	Violet
19 mm							
3/8 Reverse Cutting			GS1706	6/0 0,7	75	12	Violet
			GS2640	4/0 1,5	75	12	Violet
3/8 Reverse Cutting Premium			GS2211	5/0 1,0	45	12	Violet



Glikosorb









Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided



Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm					
1/2 Round Body, Taper Point 	GS2229	5/0 1,0	75	12	Violet
	GS2685	4/0 1,5	75	12	Violet
	GS3142	3/0 2,0	75	12	Violet
1/2 Tapercutting 	GS2673	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 	GS2250	5/0 1,0	75	12	Undyed
	GS2716	4/0 1,5	45	12	Undyed
	GS2719	4/0 1,5	75	12	Violet
	GS3164	3/0 2,0	45	12	Violet
	GS3168	3/0 2,0	75	12	Violet
	GS3663	2/0 3,0	75	12	Violet
22 mm					
1/2 Round Body, Taper Point 	GS2747	4/0 1,5	75	12	Violet
	GS3196	3/0 2,0	75	12	Violet
	GS3671	2/0 3,0	75	12	Violet
	GS4575	1 4,0	75	12	Violet
3/8 Reverse Cutting 	GS2763	4/0 1,5	75	12	Violet
	GS3214	3/0 2,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided















Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
24 mm					
1/2 Reverse Cutting  	GS5524	4/0 1,5	75	12	Violet
	GS5523	3/0 2,0	75	12	Violet
	GS3692	2/0 3,0	75	12	Violet

25 mm					
1/2 Round Body, Taper Point  	GS2788	4/0 1,5	75	12	Violet
	GS3252	3/0 2,0	75	12	Violet
	GS3746	2/0 3,0	100	12	Violet
	GS3724	2/0 3,0	75	12	Violet
1/2 Tapercutting  	GS3718	2/0 3,0	75	12	Violet
3/8 Reverse Cutting  	GS2804	4/0 1,5	75	12	Violet
	GS3757	2/0 3,0	75	12	Violet
3/8 Inside Cutting  	GS2812	4/0 1,5	45	12	Violet
	GS3284	3/0 2,0	45	12	Violet
	GS3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point  	GS2292	5/0 1,0	75	12	Violet
	GS2828	4/0 1,5	75	12	Violet
	GS3313	3/0 2,0	75	12	Violet
	GS3318	3/0 2,0	90	12	Violet
	GS3786	2/0 3,0	75	12	Violet







Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point  	GS4215	0 3,5	75	12	Violet
	GS4592	1 4,0	90	12	Violet
1/2 Reverse Cutting  	GS4230	0 3,5	75	12	Violet
3/8 Round Body, Taper Point  	GS2841	4/0 1,5	75	12	Violet
3/8 Reverse Cutting  	GS3345	3/0 2,0	75	12	Undyed
	GS3346	3/0 2,0	75	12	Violet
	GS3821	2/0 3,0	75	12	Undyed
	GS3822	2/0 3,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point  	GS3369	3/0 2,0	75	12	Violet
	GS3374	3/0 2,0	90	12	Violet
	GS3858	2/0 3,0	75	12	Violet
	GS4244	0 3,5	75	12	Violet
	GS4608	1 4,0	75	12	Violet
1/2 Reverse Cutting  	GS3882	2/0 3,0	75	12	Violet
	GS4266	0 3,5	75	12	Violet
1/2 Tapercutting  	GS3843	2/0 3,0	75	12	Violet





Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
3/8 Reverse Cutting 	GS2887	4/0 1,5	75	12	Violet
	GS3397	3/0 2,0	75	12	Violet
	GS3891	2/0 3,0	75	12	Violet
	GS4273	0 3,5	75	12	Violet
32 mm					
1/2 Round Body, Taper Point 	GS5508	3/0 2,0	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GS3411	3/0 2,0	75	12	Violet
	GS3916	2/0 3,0	75	12	Violet
	GS3922	2/0 3,0	90	12	Violet
	GS4300	0 3,5	75	12	Violet
	GS4644	1 4,0	75	12	Violet
3/8 Reverse Cutting 	GS3423	3/0 2,0	75	12	Violet
	GS3933	2/0 3,0	75	12	Violet
36 mm					
1/2 Round Body, Taper Point 	GS5509	3/0 2,0	75	12	Violet
1/2 Reverse Cutting 	GS4675	1 4,0	75	12	Violet









Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
37 mm					
1/2 Round Body, Taper Point 	GS3950	2/0 3,0	75	12	Violet
	GS4337	0 3,5	75	12	Violet
	GS4681	1 4,0	75	12	Violet
1/2 Reverse Cutting 	GS5507	2/0 3,0	75	12	Violet
	GS4354	0 3,5	90	12	Violet
	GS4994	2 5,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	GS3982	2/0 3,0	75	12	Violet
	GS4368	0 3,5	75	12	Violet
	GS4375	0 3,5	90	12	Violet
	GS4713	1 4,0	75	12	Violet
	GS4719	1 4,0	90	12	Violet
	GS5003	2 5,0	75	12	Violet
	GS5006	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GS4397	0 3,5	75	12	Violet
	GS4745	1 4,0	75	12	Violet
	GS4749	1 4,0	90	12	Violet
	GS5019	2 5,0	75	12	Violet






Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Reverse Cutting 	GS5023	2 5,0	90	12	Violet
1/2 Tapercutting 	GS4363	0 3,5	75	12	Violet
3/8 Reverse Cutting 	GS4756	1 4,0	75	12	Violet
	GS4408	0 3,5	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GS4422	0 3,5	75	12	Violet
	GS4772	1 4,0	75	12	Violet
	GS5035	2 5,0	75	12	Violet
	GS5040	2 5,0	90	12	Violet
1/2 Tapercutting 	GS4766	1 4,0	75	12	Violet
	GS5738	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GS4432	0 3,5	75	12	Violet
	GS4787	1 4,0	75	12	Violet
Straight Reverse Cutting 	GS2909	4/0 1,5	45	12	Violet
48 mm					
1/2 Round Body, Taper Point 	GS4796	1 4,0	75	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GS4454	0 3,5	75	12	Violet
	GS4459	0 3,5	90	12	Violet
	GS4812	1 4,0	75	12	Violet
	GS4818	1 4,0	90	12	Violet
	GS5066	2 5,0	75	12	Violet
	GS5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GS4471	0 3,5	75	12	Violet
60 mm					
Straight Reverse Cutting 	GS3492	3/0 2,0	75	12	Violet
	GS4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GS4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GS4500	0 3,5	90	12	Violet
	GS4861	1 4,0	90	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOSORB Non-Needle Suture	GS2950	3/0 2,0	6x45	12	Violet
	GS3543	2/0 3,0	6x45	12	Violet
	GS4137	0 3,5	6x45	12	Violet
	GS4563	1 4,0	6x45	12	Violet
10x50 cm					
GLIKOSORB Non-Needle Suture	GS2954	3/0 2,0	10x50	12	Violet
	GS3548	2/0 3,0	10x75	12	Violet
12x45 cm					
GLIKOSORB Non-Needle Suture	GS2344	4/0 1,5	12x45	12	Violet
	GS2958	3/0 2,0	12x45	12	Violet
	GS3552	2/0 3,0	12x45	12	Violet
	GS4145	0 3,5	12x45	12	Violet
	GS4554	1 4,0	12x45	12	Violet
17x45 cm					
GLIKOSORB Non-Needle Suture	GS2966	3/0 2,0	17x45	12	Violet
	GS3558	2/0 3,0	17x45	12	Violet
150 cm					
GLIKOSORB Non-Needle Suture	GS1737	5/0 1,0	150	12	Violet
	GS2326	4/0 1,5	150	12	Violet
	GS2935	3/0 2,0	150	12	Violet

Glikosorb

Surgical Suture, Polyglycolicacid (PGA) Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
GLIKOSORB Non-Needle Suture	GS3511	2/0 3,0	150	12	Violet
	GS4114	0 3,5	150	12	Violet
	GS4532	1 4,0	150	12	Violet

250 cm					
GLIKOSORB Non-Needle Suture	GS2332	4/0 1,5	250	12	Violet
	GS2942	3/0 2,0	250	12	Violet
	GS3518	2/0 3,0	250	12	Violet
	GS4121	0 3,5	250	12	Violet
	GS4541	1 4,0	250	12	Violet