

GLIKOLAK

Absorbable Surgical Suture | Emilebilir Cerrahi Sütür



PGLA

Poly [glycolide(90%)-co-lactide(10%)]



Material

Poly [glycolide(90%)-co-lactide(10%)], Synthetic

Colour

Violet Purple and Undyed

Structure

Braided-Multifilament

Coating

Calcium Stearate, %30 Glycolide, %70 Lactide

Tensile Strength

2nd week ~ %75 wound support

Absorbtion

~ 55 - 70 days

USP - EP Range

USP : 8/0 - 2 EP : 0.4 - 5

Sterilization

EO (Ethylene Oxide)

About of Glikolak Surgical Sutures

Glikolak sutures are synthetic, stripped, absorbable, sterile surgical sutures made from polyglycolide and lactide. Stripped sutures have special coating material. Coating of the sutures were made with Calcium Stearate and Poly (Glycolide Lactide)(30/70) mixture.

Glikolak sutures have violet colour to make surgical operations easier. Natural colouring agents were used in colouring process. Glikolak sutures are manufactured with or without stainless steel needle, at 8/0 range, and with various lengths between 30cm to 150cm.

Main Indication



Gastrointestinal Surgery



Gynaecology



Ophthalmic Surgery



Orthopedics



Urology

































Skin Closure



Ligatures





Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
6 mm							
3/8 Reverse Cutting 150 Micron			GL1165	8/0 0,4	45	12	Violet
3/8 Reverse Cutting 150 Micron			GL1256	7/0 0,5	45	12	Violet
3/8 Spatula 230 Micron			GL1155	8/0 0,4	30	12	Violet
3/8 Spatula 230 Micron			GL1253	7/0 0,5	45	12	Violet
8 mm							
1/2 Spatula 330 Micron			GL1394	6/0 0,7	30	12	Violet
1/4 Spatula 340 Micron			GL1405	6/0 0,7	45	12	Violet
1/4 Spatula 340 Micron			GL1404	6/0 0,7	30	12	Violet
			GL1764	5/0 1,0	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1299	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1432	6/0 0,7	75	12	Violet
3/8 Spatula 230 Micron			GL5827	8/0 0,4	45	12	Violet
3/8 Spatula 280 Micron			GL1413	6/0 0,7	45	12	Violet
10 mm							
3/8 Round Body, Taper Point 340 Micron			GL1803	5/0 1,0	45	12	Violet
3/8 Reverse Cutting 330 Micron			GL1496	6/0 0,7	45	12	Violet
3/8 Spatula 280 Micron			GL1450	6/0 0,7	45	12	Violet
3/8 Spatula 330 Micron			GL1790	5/0 1,0	75	12	Violet













Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
12 mm						
1/2 Round Body, Taper Point 430 Micron	 	GL1836	5/0 1,0	45	12	Violet
13 mm						
1/2 Round Body, Taper Point 330 Micron	 	GL1560	6/0 0,7	45	12	Violet
		GL1564	6/0 0,7	75	12	Violet
1/2 Round Body, Taper Point 430 Micron	 	GL1894	5/0 1,0	75	12	Violet
		GL1900	5/0 1,0	90	12	Violet
3/8 Round Body, Taper Point 330 Micron	 	GL1580	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 380 Micron	 	GL1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron	 	GL2434	4/0 1,5	60	12	Violet
3/8 Reverse Cutting Premium 330 Micron	 	GL1626	6/0 0,7	45	12	Violet
16 mm						
1/2 Round Body, Taper Point	 	GL1990	5/0 1,0	75	12	Violet
		GL1994	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point	 	GL2504	4/0 1,5	75	12	Violet
		GL3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting	 	GL2021	5/0 1,0	75	12	Violet
		GL2540	4/0 1,5	75	12	Violet
		GL3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting	 	GL2144	5/0 1,0	75	12	Violet


Glikolak



Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided



Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm					
3/8 Inside Cutting  	GL2557	4/0 1,5	75	12	Undyed
	GL2553	4/0 1,5	45	12	Violet
	GL3039	3/0 2,0	75	12	Violet
17 mm					
1/2 Round Body, Taper Point  	GL2165	5/0 1,0	75	12	Violet
18 mm					
1/2 Round Body, Taper Point  	GL2613	4/0 1,5	75	12	Violet
19 mm					
3/8 Reverse Cutting  	GL1706	6/0 0,7	75	12	Violet
	GL2206	5/0 1,0	75	12	Violet
	GL2640	4/0 1,5	75	12	Violet
	GL3109	3/0 2,0	75	12	Violet
20 mm					
1/2 Round Body, Taper Point  	GL2229	5/0 1,0	75	12	Violet
	GL2685	4/0 1,5	75	12	Violet
	GL2690	4/0 1,5	90	12	Violet
	GL3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting  	GL2716	4/0 1,5	45	12	Undyed
	GL3164	3/0 2,0	45	12	Violet
	GL3663	2/0 3,0	75	12	Violet


Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	GL2788	4/0 1,5	75	12	Violet
	GL3252	3/0 2,0	75	12	Violet
	GL3746	2/0 3,0	100	12	Violet
	GL3724	2/0 3,0	75	12	Violet






25 mm					
3/8 Reverse Cutting 	GL2804	4/0 1,5	75	12	Violet
	GL2812	4/0 1,5	45	12	Violet
3/8 Inside Cutting 	GL3284	3/0 2,0	45	12	Violet
	GL3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point 	GL3313	3/0 2,0	75	12	Violet
	GL4220	0 3,5	90	12	Violet
3/8 Reverse Cutting 	GL2852	4/0 1,5	75	12	Violet
	GL3346	3/0 2,0	75	12	Violet
	GL3822	2/0 3,0	75	12	Violet

30 mm					
1/2 Round Body, Taper Point 	GL3369	3/0 2,0	75	12	Violet
	GL3858	2/0 3,0	75	12	Violet
	GL3865	2/0 3,0	90	12	Violet







Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	GL4244	0 3,5	75	12	Violet
	GL4250	0 3,5	90	12	Violet
	GL4608	1 4,0	75	12	Violet
	GL4614	1 4,0	90	12	Violet
1/2 Round Body, Taper Point 	GL3854	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GL3397	3/0 2,0	75	12	Violet
	GL3891	2/0 3,0	75	12	Violet
	GL3895	2/0 3,0	90	12	Violet
	GL4273	0 3,5	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GL3916	2/0 3,0	75	12	Violet
	GL4300	0 3,5	75	12	Violet
	GL4644	1 4,0	75	12	Violet
	GL4649	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GL3423	3/0 2,0	75	12	Violet
	GL3933	2/0 3,0	75	12	Violet
	GL5881	1 4,0	75	12	Violet






Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	GL3982	2/0 3,0	75	12	Violet
	GL4368	0 3,5	75	12	Violet
	GL4375	0 3,5	90	12	Violet
	GL4713	1 4,0	75	12	Violet
	GL4719	1 4,0	90	12	Violet
	GL5003	2 5,0	75	12	Violet
1/2 Reverse Cutting 	GL4399	0 3,5	90	12	Violet
	GL4749	1 4,0	90	12	Violet
	GL5019	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4408	0 3,5	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GL4772	1 4,0	75	12	Violet
	GL4775	1 4,0	90	12	Violet
	GL5035	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4022	2/0 3,0	75	12	Violet
	GL4432	0 3,5	75	12	Violet
	GL4787	1 4,0	75	12	Violet
Straight Reverse Cutting 	GL2909	4/0 1,5	45	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GL4454	0 3,5	75	12	Violet
	GL4459	0 3,5	90	12	Violet
	GL4812	1 4,0	75	12	Violet
	GL4818	1 4,0	90	12	Violet
	GL5066	2 5,0	75	12	Violet
	GL5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GL5078	2 5,0	75	12	Violet
60 mm					
Straight Reverse Cutting 	GL2923	4/0 1,5	75	12	Violet
	GL3492	3/0 2,0	75	12	Violet
	GL4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GL4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GL4861	1 4,0	90	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOLAK Non-Needle Suture	GL2950	3/0 2,0	6x45	12	Violet
	GL4137	0 3,5	6x45	12	Violet
	GL4563	1 4,0	6x45	12	Violet

12x45 cm					
GLIKOLAK Non-Needle Suture	GL2344	4/0 1,5	12x45	12	Violet
	GL2958	3/0 2,0	12x45	12	Violet
	GL3552	2/0 3,0	12x45	12	Violet
	GL4145	0 3,5	12x45	12	Violet
	GL4554	1 4,0	12x45	12	Violet

17x45 cm					
GLIKOLAK Non-Needle Suture	GL2966	3/0 2,0	17x45	12	Violet
	GL4559	1 4,0	17x45	12	Violet

10x50 cm					
GLIKOLAK Non-Needle Suture	GL4140	0 3,5	10x50	12	Violet
	GL2955	3/0 2,0	10x75	12	Violet

10x75 cm					
GLIKOLAK Non-Needle Suture	GL4550	1 4,0	10x75	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
GLIKOLAK Non-Needle Suture	GL1737	5/0 1,0	150	12	Violet
	GL2326	4/0 1,5	150	12	Violet
	GL2935	3/0 2,0	150	12	Violet
	GL3511	2/0 3,0	150	12	Violet
	GL4114	0 3,5	150	12	Violet
	GL4532	1 4,0	150	12	Violet