

REOXCEL®

ABSORBABLE HAEMOSTAT, OXIDIZED REGENERATED CELLULOSE (ORC)

 **Bactericidal Effect Against 40 Kind of Microorganisms**

 **Clot Formation in 3-4 Minutes**

 **Completely Absorbed in 7-14 Days**


 **Absorption in implantation areas without tissue reaction**


Reoxcel absorbable hemostat is a sterile hemostatic prepare with plain textile structure made from oxidised regenerated cellulose (ORC) (polyanhydro glucuronic acid). Reoxcel absorbable hemostat is compliant with USP American Pharmacopeia "Oxidised Regenerated Cellulose" standard requirements.


Reoxcel absorbable hemostat has no animal or chalcogen additive substances. This hemostat can be sutured and cut without any wear. It has durable structure and should be stored at controlled room temperature. The textile has light yellow colour and slight caramel-like odour.


Wide Usage Area


Reoxcel absorbable hemostat is designed to stop bleeding in capillary bleeding, parenchymatous organ bleeding, and bleeding in reaction areas during surgical operations. This hemostat is used when other hemostatic procedures are non-applicable. It has wide usage area. It can be used in important surgical operations.


General Surgical



Digestive Surgery



Neurosurgery Operations



Plastic Surgery



Orthopedy Operations



Gynecology Operations



Urology Operations



Stomatology Operations



Throat Surgeries



Nose surgeries



Traumatology Surgical


Cardiovascular Surgery


Face and Jaw Surgery


Lung Operations


Haemorrhoidectomy Operations


Biopsi Operations


These hemostats are suitable for all other surgical fields including general surgery, intestine system surgery, neurosurgery (especially brain operations), plastic surgery, orthopaedics, gynaecology, urology, stomatology, otorhinolaryngology operations, traumatology and cardiovascular operations, vascular prosthetics operations, haemorrhoidectomy, biopsy, liver and gall bladder operations, gastric resections, thoracic and abdominal sympathectomy. Reoxcel absorbable hemostat can be applied to cavities (after tumour scraping) as well as endoscopic operations and dental applications.

Product Features

Product Code

Description of Product

Unit of Sale

RX0502

1,25cm x 5cm

12 Sachets per box

RX1001

2,5cm x 2,5cm

12 Sachets per box

RX1002

2,5cm x 5cm

12 Sachets per box

RX1003

2,5cm x 7,5cm

12 Sachets per box

RX1063

2,5cm x 9cm

12 Sachets per box

RX2003

5cm x 7,5cm

12 Sachets per box

RX2004

5cm x 10cm

12 Sachets per box

RX2014

5cm x 35cm

12 Sachets per box

RX3004

7,5cm x 10cm

12 Sachets per box

RX4004

10cm x 10cm

12 Sachets per box

RX4008

10cm x 20cm

12 Sachets per box

RX6009

15cm x 23cm

12 Sachets per box

RX0502X

1,25cm x 5cm

10 Sachets per box

RX1001X

2,5cm x 2,5cm

10 Sachets per box

RX1002X

2,5cm x 5cm

10 Sachets per box

RX2003X

5cm x 7,5cm

10 Sachets per box

RX2004X

5cm x 10cm

10 Sachets per box

RX2014X

5cm x 35cm

10 Sachets per box

RX3004X

7,5cm x 10cm

10 Sachets per box

RX4004X

10cm x 10cm

10 Sachets per box

RX4008X

10cm x 20cm

10 Sachets per box

RX6009X

15cm x 23cm

10 Sachets per box



Reliable Performance

Safe Haemostasis

Effective Against Bacteria

- ▶ Light knit formation covers uneven tissue
- ▶ Applicable on wide range of operations

