

ZEOCEL®

ABSORBABLE HAEMOSTAT, OXIDIZED CELLULOSE (OC)

✓ Bactericidal effects against 40 types of microorganism

⌚ Clot Formation in 3-4 Minutes

📅 Completely Absorbed in 7-14 Days

➡ Absorption in implantation areas without tissue reaction


Zeocel absorbable hemostat is a weaved, sterile hemostatic agent with Oxidised Cellulose (OC) (polyanhydro glucuronic acid) form. This hemostat is manufactured from alpha-graded cotton. Agent is light yellow materials and has caramel-like odour. Since this hemostat is stored at controlled room temperature, the colour may change over time. But such colour change will have no impact on material performance.

Zeocel absorbable hemostat is compliant with USP American Pharmacopeia “Oxidised Cellulose” standard requirements. This hemostat has no animal or chalcogen additive substances. Zeocel can be applied to bleeding surfaces in various ways. Material can be wrapped around, pressed, or sutured.


Zeocel Usage Areas

Zeocel absorbable hemostat is designed to stop other different type capillary bleeding and control bleeding leakage based on tissue resection in any type of surgical operation. When other hemostatic procedures are unsuitable, Zeocel absorbable hemostat can be used.


General Surgical




Digestive Surgery




Plastic Surgery




Orthopedy Operations




Gynecology Operations




Urology Operations



Stomatology Operations



Traumatology Surgical



Zeocel absorbable hemostat can be used in various surgery types including general surgery, digestion system surgery, plastic surgery, orthopaedic, gynaecology, urology, stomatology, traumatology and many other surgical branches. Hemostats can be used for filling the cavities due to tumour removal in addition to endoscopic operations and dental applications.

Product Features



Product Code	Description of Product	Unit of Sale
ZC0502	1,25cm x 5cm	12 Sachets per box
ZC1001	2,5cm x 2,5cm	12 Sachets per box
ZC1002	2,5cm x 5cm	12 Sachets per box
ZC2003	5cm x 7,5cm	12 Sachets per box
ZC2004	5cm x 7cm	12 Sachets per box
ZC2014	5cm x 35cm	12 Sachets per box
ZC3004	7,5cm x 10cm	12 Sachets per box
ZC4004	10cm x 10cm	12 Sachets per box
ZC4008	10cm x 20cm	12 Sachets per box
ZC6009	15cm x 23cm	12 Sachets per box
ZC0502X	1,25cm x 5cm	10 Sachets per box
ZC1001X	2,5cm x 2,5cm	10 Sachets per box
ZC1002X	2,5cm x 5cm	10 Sachets per box
ZC2003X	5cm x 7,5cm	10 Sachets per box
ZC2004X	5cm x 7cm	10 Sachets per box
ZC2014X	5cm x 35cm	10 Sachets per box
ZC3004X	7,5cm x 10cm	10 Sachets per box
ZC4004X	10cm x 10cm	10 Sachets per box
ZC4008X	10cm x 20cm	10 Sachets per box
ZC6009X	15cm x 23cm	10 Sachets per box

Absorbable **Haemostat**

Safe Haemostasis
Effective Against Bacteria

