

Chromic Catgut

Absorbable Surgical Suture



Definition

Collagen based biomaterial, absorbable.

Material

Collagen tissue derived from purified sheep intestine

Coating

No coating

Structure

Monofilament

Colour

Brown

Size

EP : 1 to 8

USP : 6/0 to 4

Sterilization Method

Ethylene Oxide Gas

Gamma Irradiation


Main Indication

Ligatures

Ophthalmic Surgery

Chromic Catgut


Absorbable Surgical Suture

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 430 Micron With Alchol			CK1892A	5/0 1,5	75	12	Brown
3/8 Reverse Cutting 430 Micron			CK1933	5/0 1,5	75	12	Brown
16 mm							
1/2 Round Body, Taper Point			CK1989	5/0 1,5	75	12	Brown
			CK2503	4/0 2,0	75	12	Brown
			CK3004	3/0 3,0	75	12	Brown
1/2 Round Body, Taper Point With Alchol			CK1989A	5/0 1,5	75	12	Brown
			CK2503A	4/0 2,0	75	12	Brown
			CK3004A	3/0 3,0	75	12	Brown
3/8 Reverse Cutting			CK2019	5/0 1,5	75	12	Brown
			CK2539	4/0 2,0	75	12	Brown
			CK3024	3/0 3,0	75	12	Brown
3/8 Reverse Cutting With Alchol			CK1670A	6/0 1,0	75	12	Brown
			CK2539A	4/0 2,0	75	12	Brown
			CK3024A	3/0 3,0	75	12	Brown
19 mm							
3/8 Reverse Cutting			CK2639	4/0 2,0	75	12	Brown
3/8 Reverse Cutting With Alchol			CK2205A	5/0 1,5	75	12	Brown
			CK2639A	4/0 2,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut Absorbable Surgical Suture


Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			CK2228	5/0 1,5	75	12	Brown
			CK2683	4/0 2,0	75	12	Brown
			CK3141	3/0 3,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol			CK2683A	4/0 2,0	75	12	Brown
			CK3141A	3/0 3,0	75	12	Brown
3/8 Reverse Cutting			CK2250	5/0 1,5	75	12	Brown
			CK2718	4/0 2,0	75	12	Brown
			CK3167	3/0 3,0	75	12	Brown
3/8 Reverse Cutting With Alcohol			CK2250A	5/0 1,5	75	12	Brown
			CK2718A	4/0 2,0	75	12	Brown
			CK3167A	3/0 3,0	75	12	Brown
25 mm							
3/8 Reverse Cutting			CK2800	4/0 2,0	75	12	Brown
3/8 Reverse Cutting With Alcohol			CK2800A	4/0 2,0	75	12	Brown
26 mm							
1/2 Round Body, Taper Point			CK2827	4/0 2,0	75	12	Brown
			CK3312	3/0 3,0	75	12	Brown
			CK3785	2/0 3,5	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point  	CK4214	0 4,0	75	12	Brown
	CK4589	1 5,0	75	12	Brown
	CK4946	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol  	CK2827A	4/0 2,0	75	12	Brown
	CK3312A	3/0 3,0	75	12	Brown
	CK3785A	2/0 3,5	75	12	Brown
	CK4214A	0 4,0	75	12	Brown
	CK4589A	1 5,0	75	12	Brown
3/8 Reverse Cutting  	CK2851	4/0 2,0	75	12	Brown
	CK3345	3/0 3,0	75	12	Brown
	CK3821	2/0 3,5	75	12	Brown
	CK4236	0 4,0	75	12	Brown
3/8 Reverse Cutting With Alcohol  	CK2851A	4/0 2,0	75	12	Brown
	CK3345A	3/0 3,0	75	12	Brown
	CK3821A	2/0 3,5	75	12	Brown
	CK4236A	0 4,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	CK3368	3/0 3,0	75	12	Brown
	CK3857	2/0 3,5	75	12	Brown
	CK4243	0 4,0	75	12	Brown
	CK4954	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK3857A	2/0 3,5	75	12	Brown
	CK4243A	0 4,0	75	12	Brown
	CK4607A	1 5,0	75	12	Brown
	CK4954A	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK3396	3/0 3,0	75	12	Brown
	CK3890	2/0 3,5	75	12	Brown
	CK4628	1 5,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK3396A	3/0 3,0	75	12	Brown
	CK3890A	2/0 3,5	75	12	Brown
	CK4272A	0 4,0	75	12	Brown
35 mm					
1/2 Round Body, Taper Point 	CK3410	3/0 3,0	75	12	Brown
	CK3915	2/0 3,5	75	12	Brown
	CK4299	0 4,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotates the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
1/2 Round Body, Taper Point 	CK4643	1 5,0	75	12	Brown
	CK4969	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4299A	0 4,0	75	12	Brown
	CK4643A	1 5,0	75	12	Brown
1/2 Reverse Cutting 	CK4976	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK3422	3/0 3,0	75	12	Brown
	CK4319	0 4,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4319A	0 4,0	75	12	Brown
40 mm					
1/2 Round Body, Taper Point 	CK4712	1 5,0	75	12	Brown
	CK5001	2 6,0	75	12	Brown
	CK5118	3 7,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4367A	0 4,0	75	12	Brown
	CK5001A	2 6,0	75	12	Brown
	CK5118A	3 7,0	75	12	Brown
1/2 Reverse Cutting 	CK5018	2 6,0	75	12	Brown
	CK5120	3 7,0	75	12	Brown
1/2 Reverse Cutting With Alcohol 	CK5120A	3 7,0	75	12	Brown



*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.


*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
3/8 Reverse Cutting 	CK4407	0 4,0	75	12	Brown
	CK4755	1 5,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4407A	0 4,0	75	12	Brown
	CK4755A	1 5,0	75	12	Brown

45 mm					
1/2 Round Body, Taper Point 	CK4421	0 4,0	75	12	Brown
	CK5034	2 6,0	75	12	Brown
1/2 Round Body, Taper Point With Alcohol 	CK4421A	0 4,0	75	12	Brown
	CK4771A	1 5,0	75	12	Brown
	CK5034A	2 6,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4431A	0 4,0	75	12	Brown
Straight Reverse Cutting 	CK2911	4/0 2,0	75	12	Brown
Straight Reverse Cutting With Alcohol 	CK2911A	4/0 2,0	75	12	Brown





50 mm					
1/2 Round Body, Taper Point 	CK4811	1 5,0	75	12	Brown
	CK5065	2 6,0	75	12	Brown
	CK5127	3 7,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point With Alcohol 	CK4811A	1 5,0	75	12	Brown
	CK5065A	2 6,0	75	12	Brown
	CK5127A	3 7,0	75	12	Brown
1/2 Reverse Cutting 	CK5077	2 6,0	75	12	Brown
3/8 Reverse Cutting 	CK4474	0 4,0	75	12	Brown
	CK4840	1 5,0	75	12	Brown
	CK5082	2 6,0	75	12	Brown
	CK5130	3 7,0	75	12	Brown
3/8 Reverse Cutting With Alcohol 	CK4840A	1 5,0	75	12	Brown
	CK5082A	2 6,0	75	12	Brown
	CK5130A	3 7,0	75	12	Brown
60 mm					
Straight Reverse Cutting 	CK3491	3/0 3,0	75	12	Brown
	CK4073	2/0 3,5	75	12	Brown
Straight Reverse Cutting With Alcohol 	CK4073A	2/0 3,5	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denetos the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
CHROMIC CATGUT Non-Needle Suture	CK1373	6/0 1,0	150	12	Brown
	CK1736	5/0 1,5	150	12	Brown
	CK2325	4/0 2,0	150	12	Brown
	CK2934	3/0 3,0	150	12	Brown
	CK3510	2/0 3,5	150	12	Brown
	CK4113	0 4,0	150	12	Brown
	CK4531	1 5,0	150	12	Brown
	CK4907	2 6,0	150	12	Brown
	CK5106	3 7,0	150	12	Brown
	CK5134	4 8,0	150	12	Brown

150 cm With Alcohol					
CHROMIC CATGUT Non-Needle Suture With Alcohol	CK1373A	6/0 1,0	150	12	Brown
	CK1736A	5/0 1,5	150	12	Brown
	CK2325A	4/0 2,0	150	12	Brown
	CK2934A	3/0 3,0	150	12	Brown
	CK3510A	2/0 3,5	150	12	Brown
	CK4113A	0 4,0	150	12	Brown
	CK4531A	1 5,0	150	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denetos the product is "wet" e.g.

*After code number CK####A

Chromic Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm With Alchol					
CHROMIC CATGUT Non-Needle Suture With Alchol	CK4907A	2 6,0	150	12	Brown
	CK5106A	3 7,0	150	12	Brown
	CK5134A	4 8,0	150	12	Brown

75 cm					
CHROMIC CATGUT Non-Needle Suture	CK1748	5/0 1,5	75	12	Brown
	CK2316	4/0 2,0	75	12	Brown
	CK2925	3/0 3,0	75	12	Brown
	CK3530	2/0 3,5	75	12	Brown
	CK4104	0 4,0	75	12	Brown
	CK4521	1 5,0	75	12	Brown
	CK4898	2 6,0	75	12	Brown

75 cm With Alchol					
CHROMIC CATGUT Non-Needle Suture With Alchol	CK1748A	5/0 1,5	75	12	Brown
	CK2316A	4/0 2,0	75	12	Brown
	CK2925A	3/0 3,0	75	12	Brown
	CK3530A	2/0 3,5	75	12	Brown
	CK4104A	0 4,0	75	12	Brown
	CK4521A	1 5,0	75	12	Brown
	CK4898A	2 6,0	75	12	Brown

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denetos the product is "wet" e.g.

*After code number CK####A