


Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)

Fast Absorbable, Synthetic, Braided, Undyed

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Coating

Poly(glycolide-co-lactide) (30/70) & calcium stearate

Structure

Braided

Colour

Beige (undyed)

Size

EP : 0.4 to 5

USP : 8/0 to 2

Tensile Strength Retention

7 Days 50 % - 10-14 Days 0 %

Mass Absorption

42 days

Sterilization Method


Gamma Irradiation

Main Indication

Gynaecology (episiotomies)
Urology
Ophtalmic surgery (conjunctive suturing)
Odontology (oral mucosa)
Pediatric surgery
Skin closer - intracutaneous, subcutaneous.
Ligatures

















Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
3/8 Reverse Cutting 330 Micron			FL1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 380 Micron			FL1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FL1926	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FL1989	5/0 1,0	75	12	Undyed
			FL2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FL1673	6/0 0,7	75	12	Undyed
			FL2019	5/0 1,0	75	12	Undyed
			FL2539	4/0 1,5	75	12	Undyed
			FL3024	3/0 2,0	75	12	Undyed
3/8 Inside Cutting			FL2557	4/0 1,5	75	12	Undyed
			FL3038	3/0 2,0	75	12	Undyed
17 mm							
1/2 Round Body, Taper Point			FL1687	6/0 0,7	75	12	Undyed
18 mm							
1/2 Round Body, Taper Point			FL3080	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FL2639	4/0 1,5	75	12	Undyed
			FL3108	3/0 2,0	75	12	Undyed













Fastlak






Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm						
1/2 Round Body, Taper Point  		FL2228	5/0 1,0	75	12	Undyed
		FL2683	4/0 1,5	75	12	Undyed
		FL3141	3/0 2,0	75	12	Undyed
		FL3150	3/0 2,0	90	12	Undyed
3/8 Reverse Cutting  		FL2250	5/0 1,0	75	12	Undyed
22 mm						
1/2 Round Body, Taper Point  		FL2745	4/0 1,5	75	12	Undyed
		FL3195	3/0 2,0	75	12	Undyed
1/2 Inside Cutting  		FL3216	3/0 2,0	45	12	Undyed
3/8 Reverse Cutting  		FL3210	3/0 2,0	75	12	Undyed
25 mm						
1/2 Round Body, Taper Point  		FL2787	4/0 1,5	75	12	Undyed
		FL3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting  		FL2805	4/0 1,5	75	12	Undyed
26 mm						
1/2 Round Body, Taper Point  		FL3312	3/0 2,0	75	12	Undyed
		FL3785	2/0 3,0	75	12	Undyed
3/8 Reverse Cutting  		FL3345	3/0 2,0	75	12	Undyed
		FL3821	2/0 3,0	75	12	Undyed

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
30 mm							
1/2 Round Body, Taper Point			FL3368	3/0 2,0	75	12	Undyed
			FL3857	2/0 3,0	75	12	Undyed
			FL4243	0 3,5	75	12	Undyed
			FL4607	1 4,0	75	12	Undyed
3/8 Reverse Cutting			FL2885	4/0 1,5	75	12	Undyed
			FL3396	3/0 2,0	75	12	Undyed
			FL3890	2/0 3,0	75	12	Undyed
35 mm							
1/2 Round Body, Taper Point			FL3915	2/0 3,0	75	12	Undyed
			FL4299	0 3,5	75	12	Undyed
			FL4305	0 3,5	90	12	Undyed
			FL4643	1 4,0	75	12	Undyed
3/8 Reverse Cutting			FL3422	3/0 2,0	75	12	Undyed
			FL3932	2/0 3,0	75	12	Undyed
36 mm							
1/2 Round Body, Taper Point			FL4667	1 4,0	75	12	Undyed
40 mm							
1/2 Round Body, Taper Point			FL3981	2/0 3,0	75	12	Undyed
			FL3988	2/0 3,0	90	12	Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	FL4367	0 3,5	75	12	Undyed
	FL4374	0 3,5	90	12	Undyed
	FL4712	1 4,0	75	12	Undyed
	FL4718	1 4,0	90	12	Undyed
	FL5001	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	FL4005	2/0 3,0	75	12	Undyed
	FL4755	1 4,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	FL4424	0 3,5	90	12	Undyed
	FL4771	1 4,0	75	12	Undyed
	FL5034	2 5,0	75	12	Undyed
Straight Reverse Cutting 	FL2911	4/0 1,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	FL4037	2/0 3,0	75	12	Undyed
	FL4453	0 3,5	75	12	Undyed
	FL4811	1 4,0	75	12	Undyed
	FL5065	2 5,0	75	12	Undyed

Fastlak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
FASTLAK Non-Needle Suture	FL1736	5/0 1,0	150	12	Undyed
	FL2325	4/0 1,5	150	12	Undyed
	FL2934	3/0 2,0	150	12	Undyed
	FL3510	2/0 3,0	150	12	Undyed
	FL4113	0 3,5	150	12	Undyed
	FL4531	1 4,0	150	12	Undyed
	FL4907	2 5,0	150	12	Undyed