

# Fastsorb

Surgical Suture, Polyglycolicacid (PGA)  
Fast Absorbable, Synthetic, Braided, Undyed

**Surgical Sutures**  
**Absorbable**



**Boz®**

## Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

## Material

Polyglycolic Acid (PGA)

## Coating

Calcium Stearate and polycaprolactone

## Structure

Braided

## Colour

Beige (undyed)

## Size

EP : 0.4 to 8

USP : 8/0 to 6

## Tensile Strength Retention

7 Days 75 %

## Mass Absorption

42 days

## Sterilization Method

Gamma Irradiation

## Main Indication

Gynaecology (episiotomies)

Urology

Ophtalmic surgery ( conjunctive suturing)

Odontology (oral mucosa)






















Pediatric surgery

Skin closer - intracutaneous, subcutaneous.

Ligatures

















# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
3/8 Reverse Cutting 330 Micron			FS1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 330 Micron			FS1563	6/0 0,7	75	12	Undyed
1/2 Round Body, Taper Point 380 Micron			FS1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 330 Micron			FS5912	6/0 0,7	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FS1929	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting Premium 380 Micron			FS1941	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FS1989	5/0 1,0	75	12	Undyed
			FS2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FS1673	6/0 0,7	75	12	Undyed
			FS2019	5/0 1,0	75	12	Undyed
			FS2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			FS2557	4/0 1,5	75	12	Undyed
			FS3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FS2639	4/0 1,5	75	12	Undyed
			FS3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			FS2228	5/0 1,0	75	12	Undyed






# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			FS2683	4/0 1,5	75	12	Undyed
			FS3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS2250	5/0 1,0	75	12	Undyed
22 mm							
1/2 Round Body, Taper Point			FS2745	4/0 1,5	75	12	Undyed
25 mm							
1/2 Round Body, Taper Point			FS2787	4/0 1,5	75	12	Undyed
			FS3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS3273	3/0 2,0	75	12	Undyed
			FS3756	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body, Taper Point			FS3785	2/0 3,0	75	12	Undyed
			FS4214	0 3,5	75	12	Undyed
3/8 Reverse Cutting			FS3821	2/0 3,0	75	12	Undyed
			FS4236	0 3,5	75	12	Undyed
			FS3345	3/0 2,0	75	12	Undyed
30 mm							
1/2 Round Body, Taper Point			FS3368	3/0 2,0	75	12	Undyed
			FS3857	2/0 3,0	75	12	Undyed
			FS4243	0 3,5	75	12	Undyed
			FS4607	1 4,0	75	12	Undyed







# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting		FS2885	4/0 1,5	75	12	Undyed
		FS3396	3/0 2,0	75	12	Undyed
		FS4628	1 4,0	75	12	Undyed
35 mm						
1/2 Round Body, Taper Point		FS3915	2/0 3,0	75	12	Undyed
		FS4299	0 3,5	75	12	Undyed
		FS4305	0 3,5	90	12	Undyed
		FS4643	1 4,0	75	12	Undyed
		FS4648	1 4,0	90	12	Undyed
3/8 Reverse Cutting		FS3422	3/0 2,0	75	12	Undyed
		FS3932	2/0 3,0	75	12	Undyed
37 mm						
1/2 Round Body, Taper Point		FS3954	2/0 3,0	90	12	Undyed
40 mm						
1/2 Round Body, Taper Point		FS3981	2/0 3,0	75	12	Undyed
		FS3988	2/0 3,0	90	12	Undyed
		FS4367	0 3,5	75	12	Undyed
		FS4374	0 3,5	90	12	Undyed

# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>40 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4712</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS4718</b>	<b>1</b> 4,0	90	12	Undyed
	<b>FS5001</b>	<b>2</b> 5,0	75	12	Undyed
<b>3/8</b> Reverse Cutting 	<b>FS4005</b>	<b>2/0</b> 3,0	75	12	Undyed
<b>45 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4421</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4771</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS5034</b>	<b>2</b> 5,0	75	12	Undyed
<b>Straight</b> Reverse Cutting 	<b>FS2911</b>	<b>4/0</b> 1,5	75	12	Undyed
<b>50 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4037</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>FS4453</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4811</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS5065</b>	<b>2</b> 5,0	75	12	Undyed
<b>65 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4864</b>	<b>1</b> 4,0	100	12	Undyed

# Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
<b>FASTSORB</b> Non-Needle Suture	<b>FS1373</b>	<b>6/0</b> 0,7	150	12	Undyed
	<b>FS2325</b>	<b>4/0</b> 1,5	150	12	Undyed
	<b>FS2934</b>	<b>3/0</b> 2,0	150	12	Undyed
	<b>FS3510</b>	<b>2/0</b> 3,0	150	12	Undyed
	<b>FS4113</b>	<b>0</b> 3,5	150	12	Undyed
	<b>FS4531</b>	<b>1</b> 4,0	150	12	Undyed
	<b>FS4907</b>	<b>2</b> 5,0	150	12	Undyed