

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, coated, mid term wound support and mid term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Coating

Poly(glycolide-co-lactide) (30/70) & calcium stearate

Structure

Braided-Multifilament

Colour

Beige (undyed) or Violet (dyed)

Size

EP : 0.4 to 5

USP : 8/0 to 2

Tensile Strength Retention

2. Week 75%

Mass Absorption

55-70 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery
Gynaecology
Ophthalmic Surgery
Orthopedics
Urology
Skin closer - intracutaneous, subcutaneous.
Ligatures







Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
6 mm							
3/8 Reverse Cutting 150 Micron			GL1165	8/0 0,4	45	12	Violet
3/8 Reverse Cutting 150 Micron			GL1256	7/0 0,5	45	12	Violet
3/8 Spatula 230 Micron			GL1155	8/0 0,4	30	12	Violet
3/8 Spatula 230 Micron			GL1253	7/0 0,5	45	12	Violet
8 mm							
1/2 Spatula 330 Micron			GL1394	6/0 0,7	30	12	Violet
1/4 Spatula 340 Micron			GL1405	6/0 0,7	45	12	Violet
1/4 Spatula 340 Micron			GL1404	6/0 0,7	30	12	Violet
			GL1764	5/0 1,0	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1299	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			GL1432	6/0 0,7	75	12	Violet
3/8 Spatula 230 Micron			GL5827	8/0 0,4	45	12	Violet
3/8 Spatula 280 Micron			GL1413	6/0 0,7	45	12	Violet
10 mm							
3/8 Round Body, Taper Point 340 Micron			GL1803	5/0 1,0	45	12	Violet
3/8 Reverse Cutting 330 Micron			GL1496	6/0 0,7	45	12	Violet
3/8 Spatula 280 Micron			GL1450	6/0 0,7	45	12	Violet
3/8 Spatula 330 Micron			GL1790	5/0 1,0	75	12	Violet


Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
1/2 Round Body, Taper Point 430 Micron			GL1836	5/0 1,0	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 330 Micron			GL1560	6/0 0,7	45	12	Violet
			GL1564	6/0 0,7	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			GL1894	5/0 1,0	75	12	Violet
			GL1900	5/0 1,0	90	12	Violet
3/8 Round Body, Taper Point 330 Micron			GL1580	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 380 Micron			GL1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GL2434	4/0 1,5	60	12	Violet
3/8 Reverse Cutting Premium 330 Micron			GL1626	6/0 0,7	45	12	Violet
16 mm							
1/2 Round Body, Taper Point			GL1990	5/0 1,0	75	12	Violet
			GL1994	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point			GL2504	4/0 1,5	75	12	Violet
			GL3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GL2021	5/0 1,0	75	12	Violet
			GL2540	4/0 1,5	75	12	Violet
			GL3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GL2144	5/0 1,0	75	12	Violet



Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm					
3/8 Inside Cutting  	GL2557	4/0 1,5	75	12	Undyed
	GL2553	4/0 1,5	45	12	Violet
	GL3039	3/0 2,0	75	12	Violet
17 mm					
1/2 Round Body, Taper Point  	GL2165	5/0 1,0	75	12	Violet
18 mm					
1/2 Round Body, Taper Point  	GL2613	4/0 1,5	75	12	Violet
19 mm					
3/8 Reverse Cutting  	GL1706	6/0 0,7	75	12	Violet
	GL2206	5/0 1,0	75	12	Violet
	GL2640	4/0 1,5	75	12	Violet
	GL3109	3/0 2,0	75	12	Violet
20 mm					
1/2 Round Body, Taper Point  	GL2229	5/0 1,0	75	12	Violet
	GL2685	4/0 1,5	75	12	Violet
	GL2690	4/0 1,5	90	12	Violet
	GL3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting  	GL2716	4/0 1,5	45	12	Undyed
	GL3164	3/0 2,0	45	12	Violet
	GL3663	2/0 3,0	75	12	Violet


Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color	
25 mm							
1/2 Round Body, Taper Point			GL2788	4/0 1,5	75	12	Violet
			GL3252	3/0 2,0	75	12	Violet
			GL3746	2/0 3,0	100	12	Violet
			GL3724	2/0 3,0	75	12	Violet
25 mm							
3/8 Reverse Cutting			GL2804	4/0 1,5	75	12	Violet
3/8 Inside Cutting			GL2812	4/0 1,5	45	12	Violet
			GL3284	3/0 2,0	45	12	Violet
			GL3763	2/0 3,0	75	12	Violet
26 mm							
1/2 Round Body, Taper Point			GL3313	3/0 2,0	75	12	Violet
			GL4220	0 3,5	90	12	Violet
3/8 Reverse Cutting			GL2852	4/0 1,5	75	12	Violet
			GL3346	3/0 2,0	75	12	Violet
			GL3822	2/0 3,0	75	12	Violet
30 mm							
1/2 Round Body, Taper Point			GL3369	3/0 2,0	75	12	Violet
			GL3858	2/0 3,0	75	12	Violet
			GL3865	2/0 3,0	90	12	Violet



Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	GL4244	0 3,5	75	12	Violet
	GL4250	0 3,5	90	12	Violet
	GL4608	1 4,0	75	12	Violet
	GL4614	1 4,0	90	12	Violet
1/2 Round Body, Taper Point 	GL3854	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GL3397	3/0 2,0	75	12	Violet
	GL3891	2/0 3,0	75	12	Violet
	GL3895	2/0 3,0	90	12	Violet
	GL4273	0 3,5	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GL3916	2/0 3,0	75	12	Violet
	GL4300	0 3,5	75	12	Violet
	GL4644	1 4,0	75	12	Violet
	GL4649	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GL3423	3/0 2,0	75	12	Violet
	GL3933	2/0 3,0	75	12	Violet
	GL5881	1 4,0	75	12	Violet




Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	GL3982	2/0 3,0	75	12	Violet
	GL4368	0 3,5	75	12	Violet
	GL4375	0 3,5	90	12	Violet
	GL4713	1 4,0	75	12	Violet
	GL4719	1 4,0	90	12	Violet
	GL5003	2 5,0	75	12	Violet
1/2 Reverse Cutting 	GL4399	0 3,5	90	12	Violet
	GL4749	1 4,0	90	12	Violet
	GL5019	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4408	0 3,5	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GL4772	1 4,0	75	12	Violet
	GL4775	1 4,0	90	12	Violet
	GL5035	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GL4022	2/0 3,0	75	12	Violet
	GL4432	0 3,5	75	12	Violet
	GL4787	1 4,0	75	12	Violet
Straight Reverse Cutting 	GL2909	4/0 1,5	45	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GL4454	0 3,5	75	12	Violet
	GL4459	0 3,5	90	12	Violet
	GL4812	1 4,0	75	12	Violet
	GL4818	1 4,0	90	12	Violet
	GL5066	2 5,0	75	12	Violet
	GL5070	2 5,0	90	12	Violet
1/2 Reverse Cutting 	GL5078	2 5,0	75	12	Violet
60 mm					
Straight Reverse Cutting 	GL2923	4/0 1,5	75	12	Violet
	GL3492	3/0 2,0	75	12	Violet
	GL4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GL4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GL4861	1 4,0	90	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOLAK Non-Needle Suture	GL2950	3/0 2,0	6x45	12	Violet
	GL4137	0 3,5	6x45	12	Violet
	GL4563	1 4,0	6x45	12	Violet

12x45 cm					
GLIKOLAK Non-Needle Suture	GL2344	4/0 1,5	12x45	12	Violet
	GL2958	3/0 2,0	12x45	12	Violet
	GL3552	2/0 3,0	12x45	12	Violet
	GL4145	0 3,5	12x45	12	Violet
	GL4554	1 4,0	12x45	12	Violet

17x45 cm					
GLIKOLAK Non-Needle Suture	GL2966	3/0 2,0	17x45	12	Violet
	GL4559	1 4,0	17x45	12	Violet

10x50 cm					
GLIKOLAK Non-Needle Suture	GL4140	0 3,5	10x50	12	Violet
	GL2955	3/0 2,0	10x75	12	Violet

10x75 cm					
GLIKOLAK Non-Needle Suture	GL4550	1 4,0	10x75	12	Violet

Glikolak

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
GLIKOLAK Non-Needle Suture	GL1737	5/0 1,0	150	12	Violet
	GL2326	4/0 1,5	150	12	Violet
	GL2935	3/0 2,0	150	12	Violet
	GL3511	2/0 3,0	150	12	Violet
	GL4114	0 3,5	150	12	Violet
	GL4532	1 4,0	150	12	Violet