

Antibacterial Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, coated, antibacterial mid term wound support and mid term mass absorption.

Material

Poly(glycolide-co-lactide)
glycolide 90% -lactide 10 % (PGLA)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

Poly(glycolide-co-lactide) & calcium stearate

Structure

Braided-Multifilament

Colour

Beige (undyed) or Violet (dyed)

Size

EP : 1 to 5
USP : 5/0 to 2

Tensile Strenght Retention

2. Week 75%

Mass Absorption

55-70 days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gastrointestinal Surgery
Gynaecology
Orthopedics
Urology
Skin closer - intracutaneous, subcutaneous.
Ligatures

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 340 Micron			GLP1764	5/0 1,0	45	12	Violet
3/8 Spatula 230 Micron			GLP1768	5/0 1,0	45	12	Violet
3/8 Round Body, Taper Point 380 Micron			GLP1778	5/0 1,0	75	12	Violet

13 mm							
1/2 Round Body, Taper Point 380 Micron			GLP1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			GLP1897	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			GLP1914	5/0 1,0	75	12	Violet
3/8 Round Body, Taper Point 430 Micron			GLP2434	4/0 1,5	60	12	Violet
3/8 Round Body, Taper Point 430 Micron			GLP2436	4/0 1,5	75	12	Violet
3/8 Reverse Cutting 380 Micron			GLP1928	5/0 1,0	75	12	Violet

16 mm							
1/2 Round Body, Taper Point			GLP1990	5/0 1,0	75	12	Violet
			GLP2504	4/0 1,5	75	12	Violet
			GLP3005	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GLP2021	5/0 1,0	75	12	Violet
			GLP2540	4/0 1,5	75	12	Violet
			GLP3025	3/0 2,0	75	12	Violet
3/8 Inside Cutting			GLP2144	5/0 1,0	75	12	Violet
			GLP2553	4/0 1,5	45	12	Violet

Antibacterial**Glikolak Plus**Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
3/8 Inside Cutting			GLP2557	4/0 1,5	75	12	Undyed
			GLP3039	3/0 2,0	75	12	Violet
18 mm							
1/2 Round Body, Taper Point			GLP2613	4/0 1,5	75	12	Violet
19 mm							
3/8 Reverse Cutting			GLP2206	5/0 1,0	75	12	Violet
			GLP2640	4/0 1,5	75	12	Violet
20 mm							
1/2 Round Body, Taper Point			GLP2229	5/0 1,0	75	12	Violet
			GLP2685	4/0 1,5	75	12	Violet
			GLP2690	4/0 1,5	90	12	Violet
			GLP3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			GLP2716	4/0 1,5	45	12	Undyed
			GLP3164	3/0 2,0	45	12	Violet
			GLP3168	3/0 2,0	75	12	Violet
22 mm							
3/8 Reverse Cutting			GLP2763	4/0 1,5	75	12	Violet
25 mm							
1/2 Round Body, Taper Point			GLP2788	4/0 1,5	75	12	Violet
			GLP3252	3/0 2,0	75	12	Violet

Antibacterial

Glikolak Plus


Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	GLP3746	2/0 3,0	100	12	Violet
	GLP3724	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GLP2804	4/0 1,5	75	12	Violet
3/8 Inside Cutting 	GLP2812	4/0 1,5	45	12	Violet
	GLP3284	3/0 2,0	45	12	Violet
	GLP3763	2/0 3,0	75	12	Violet

26 mm					
1/2 Round Body, Taper Point 	GLP2828	4/0 1,5	75	12	Violet
	GLP3313	3/0 2,0	75	12	Violet
	GLP3786	2/0 3,0	75	12	Violet
	GLP4220	0 3,5	90	12	Violet
3/8 Reverse Cutting 	GLP2852	4/0 1,5	75	12	Violet
	GLP3346	3/0 2,0	75	12	Violet
	GLP3822	2/0 3,0	75	12	Violet

30 mm					
1/2 Round Body, Taper Point 	GLP3369	3/0 2,0	75	12	Violet
	GLP3858	2/0 3,0	75	12	Violet
	GLP3865	2/0 3,0	90	12	Violet
	GLP4244	0 3,5	75	12	Violet




Antibacterial**Glikolak Plus**Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point 	GLP4250	0 3,5	90	12	Violet
	GLP4608	1 4,0	75	12	Violet
	GLP4614	1 4,0	90	12	Violet
1/2 Round Body, Taper Point 	GLP3854	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	GLP3397	3/0 2,0	75	12	Violet
	GLP3891	2/0 3,0	75	12	Violet
	GLP3895	2/0 3,0	90	12	Violet
	GLP4273	0 3,5	75	12	Violet
35 mm					
1/2 Round Body, Taper Point 	GLP3411	3/0 2,0	75	12	Violet
	GLP3916	2/0 3,0	75	12	Violet
	GLP4300	0 3,5	75	12	Violet
	GLP4644	1 4,0	75	12	Violet
	GLP4649	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GLP3423	3/0 2,0	75	12	Violet
	GLP3933	2/0 3,0	75	12	Violet
3/8 Inside Cutting 	GLP3938	2/0 3,0	75	12	Violet

Antibacterial

Glikolak Plus




Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	GLP3982	2/0 3,0	75	12	Violet
	GLP4368	0 3,5	75	12	Violet
	GLP4375	0 3,5	90	12	Violet
	GLP4713	1 4,0	75	12	Violet
	GLP4719	1 4,0	90	12	Violet
	GLP5003	2 5,0	75	12	Violet
1/2 Reverse Cutting 	GLP4399	0 3,5	90	12	Violet
	GLP4749	1 4,0	90	12	Violet
3/8 Reverse Cutting 	GLP4756	1 4,0	75	12	Violet
45 mm					
1/2 Round Body, Taper Point 	GLP4422	0 3,5	75	12	Violet
	GLP4772	1 4,0	75	12	Violet
	GLP4775	1 4,0	90	12	Violet
	GLP5035	2 5,0	75	12	Violet
3/8 Reverse Cutting 	GLP4787	1 4,0	75	12	Violet
50 mm					
1/2 Round Body, Taper Point 	GLP4454	0 3,5	75	12	Violet
	GLP4459	0 3,5	90	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	GLP4812	1 4,0	75	12	Violet
	GLP4818	1 4,0	90	12	Violet
	GLP5066	2 5,0	75	12	Violet
	GLP5070	2 5,0	90	12	Violet
60 mm					
Straight Reverse Cutting 	GLP2923	4/0 1,5	75	12	Violet
	GLP3492	3/0 2,0	75	12	Violet
	GLP3498	3/0 2,0	90	12	Violet
	GLP4074	2/0 3,0	75	12	Violet
63 mm					
3/8 Blunt 	GLP4085	2/0 3,0	75	12	Violet
	GLP4858	1 4,0	75	12	Violet
64 mm					
3/8 Blunt 	GLP4861	1 4,0	90	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
6x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2950	3/0 2,0	6x45	12	Violet
	GLP4137	0 3,5	6x45	12	Violet
12x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2344	4/0 1,5	12x45	12	Violet
	GLP2958	3/0 2,0	12x45	12	Violet
	GLP3552	2/0 3,0	12x45	12	Violet
	GLP4145	0 3,5	12x45	12	Violet
17x45 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2349	4/0 1,5	17x45	12	Violet
	GLP2966	3/0 2,0	17x45	12	Violet
	GLP3558	2/0 3,0	17x45	12	Violet
	GLP4149	0 3,5	17x45	12	Violet
10x50 cm					
GLIKOLAK Plus Non-Needle Suture	GLP3547	2/0 3,0	10x50	12	Violet
	GLP4140	0 3,5	10x50	12	Violet
	GLP4549	1 4,0	10x50	12	Violet
7x75 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2952	3/0 2,0	7x75	12	Violet
	GLP3545	2/0 3,0	7x75	12	Violet
	GLP4139	0 3,5	7x75	12	Violet

Antibacterial

Glikolak Plus

Surgical Suture, Poly[glycolide(90%)-co-lactide(10%)](PGLA)
Absorbable, Synthetic, Braided

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
10x75 cm					
GLIKOLAK Plus Non-Needle Suture	GLP2341	4/0 1,5	10x75	12	Violet
	GLP4141	0 3,5	10x75	12	Violet

150 cm					
GLIKOLAK Plus Non-Needle Suture	GLP1737	5/0 1,0	150	12	Violet
	GLP2326	4/0 1,5	150	12	Violet
	GLP2934	3/0 2,0	150	12	Undyed
	GLP2935	3/0 2,0	150	12	Violet
	GLP3511	2/0 3,0	150	12	Violet
	GLP4114	0 3,5	150	12	Violet