


# Monamid

Surgical Suture, Monofilament Polyamide (PA)  
Nonabsorbable, Synthetic

**BOZ**<sup>®</sup>



## Definition

Synthetic Non-Absorbable

## Material

Nylon 6 and Nylon 6.6 (PA)

## Coating

Uncoated

## Structure

Monofilament

## Colour

Blue, Black, Undyed

## Size

EP : 0.2 to 5

USP : 10/0 to 2

## Tensile Strength Retention

Moderate Loss

## Mass Absorption

Permanent

## Sterilization Method

Ethylene Oxide Gas

## Main Indication

General Soft tissue approximation or ligation

Cardiovascular

Ophthalmic Surgery





Neurosurgical Surgery

Plastic Surgery



















# Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color		
<b>4 mm</b>							
<b>3/8</b> Round Body, Taper Point <b>70 Micron</b>			<b>MA1001</b>	<b>11/0</b> 0,1	10	12	Black
<b>3/8</b> Round Body, Taper Point <b>100 Micron</b>			<b>MA1023</b>	<b>10/0</b> 0,2	15	12	Black
<b>3/8</b> Round Body, Taper Point <b>150 Micron</b>			<b>MA1099</b>	<b>9/0</b> 0,3	30	12	Black
			<b>MA1129</b>	<b>8/0</b> 0,4	13	12	Black
<b>5 mm</b>							
<b>3/8</b> Round Body, Taper Point <b>100 Micron</b>			<b>MA1045</b>	<b>10/0</b> 0,2	13	12	Black
			<b>MA1107</b>	<b>9/0</b> 0,3	13	12	Black
<b>3/8</b> Round Body, Taper Point <b>130 Micron</b>			<b>MA1133</b>	<b>8/0</b> 0,4	13	12	Black
<b>6 mm</b>							
<b>3/8</b> Round Body, Taper Point <b>140 Micron</b>			<b>MA1162</b>	<b>8/0</b> 0,4	15	12	Black
<b>3/8</b> Spatula <b>150 Micron Micropoint</b>			<b>MA1065</b>	<b>10/0</b> 0,2	30	12	Black
<b>3/8</b> Spatula <b>150 Micron</b>			<b>MA1066</b>	<b>10/0</b> 0,2	30	12	Black
			<b>MA1068</b>	<b>10/0</b> 0,2	45	12	Black
			<b>MA1118</b>	<b>9/0</b> 0,3	30	12	Black
<b>6,6 mm</b>							
<b>1/2</b> Cobra <b>150 Micron</b>			<b>MA5520</b>	<b>10/0</b> 0,2	60	12	Blue
			<b>MA5521</b>	<b>9/0</b> 0,3	60	12	Blue
<b>10 mm</b>							
<b>1/2</b> Reverse Cutting <b>330 Micron</b>			<b>MA1488</b>	<b>6/0</b> 0,7	45	12	Black
<b>3/8</b> Reverse Cutting <b>330 Micron</b>			<b>MA1495</b>	<b>6/0</b> 0,7	45	12	Black







# Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>13 mm</b>					
<b>3/8</b> Reverse Cutting 330 Micron 	<b>MA1612</b>	<b>6/0</b> 0,7	45	12	Blue
	<b>MA1619</b>	<b>6/0</b> 0,7	75	12	Blue
<b>16 mm</b>					
<b>1/2</b> Reverse Cutting       <b>3/8</b> Reverse Cutting 	<b>MA2012</b>	<b>5/0</b> 1,0	75	12	Blue
	<b>MA1672</b>	<b>6/0</b> 0,7	75	12	Blue
	<b>MA1667</b>	<b>6/0</b> 0,7	45	12	Black
	<b>MA2017</b>	<b>5/0</b> 1,0	45	12	Blue
	<b>MA2536</b>	<b>4/0</b> 1,5	45	12	Blue
	<b>MA2542</b>	<b>4/0</b> 1,5	75	12	Black
	<b>MA3021</b>	<b>3/0</b> 2,0	45	12	Black
	<b>MA3027</b>	<b>3/0</b> 2,0	75	12	Blue
<b>19 mm</b>					
<b>3/8</b> Reverse Cutting 	<b>MA2201</b>	<b>5/0</b> 1,0	45	12	Blue
	<b>MA2207</b>	<b>5/0</b> 1,0	75	12	Blue
	<b>MA2638</b>	<b>4/0</b> 1,5	45	12	Blue
	<b>MA2642</b>	<b>4/0</b> 1,5	75	12	Blue
<b>3/8</b> Inside Cutting 	<b>MA3120</b>	<b>3/0</b> 2,0	45	12	Blue
	<b>MA3121</b>	<b>3/0</b> 2,0	45	12	Black

# Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic








Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>19 mm</b>					
<b>3/8</b> Inside Cutting  	<b>MA3125</b>	<b>3/0</b> 2,0	75	12	Blue
	<b>MA2659</b>	<b>4/0</b> 1,5	45	12	Blue
	<b>MA2663</b>	<b>4/0</b> 1,5	75	12	Blue
<b>20 mm</b>					
<b>1/2</b> Reverse Cutting  	<b>MA3162</b>	<b>3/0</b> 2,0	75	12	Blue
<b>3/8</b> Reverse Cutting  	<b>MA3165</b>	<b>3/0</b> 2,0	45	12	Black
<b>24 mm</b>					
<b>3/8</b> Inside Cutting  	<b>MA3235</b>	<b>3/0</b> 2,0	45	12	Blue
<b>25 mm</b>					
<b>1/2</b> Tapercutting  	<b>MA5513</b>	<b>3/0</b> 2,0	45	12	Blue
<b>3/8</b> Reverse Cutting  	<b>MA2802</b>	<b>4/0</b> 1,5	45	12	Blue
	<b>MA3758</b>	<b>2/0</b> 3,0	75	12	Blue
<b>3/8</b> Inside Cutting  	<b>MA2811</b>	<b>4/0</b> 1,5	45	12	Blue
<b>26 mm</b>					
<b>1/2</b> Reverse Cutting  	<b>MA3812</b>	<b>2/0</b> 3,0	75	12	Blue
<b>3/8</b> Reverse Cutting  	<b>MA3344</b>	<b>3/0</b> 2,0	45	12	Blue
	<b>MA3347</b>	<b>3/0</b> 2,0	75	12	Black
	<b>MA3348</b>	<b>3/0</b> 2,0	75	12	Blue

# Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic








Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>26 mm</b>					
<b>3/8</b> Reverse Cutting 	<b>MA3820</b>	<b>2/0</b> 3,0	45	12	Blue
	<b>MA3823</b>	<b>2/0</b> 3,0	75	12	Blue
<b>3/8</b> Inside Cutting 	<b>MA3358</b>	<b>3/0</b> 2,0	45	12	Blue
<b>30 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>MA4254</b>	<b>0</b> 3,5	100	12	Blue
	<b>MA4255</b>	<b>0</b> 3,5	100	12	Black
	<b>MA4615</b>	<b>1</b> 4,0	100	12	Blue
	<b>MA4610</b>	<b>1</b> 4,0	75	12	Blue
<b>1/2</b> Reverse Cutting 	<b>MA3392</b>	<b>3/0</b> 2,0	90	12	Blue
	<b>MA3884</b>	<b>2/0</b> 3,0	75	12	Blue
<b>3/8</b> Reverse Cutting 	<b>MA3398</b>	<b>3/0</b> 2,0	75	12	Blue
	<b>MA3399</b>	<b>3/0</b> 2,0	75	12	Black
	<b>MA3893</b>	<b>2/0</b> 3,0	75	12	Black
	<b>MA3894</b>	<b>2/0</b> 3,0	75	12	Blue
	<b>MA4276</b>	<b>0</b> 3,5	75	12	Blue
	<b>MA4631</b>	<b>1</b> 4,0	75	12	Blue
<b>3/8</b> Inside Cutting 	<b>MA3902</b>	<b>2/0</b> 3,0	75	12	Blue

# Monamid

Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>35 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>MA4309</b>	<b>0</b> 3,5	100	12	Black
	<b>MA4298</b>	<b>0</b> 3,5	45	12	Blue
	<b>MA4652</b>	<b>1</b> 4,0	100	12	Black
<b>1/2</b> Reverse Cutting 	<b>MA4317</b>	<b>0</b> 3,5	75	12	Black
<b>3/8</b> Reverse Cutting 	<b>MA3936</b>	<b>2/0</b> 3,0	75	12	Blue
	<b>MA4664</b>	<b>1</b> 4,0	100	12	Blue
<b>38 mm</b>					
<b>1/2</b> Reverse Cutting 	<b>MA5512</b>	<b>2/0</b> 3,0	75	12	Black
	<b>MA5510</b>	<b>0</b> 3,5	75	12	Black
<b>40 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>MA4377</b>	<b>0</b> 3,5	100	12	Blue
	<b>MA4370</b>	<b>0</b> 3,5	75	12	Black
	<b>MA4725</b>	<b>1</b> 4,0	100	12	Black
	<b>MA4726</b>	<b>1</b> 4,0	100	12	Blue
<b>1/2</b> Round Body, Taper Point <b>LOOP</b> 	<b>MA4387</b>	<b>0</b> 3,5	150	12	Blue
<b>1/2</b> Reverse Cutting 	<b>MA4401</b>	<b>0</b> 3,5	100	12	Black
	<b>MA4402</b>	<b>0</b> 3,5	100	12	Blue

# Monamid Surgical Suture, Monofilament Polyamide (PA) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>40 mm</b>					
<b>1/2</b> Reverse Cutting Heavy 	<b>MA4405</b>	<b>0</b> 3,5	100	12	Black
<b>3/8</b> Reverse Cutting 	<b>MA4758</b>	<b>1</b> 4,0	75	12	Blue
<b>45 mm</b>					
<b>1/2</b> Round Body, Taper Point, LOOP 	<b>MA4781</b>	<b>1</b> 4,0	150	12	Blue
<b>3/8</b> Reverse Cutting 	<b>MA4789</b>	<b>1</b> 4,0	75	12	Blue
<b>50 mm</b>					
<b>1/2</b> Round Body, Taper Point LOOP 	<b>MA4464</b>	<b>0</b> 3,5	150	12	Black
	<b>MA4465</b>	<b>0</b> 3,5	150	12	Blue
	<b>MA4828</b>	<b>1</b> 4,0	150	12	Black
	<b>MA4829</b>	<b>1</b> 4,0	150	12	Blue
<b>60 mm</b>					
<b>Straight</b> Reverse Cutting 	<b>MA4077</b>	<b>2/0</b> 3,0	75	12	Blue
<b>61 mm</b>					
<b>Straight</b> Reverse Cutting 	<b>MA4083</b>	<b>2/0</b> 3,0	100	12	Blue