

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Surgical Sutures
Absorbable



Boz®

Definition

Short-term synthetic absorbable, closure of superficial wounds

Material

Poly(glycolide (75%)- co-caprolactone (25%)) (PGCL)

Coating

No coating

Structure

Monofilament

Colour

Undyed Beige (undyed) or Violet (dyed)

Size

EP : 0.7 to 5

USP : 6/0 to 2

Tensile Strength Retention

1. Week 45 % - 70 %
2. Week 20 %

Mass Absorption

90 - 120 Days

Sterilization Method



















Ethylene Oxide Gas

Main Indication

Gynaecology
Obstetrics
Urology
Plastic Surgery
Odontology (oral mucosa)
Skin closer - intracutaneous, subcutaneous.
Reconstructive Surgery
Gastrointestinal Surgery



















Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 330 Micron			MK1398	6/0 0,7	45	12	Undyed
			MK1759	5/0 1,0	45	12	Undyed
10 mm							
3/8 Reverse Cutting 330 Micron			MK1494	6/0 0,7	45	12	Undyed
13 mm							
3/8 Round Body, Taper Point 330 Micron			MK1578	6/0 0,7	45	12	Undyed
3/8 Round Body, Taper Point 330 Micron			MK1905	5/0 1,0	50	12	Undyed
16 mm							
1/2 Round Body, Taper Point			MK1989	5/0 1,0	75	12	Undyed
			MK2503	4/0 1,5	75	12	Undyed
			MK3004	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			MK2139	5/0 1,0	45	12	Undyed
			MK2557	4/0 1,5	75	12	Undyed
			MK3038	3/0 2,0	75	12	Undyed
18 mm							
3/8 Reverse Cutting			MK1704	6/0 0,7	45	12	Undyed
19 mm							
3/8 Reverse Cutting			MK2639	4/0 1,5	75	12	Undyed
			MK3108	3/0 2,0	75	12	Undyed













Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body,Taper Point			MK2683	4/0 1,5	75	12	Undyed
			MK3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2718	4/0 1,5	75	12	Undyed
			MK3167	3/0 2,0	75	12	Undyed
24 mm							
3/8 Inside Cutting			MK3236	3/0 2,0	75	12	Undyed
25 mm							
1/2 Round Body,Taper Point			MK2787	4/0 1,5	75	12	Undyed
			MK3251	3/0 2,0	75	12	Undyed
			MK3252	3/0 2,0	75	12	Violet
			MK3723	2/0 3,0	75	12	Undyed
3/8 Reverse Cutting			MK2805	4/0 1,5	75	12	Undyed
			MK3273	3/0 2,0	75	12	Undyed
3/8 Inside Cutting			MK3762	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body,Taper Point			MK3312	3/0 2,0	75	12	Undyed
			MK3313	3/0 2,0	75	12	Violet
			MK3785	2/0 3,0	75	12	Undyed
1/2 Round Body,Taper Point			MK5846	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting			MK3335	3/0 2,0	75	12	Undyed





Monokaprol




Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament


Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm						
<div>1/2</div> <div>Reverse Cutting</div>	 	MK4229	0 3,5	75	12	Undyed
		MK2851	4/0 1,5	75	12	Undyed
<div>3/8</div> <div>Reverse Cutting</div>	 	MK3345	3/0 2,0	75	12	Undyed
		MK3821	2/0 3,0	75	12	Undyed
30 mm						
<div>1/2</div> <div>Round Body,Taper Point</div>	 	MK3368	3/0 2,0	75	12	Undyed
		MK3369	3/0 2,0	75	12	Violet
		MK3857	2/0 3,0	75	12	Undyed
		MK3858	2/0 3,0	75	12	Violet
		MK4243	0 3,5	75	12	Undyed
		MK4607	1 4,0	75	12	Undyed
<div>3/8</div> <div>Round Body,Taper Point</div>	 	MK4263	0 3,5	75	12	Undyed
<div>3/8</div> <div>Reverse Cutting</div>	 	MK3890	2/0 3,0	75	12	Undyed
		MK4272	0 3,5	75	12	Undyed
		MK4628	1 4,0	75	12	Undyed
35 mm						
<div>1/2</div> <div>Round Body,Taper Point</div>	 	MK3915	2/0 3,0	75	12	Undyed
		MK4299	0 3,5	75	12	Undyed
		MK4643	1 4,0	75	12	Undyed


Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
35 mm					
<div>1/2</div> <div>Reverse Cutting</div> <div></div> <div></div>	MK4658	<div>1</div> <div>4,0</div>	75	12	Undyed
<div>3/8</div> <div>Reverse Cutting</div> <div></div> <div></div>	MK3932	<div>2/0</div> <div>3,0</div>	75	12	Undyed




40 mm						
<div>1/2</div> <div>Round Body, Taper Point</div> <div></div>	MK3981	2/0 3,0	75	12	Undyed	
	MK4367	0 3,5	75	12	Undyed	
	MK4374	0 3,5	90	12	Undyed	
	MK4712	1 4,0	75	12	Undyed	
	MK4713	1 4,0	75	12	Violet	
	MK4718	1 4,0	90	12	Undyed	
	MK5001	2 5,0	75	12	Undyed	
<div>1/2</div> <div>Reverse Cutting</div> <div></div>	MK4744	1 4,0	75	12	Undyed	
	MK5018	2 5,0	75	12	Undyed	
<div>3/8</div> <div>Reverse Cutting</div> <div></div>	MK4407	0 3,5	75	12	Undyed	

45 mm						
<div>1/2</div> <div>Round Body, Taper Point</div> <div></div>	MK4771	1 4,0	75	12	Undyed	
	MK5034	2 5,0	75	12	Undyed	

50 mm					
<div><div>1/2</div><div>Round Body,Taper Point</div><div></div></div>	MK4037	<div>2/0</div> <div>3,0</div>	75	12	Undyed
	MK4811	<div>1</div> <div>4,0</div>	75	12	Undyed

Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	MK4817	1 4,0	90	12	Undyed
	MK5065	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	MK4474	0 3,5	75	12	Undyed
	MK4840	1 4,0	75	12	Undyed
63 mm					
3/8 Blunt 	MK4859	1 4,0	100	12	Undyed
	MK4857	1 4,0	75	12	Undyed
Non-Needle Suture	Product No	USP EP	Length (cm)	Unit	Color
150 cm					
MONOKAPROL Non-Needle Suture	MK2934	3/0 2,0	150	12	Undyed
	MK3510	2/0 3,0	150	12	Undyed
	MK4113	0 3,5	150	12	Undyed
	MK4531	1 4,0	150	12	Undyed
	MK4907	2 5,0	150	12	Undyed