

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, antibacterial, closure of superficial wounds

Material

Poly(glycolide (75%)- co-caprolactone (25%) (PGCL)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

No coating

Structure

Monofilament

Colour

Undyed Beige (undyed) or Violet (dyed)

Size

EP : 1 to 5

USP : 5/0 to 2

Tensile Strenght Retention

1. Week 70 % - 75 % - 2. Week 40 % - 50%

Mass Absorption

90 - 120 Days

Sterilization Method

Ethylene Oxide Gas

Main Indication

Gynaecology

Obstetrics

Urology


Plastic Surgery

Odontology (oral mucosa)

Skin closer - intracutaneous, subcutaneous.


Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 330 Micron			MKP1759	5/0 1,0	45	12	Undyed
13 mm							
3/8 Round Body, Taper Point 330 Micron			MKP1905	5/0 1,0	50	12	Undyed
16 mm							
1/2 Round Body, Taper Point			MKP1985	5/0 1,0	45	12	Undyed
			MKP1989	5/0 1,0	75	12	Undyed
			MKP2503	4/0 1,5	75	12	Undyed
			MKP3004	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MKP2019	5/0 1,0	75	12	Undyed
			MKP2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			MKP2139	5/0 1,0	45	12	Undyed
			MKP2557	4/0 1,5	75	12	Undyed
			MKP3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			MKP2639	4/0 1,5	75	12	Undyed
			MKP3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			MKP2683	4/0 1,5	75	12	Undyed
			MKP3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MKP2718	4/0 1,5	75	12	Undyed
			MKP3167	3/0 2,0	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	MKP2787	4/0 1,5	75	12	Undyed
	MKP3251	3/0 2,0	75	12	Undyed
	MKP3723	2/0 3,0	75	12	Undyed
3/8 Inside Cutting 	MKP3762	2/0 3,0	75	12	Undyed
26 mm					
1/2 Round Body, Taper Point 	MKP3312	3/0 2,0	75	12	Undyed
	MKP3785	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting 	MKP3335	3/0 2,0	75	12	Undyed
30 mm					
1/2 Round Body, Taper Point 	MKP3368	3/0 2,0	75	12	Undyed
	MKP3857	2/0 3,0	75	12	Undyed
	MKP4243	0 3,5	75	12	Undyed
	MKP4249	0 3,5	90	12	Undyed
	MKP4607	1 4,0	75	12	Undyed
1/2 Round Body, Taper Point 	MKP3853	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting 	MKP4265	0 3,5	75	12	Undyed
3/8 Round Body, Taper Point 	MKP4263	0 3,5	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament





Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
3/8 Reverse Cutting 	MKP3890	2/0 3,0	75	12	Undyed
	MKP4272	0 3,5	75	12	Undyed

35 mm					
1/2 Round Body, Taper Point 	MKP3915	2/0 3,0	75	12	Undyed
	MKP4299	0 3,5	75	12	Undyed
	MKP4643	1 4,0	75	12	Undyed
3/8 Reverse Cutting 	MKP3932	2/0 3,0	75	12	Undyed
	MKP4319	0 3,5	75	12	Undyed

40 mm					
1/2 Round Body, Taper Point 	MKP3981	2/0 3,0	75	12	Undyed
	MKP4367	0 3,5	75	12	Undyed
	MKP4374	0 3,5	90	12	Undyed
	MKP4723	1 4,0	100	12	Undyed
	MKP4712	1 4,0	75	12	Undyed
	MKP4718	1 4,0	90	12	Undyed
	MKP5001	2 5,0	75	12	Undyed
1/2 Reverse Cutting 	MKP4744	1 4,0	75	12	Undyed
	MKP5018	2 5,0	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
3/8 Reverse Cutting 	MKP4407	0 3,5	75	12	Undyed
	MKP4755	1 4,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	MKP4421	0 3,5	75	12	Undyed
	MKP4771	1 4,0	75	12	Undyed
3/8 Reverse Cutting 	MKP4431	0 3,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	MKP4037	2/0 3,0	75	12	Undyed
	MKP4453	0 3,5	75	12	Undyed
	MKP4811	1 4,0	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	MKP4817	1 4,0	90	12	Undyed
	MKP5065	2 5,0	75	12	Undyed
63 mm					
3/8 Blunt 	MKP4857	1 4,0	75	12	Undyed

Antibacterial Monokaprol Plus

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)
Absorbable, Synthetic, Monofilament

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
MONOKAPROL Plus Non-Needle Suture	MKP2934	3/0 2,0	150	12	Undyed
	MKP3510	2/0 3,0	150	12	Undyed
	MKP4113	0 3,5	150	12	Undyed
	MKP4531	1 4,0	150	12	Undyed
	MKP4907	2 5,0	150	12	Undyed