


Surgical Sutures

Non-Absorbable

Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament



Boz®

Definition

Synthetic Non-Absorbable

Material

(PTFE) Polytetrafluoroethylene

Coating

No Coating

Structure

Monofilament

Colour

White

Size

EP : 0.5 to 2

USP : 7/0 to 3/0

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

General Soft Tissue Approximation and/or Ligation



Cardiovascular Surgery

Dental and General Surgeries as

Well as Repair of the Dura Mater


Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/2 Round Body, Taper Point 280 Micron			MT5770	6/0 0,7	45	12	White
8.5 mm							
3/8 Reverse Cutting Premium 280 Micron			MT5768	7/0 0,5	55	12	White
10 mm							
3/8 Reverse Cutting Premium 330 Micron			MT5767	6/0 0,7	55	12	White
12 mm							
3/8 Reverse Cutting 330 Micron			MT1540	6/0 0,7	45	12	White
12.3 mm							
3/8 Reverse Cutting Premium 380 Micron			MT5763	5/0 1,0	55	12	White
3/8 Reverse Cutting Premium 430 Micron			MT5765	4/0 1,5	75	12	White
13 mm							
3/8 Round Body, Taper Point 380 Micron			MT5852	5/0 1,0	75	12	White
3/8 Round Body, Taper Point 430 Micron			MT5854	4/0 1,5	75	12	White
3/8 CC Needle 380 Micron			MT5853	5/0 1,0	75	12	White
3/8 CC Needle 430 Micron			MT5855	4/0 1,5	75	12	White
15 mm							
3/8 Round Body, Taper Point			MT5856	3/0 2,0	75	12	White




Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
15 mm							
3/8 Reverse Cutting			MT5769	6/0 0,7	45	12	White
3/8 CC Needle			MT5857	3/0 2,0	75	12	White
16 mm							
3/8 Reverse Cutting			MT5771	5/0 1,0	45	12	White
			MT5772	4/0 1,5	45	12	White
3/8 Inside Cutting			MT5755	5/0 1,0	45	12	White
16.2 mm							
3/8 Reverse Cutting Premium 430 Micron			MT5766	4/0 1,5	75	12	White
17 mm							
1/2 Round Body, Taper Point			MT5860	5/0 1,0	75	12	White
1/2 CC Needle			MT5861	5/0 1,0	75	12	White
18 mm							
3/8 Reverse Cutting			MT5829	3/0 2,0	45	12	White
			MT5558	3/0 2,0	75	12	White
18.7 mm							
3/8 Reverse Cutting 560 Micron			MT5764	3/0 2,0	75	12	White
19 mm							
3/8 Reverse Cutting			MT5773	4/0 1,5	45	12	White

Monotef

Surgical Suture, Polytetrafluoroethylene (PTFE)
Nonabsorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
3/8 Reverse Cutting			MT5566	2/0 3,0	75	12	White
25 mm							
3/8 Round Body, Taper Point			MT5858	2/0 3,0	75	12	White
3/8 CC Needle			MT5859	2/0 3,0	75	12	White
30 mm							
1/2 Round Body, Taper Point			MT5862	0 3,5	75	12	White
1/2 CC Needle			MT5863	0 3,5	75	12	White