

## Surgical Sutures

### Non - Absorbable

# Monowire

Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament



**Boz**®

#### Definition

Non-Absorbable

#### Material

316L Stainless Steel conforming to  
ISO 5832 Part-1

#### Coating

FEP or PE

#### Structure

Monofilament or Multifilament

#### Colour

Metallic, Blue, Transperant, Orange

#### Size

EP : 0.7 to 10  
USP : 6/0 to 8

#### Tensile Strenght Retention

Permanent

#### Mass Absorption

Permanent

#### Sterilization Method




Ethylene Oxide Gas

#### Main Indication

Sternum Closure  
Orthopaedics  
Abdominal wound closure  
Hernia repair



# Monowire

Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament

| Needle Length and description                    |                                                                                     |                                                                                     | Product No    | USP EP            | Length (cm) | Unit Pcs | Color |
|--------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------|-------------------|-------------|----------|-------|
| <b>30 mm</b>                                     |                                                                                     |                                                                                     |               |                   |             |          |       |
| <b>1/2</b><br>Reverse Cutting                    |    |    | <b>MW3881</b> | <b>2/0</b><br>3,0 | 75          | 12       | Mono  |
| <b>40 mm</b>                                     |                                                                                     |                                                                                     |               |                   |             |          |       |
| <b>1/2</b><br>Reverse Cutting                    |    |    | <b>MW4743</b> | <b>1</b><br>4,0   | 4 X 45      | 12       | Mono  |
| <b>1/2</b><br>Tapercutting                       |    |    | <b>MW4707</b> | <b>1</b><br>4,0   | 75          | 12       | Mono  |
| <b>48 mm</b>                                     |                                                                                     |                                                                                     |               |                   |             |          |       |
| <b>1/2</b><br>Reverse Cutting                    |    |    | <b>MW5056</b> | <b>2</b><br>5,0   | 4 X 45      | 12       | Mono  |
|                                                  |                                                                                     |                                                                                     | <b>MW5057</b> | <b>2</b><br>5,0   | 45          | 12       | Mono  |
| <b>1/2</b><br>Reverse Cutting<br>Heavy           |  |  | <b>MW5126</b> | <b>3</b><br>6,0   | 75          | 12       | Mono  |
|                                                  |                                                                                     |                                                                                     | <b>MW5153</b> | <b>4</b><br>6,0   | 75          | 12       | Mono  |
| <b>1/2</b><br>Tapercutting<br>Heavy              |  |  | <b>MW5163</b> | <b>5</b><br>7,0   | 45          | 12       | Mono  |
|                                                  |                                                                                     |                                                                                     | <b>MW5198</b> | <b>6</b><br>8,0   | 45          | 12       | Mono  |
| <b>1/2</b><br>Tapercutting<br>Heavy, Rotating    |  |  | <b>MW5146</b> | <b>4</b><br>6,0   | 45          | 12       | Mono  |
| <b>55 mm</b>                                     |                                                                                     |                                                                                     |               |                   |             |          |       |
| <b>1/2</b><br>Reverse Cutting<br>Heavy           |  |  | <b>MW5131</b> | <b>3</b><br>6,0   | 75          | 12       | Mono  |
| <b>1/2</b><br>Reverse Cutting<br>Rotating        |  |  | <b>MW5087</b> | <b>2</b><br>5,0   | 75          | 12       | Mono  |
| <b>1/2</b><br>Reverse Cutting<br>Heavy, Rotating |  |  | <b>MW5156</b> | <b>4</b><br>6,0   | 2 X 75      | 12       | Mono  |
|                                                  |                                                                                     |                                                                                     | <b>MW5187</b> | <b>5</b><br>7,0   | 45          | 12       | Mono  |

# Monowire

Surgical Stainless Steel Wire Nonabsorbable, Monofilament or Multifilament

| Needle Length and description                                                                                                        | Product No | USP EP   | Length (cm) | Unit Pcs | Color |
|--------------------------------------------------------------------------------------------------------------------------------------|------------|----------|-------------|----------|-------|
| <b>55 mm</b>                                                                                                                         |            |          |             |          |       |
| <b>1/2</b><br>Reverse Cutting<br>Heavy, Rotating    | MW5188     | 5<br>7,0 | 2 X 45      | 12       | Mono  |
|                                                                                                                                      | MW5189     | 5<br>7,0 | 4 X 45      | 12       | Mono  |
|                                                                                                                                      | MW5190     | 5<br>7,0 | 75          | 12       | Mono  |
|                                                                                                                                      | MW5191     | 5<br>7,0 | 2 X 75      | 12       | Mono  |
|                                                                                                                                      | MW5192     | 5<br>7,0 | 4 X 75      | 12       | Mono  |
|                                                                                                                                      | MW5206     | 6<br>8,0 | 4 X 75      | 12       | Mono  |
|                                                                                                                                      | MW5205     | 6<br>8,0 | 75          | 12       | Mono  |
| <b>1/2</b><br>Reverse Cutting<br>Heavy, Rotating   | MW5193     | 5<br>7,0 | 75          | 12       | Mono  |
| <b>1/2</b><br>Trocac Point<br>Heavy, Rotating     | MW5777     | 5<br>7,0 | 4 X 45      | 12       | Mono  |
|                                                                                                                                      | MW5179     | 5<br>7,0 | 75          | 12       | Mono  |
|                                                                                                                                      | MW5180     | 5<br>7,0 | 2 X 75      | 12       | Mono  |
|                                                                                                                                      | MW5203     | 6<br>8,0 | 75          | 12       | Mono  |
|                                                                                                                                      | MW5204     | 6<br>8,0 | 2 X 75      | 12       | Mono  |
|                                                                                                                                      | MW5208     | 7<br>9,0 | 75          | 12       | Mono  |
|                                                                                                                                      | MW5209     | 7<br>9,0 | 2 X 75      | 12       | Mono  |
| <b>57 mm</b>                                                                                                                         |            |          |             |          |       |
| <b>1/2</b><br>Reverse Cutting<br>Heavy, Rotating  | MW5195     | 5<br>7,0 | 75          | 12       | Mono  |
|                                                                                                                                      | MW5196     | 5<br>7,0 | 2 X 75      | 12       | Mono  |