

Surgical Sutures Non-Absorbable

Multamid

Surgical Suture, Multifilament Polyamide (PA) Silicone Coated
Nonabsorbable, Synthetic, Black

Boz®



Definition

Synthetic Non-Absorbable

Material

Nylon 6 and Nylon 6.6 (PA)

Coating

Silicone, Paraffin Wax Uncoated

Structure

Braided

Colour

Black, Blue, Undyed

Size

EP : 0.01 to 7

USP : 12/0 to 5

Tensile Strenght Retention

Moderate Loss

Mass Absorption

Permanent

Sterilization Method




Ethylene Oxide Gas

Main Indication

General Soft tissue
approximation or ligation
Cardiovascular
Ophthalmic Surgery
Neurosurgical Surgery
Plastic Surgery


Multamid

Surgical Suture, Multifilament Polyamide (PA)
Silicone Coated Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			MU3006	3/0 2,0	75	12	Black
26 mm							
1/2 Round Body, Taper Point			MU2824	4/0 1,5	8 X 45	6	Black
			MU3315	3/0 2,0	75	12	Black
			MU3789	2/0 3,0	75	12	Black
			MU3792	2/0 3,0	8 X 75	6	Black
			MU4218	0 3,5	75	12	Black
			MU4590	1 4,0	75	12	Black
3/8 Reverse Cutting			MU3343	3/0 2,0	45	12	Black
			MU4238	0 3,5	75	12	Black
30 mm							
1/2 Round Body, Taper Point			MU3371	3/0 2,0	75	12	Black
			MU3859	2/0 3,0	75	12	Black
			MU4609	1 4,0	75	12	Black
35 mm							
1/2 Round Body, Taper Point			MU4301	0 3,5	75	12	Black
3/8 Reverse Cutting			MU3425	3/0 2,0	75	12	Black
			MU3934	2/0 3,0	75	12	Black
			MU4663	1 4,0	75	12	Black

Multamid

Surgical Suture, Multifilament Polyamide (PA)
Silicone Coated Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
90 mm					
 3/8 Inside Cutting Vinyl Tube ØOutside: 4 mm ØInside: 2.5 mm Length: 40 mm	MU4886	1 4,0	2 X 50	12	Black

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit	Color
17x45 cm					
MULTAMID Non-Needle Suture	MU1740	5/0 1,0	17x45	12	Black
	MU2351	4/0 1,5	17x45	12	Black
	MU2964	3/0 2,0	17x45	12	Black
	MU3560	2/0 3,0	17x45	12	Black
	MU4561	1 4,0	17x45	12	Black

10x75 cm					
MULTAMID Non-Needle Suture	MU2340	4/0 1,5	10x75	12	Black
	MU2956	3/0 2,0	10x75	12	Black
	MU3549	2/0 3,0	10x75	12	Black
	MU4552	1 4,0	10x75	12	Black