

Plain Catgut

Absorbable Surgical Suture



Boz®

Definition

Collagen based biomaterial, absorbable.

Material

Collagen tissue derived from purified sheep intestine

Coating

No coating

Structure

Monofilament

Colour

Yellowish

Size

EP : 1 to 8

USP : 6/0 to 4

Sterilization Method

Ethylene Oxide Gas

Gamma Irradiation


Main Indication

Ligatures

Ophthalmic Surgery

Plain Catgut

Absorbable Surgical Suture



Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 330 Micron			CP1563	6/0 1,0	75	12	Yellowish
1/2 Round Body, Taper Point 430 Micron			CP1892	5/0 1,5	75	12	Yellowish
1/2 Round Body, Taper Point 330 Micron With Alcohol			CP1563A	6/0 1,0	75	12	Yellowish
3/8 Reverse Cutting 430 Micron			CP1933	5/0 1,5	75	12	Yellowish
3/8 Reverse Cutting 430 Micron With Alcohol			CP1933A	5/0 1,5	75	12	Yellowish
16 mm							
1/2 Round Body, Taper Point			CP1989	5/0 1,5	75	12	Yellowish
			CP2503	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol			CP1989A	5/0 1,5	75	12	Yellowish
			CP2503A	4/0 2,0	75	12	Yellowish
			CP3004A	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting			CP2019	5/0 1,5	75	12	Yellowish
			CP2539	4/0 2,0	75	12	Yellowish
			CP3024	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol			CP2019A	5/0 1,5	75	12	Yellowish
			CP2539A	4/0 2,0	75	12	Yellowish
			CP3024A	3/0 3,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point	 		CP2227	5/0 1,5	50	12	Yellowish
			CP2228	5/0 1,5	75	12	Yellowish
			CP2683	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol	 		CP2683A	4/0 2,0	75	12	Yellowish
3/8 Reverse Cutting	 		CP2718	4/0 2,0	75	12	Yellowish
			CP3167	3/0 3,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol	 		CP2718A	4/0 2,0	75	12	Yellowish
25 mm							
1/2 Round Body, Taper Point With Alcohol	 		CP2787A	4/0 2,0	75	12	Yellowish
3/8 Reverse Cutting	 		CP2276	5/0 1,5	75	12	Yellowish
26 mm							
1/2 Round Body, Taper Point	 		CP2827	4/0 2,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol	 		CP2827A	4/0 2,0	75	12	Yellowish
			CP3312A	3/0 3,0	75	12	Yellowish
			CP3785A	2/0 3,5	75	12	Yellowish
3/8 Reverse Cutting	 		CP2851	4/0 2,0	75	12	Yellowish
			CP3345	3/0 3,0	75	12	Yellowish
			CP3821	2/0 3,5	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
3/8 Reverse Cutting With Alcohol 	CP3345A	3/0 3,0	75	12	Yellowish
	CP3821A	2/0 3,5	75	12	Yellowish
	CP4236A	0 4,0	75	12	Yellowish
30 mm					
1/2 Round Body, Taper Point 	CP3368	3/0 3,0	75	12	Yellowish
	CP4243	0 4,0	75	12	Yellowish
	CP4607	1 5,0	75	12	Yellowish
	CP4954	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP2874A	4/0 2,0	75	12	Yellowish
	CP3368A	3/0 3,0	75	12	Yellowish
	CP3857A	2/0 3,5	75	12	Yellowish
	CP4243A	0 4,0	75	12	Yellowish
	CP4607A	1 5,0	75	12	Yellowish
3/8 Reverse Cutting 	CP3396	3/0 3,0	75	12	Yellowish
	CP3890	2/0 3,5	75	12	Yellowish
	CP4272	0 4,0	75	12	Yellowish
	CP4628	1 5,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture








Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
3/8 Reverse Cutting With Alcohol 	CP3396A	3/0 3,0	75	12	Yellowish
	CP3890A	2/0 3,5	75	12	Yellowish
	CP4272A	0 4,0	75	12	Yellowish
	CP4628A	1 5,0	75	12	Yellowish
35 mm					
1/2 Round Body, Taper Point 	CP3410	3/0 3,0	75	12	Yellowish
	CP3915	2/0 3,5	75	12	Yellowish
	CP4299	0 4,0	75	12	Yellowish
	CP4969	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP3410A	3/0 3,0	75	12	Yellowish
	CP4299A	0 4,0	75	12	Yellowish
	CP4643A	1 5,0	75	12	Yellowish
	CP4969A	2 6,0	75	12	Yellowish
1/2 Reverse Cutting 	CP4976	2 6,0	75	12	Yellowish
3/8 Reverse Cutting 	CP3422	3/0 3,0	75	12	Yellowish
	CP3932	2/0 3,5	75	12	Yellowish
	CP4319	0 4,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP3422A	3/0 3,0	75	12	Yellowish
	CP4319A	0 4,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture



Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	CP4367	0 4,0	75	12	Yellowish
	CP4374	0 4,0	90	12	Yellowish
	CP4712	1 5,0	75	12	Yellowish
	CP5001	2 6,0	75	12	Yellowish
	CP5118	3 7,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4367A	0 4,0	75	12	Yellowish
	CP4712A	1 5,0	75	12	Yellowish
	CP5001A	2 6,0	75	12	Yellowish
	CP5118A	3 7,0	75	12	Yellowish
1/2 Reverse Cutting With Alcohol 	CP5018A	2 6,0	75	12	Yellowish
3/8 Reverse Cutting 	CP4407	0 4,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP4407A	0 4,0	75	12	Yellowish
45 mm					
1/2 Round Body, Taper Point 	CP4421	0 4,0	75	12	Yellowish
	CP4771	1 5,0	75	12	Yellowish
	CP5034	2 6,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4771A	1 5,0	75	12	Yellowish
	CP5034A	2 6,0	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
45 mm					
3/8 Reverse Cutting 	CP4787	1 5,0	75	12	Yellowish
3/8 Reverse Cutting With Alcohol 	CP4787A	1 5,0	75	12	Yellowish
50 mm					
1/2 Round Body, Taper Point 	CP4453	0 4,0	75	12	Yellowish
	CP4811	1 5,0	75	12	Yellowish
	CP5065	2 6,0	75	12	Yellowish
	CP5127	3 7,0	75	12	Yellowish
1/2 Round Body, Taper Point With Alcohol 	CP4453A	0 4,0	75	12	Yellowish
	CP4811A	1 5,0	75	12	Yellowish
	CP5065A	2 6,0	75	12	Yellowish
	CP5127A	3 7,0	75	12	Yellowish
Straight Reverse Cutting 	CP3475	3/0 3,0	75	12	Yellowish
Straight Reverse Cutting With Alcohol 	CP3475A	3/0 3,0	75	12	Yellowish
60 mm					
Straight Reverse Cutting 	CP4073	2/0 3,5	75	12	Yellowish
Straight Reverse Cutting With Alcohol 	CP4073A	2/0 3,5	75	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotes the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
75 cm					
PLAIN CATGUT Non-Needle Suture	CP1748	5/0 1,5	75	12	Yellowish
	CP2316	4/0 2,0	75	12	Yellowish
	CP2925	3/0 3,0	75	12	Yellowish
	CP4104	0 4,0	75	12	Yellowish
	CP4521	1 5,0	75	12	Yellowish
	CP4898	2 6,0	75	12	Yellowish
75 cm With Alcohol					
PLAIN CATGUT Non-Needle Suture With Alcohol	CP1748A	5/0 1,5	75	12	Yellowish
	CP2316A	4/0 2,0	75	12	Yellowish
	CP2925A	3/0 3,0	75	12	Yellowish
	CP3530A	2/0 3,5	75	12	Yellowish
	CP4104A	0 4,0	75	12	Yellowish
	CP4898A	2 6,0	75	12	Yellowish
150 cm					
PLAIN CATGUT Non-Needle Suture	CP1736	5/0 1,5	150	12	Yellowish
	CP2325	4/0 2,0	150	12	Yellowish
	CP2934	3/0 3,0	150	12	Yellowish
	CP3510	2/0 3,5	150	12	Yellowish
	CP4113	0 4,0	150	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotates the product is "wet" e.g.

*After code number CP####A

Plain Catgut

Absorbable Surgical Suture

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
PLAIN CATGUT Non-Needle Suture	CP4531	1 5,0	150	12	Yellowish
	CP4907	2 6,0	150	12	Yellowish
	CP5106	3 7,0	150	12	Yellowish
	CP5134	4 8,0	150	12	Yellowish
150 cm With Alcohol					
PLAIN CATGUT Non-Needle Suture With Alcohol	CP1736A	5/0 4,5	150	12	Yellowish
	CP2325A	4/0 2,0	150	12	Yellowish
	CP2934A	3/0 3,0	150	12	Yellowish
	CP3510A	2/0 3,5	150	12	Yellowish
	CP4113A	0 3,5	150	12	Yellowish
	CP4531A	1 4,0	150	12	Yellowish
	CP4907A	2 6,0	150	12	Yellowish
	CP5106A	3 7,0	150	12	Yellowish
	CP5134A	4 8,0	150	12	Yellowish

*All Catgut Sutures are available in conditioning fluid with isopropyl alcohol. "A" denotesthe product is "wet" e.g.

*After code number CP####A