

Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable,
Synthetic, Monofilament, Violet

Surgical Sutures
Absorbable

Boz®



Definition

Synthetic absorbable, wound support through an extended healing period

Material

Polydioxanone (PDO)

Coating

No coating

Structure

Monofilament

Colour

Violet (dyed)

Size

EP : 0.5 to 5

USP : 7/0 to 2

Tensile Strength Retention

2. Week 70 %

Mass Absorption

180-240 days

Sterilization Method








Ethylene Oxide Gas

Main Indication

Orthopaedics
General Surgery
Abdominal Wall Closure
Paediatric cardiovascular tissues
Ophthalmic Surgery











Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
3/8 Round Body, Taper Point 220 Micron			PX1289	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX5444	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX1426	6/0 0,7	45	12	Violet
10 mm							
3/8 Reverse Cutting 330 Micron			PX5850	6/0 0,7	45	12	Violet
11 mm							
3/8 Round Body, Taper Point 220 Micron			PX1509	6/0 0,7	45	12	Violet
13 mm							
1/2 Round Body, Taper Point 330 Micron			PX1564	6/0 0,7	75	12	Violet
1/2 Round Body, Taper Point 380 Micron			PX1885	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point 430 Micron			PX1897	5/0 1,0	75	12	Violet
1/2 Spatula 280 Micron			PX1599	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 280 Micron			PX1362	7/0 0,5	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1575	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1576	6/0 0,7	45	12	Violet
3/8 Round Body, Taper Point 320 Micron			PX1128	6/0 0,7	75	12	Violet
3/8 Round Body, Taper Point 380 Micron			PX1914	5/0 0,7	75	12	Violet
3/8 Tapercutting 280 Micron			PX1609	6/0 0,7	45	12	Violet



Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body, Taper Point			PX1990	5/0 1,0	75	12	Violet
			PX2504	4/0 1,5	75	12	Violet
			PX3005	3/0 2,0	75	12	Violet
1/2 Round Body, Taper Point			PX1999	5/0 1,0	90	12	Violet
			PX2513	4/0 1,5	90	12	Violet
3/8 Round Body, Taper Point			PX2005	5/0 1,0	75	12	Violet
			PX2520	4/0 1,5	75	12	Violet
3/8 Round Body, Taper Point			PX2523	4/0 1,5	90	12	Violet
3/8 Reverse Cutting			PX2021	5/0 1,0	75	12	Violet
			PX2540	4/0 1,5	75	12	Violet
17 mm							
1/2 Round Body, Taper Point			PX5517	5/0 1,0	75	12	Violet
1/2 Round Body, Taper Point			PX5516	4/0 1,5	75	12	Violet
19 mm							
3/8 Reverse Cutting			PX2635	4/0 1,5	45	12	Violet
			PX3109	3/0 2,0	75	12	Violet
20 mm							
1/2 Round Body, Taper Point			PX2229	5/0 1,0	75	12	Violet
			PX2685	4/0 1,5	75	12	Violet


Polidiox



Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet


Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			PX3142	3/0 2,0	75	12	Violet
1/2 Inside Cutting			PX3175	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			PX2719	4/0 1,5	75	12	Violet
22 mm							
1/2 Inside Cutting			PX3217	3/0 2,0	45	12	Violet
25 mm							
1/2 Round Body, Taper Point			PX3252	3/0 2,0	75	12	Violet
			PX3724	2/0 3,0	75	12	Violet
			PX4190	0 3,5	75	12	Violet
26 mm							
1/2 Round Body, Taper Point			PX2828	4/0 1,5	75	12	Violet
			PX3313	3/0 2,0	75	12	Violet
			PX3786	2/0 3,0	75	12	Violet
1/2 Reverse Cutting			PX3336	3/0 2,0	75	12	Violet
			PX3813	2/0 3,0	75	12	Violet
3/8 Reverse Cutting			PX3341	3/0 2,0	45	12	Violet
			PX3346	3/0 2,0	75	12	Violet
3/8 Inside Cutting			PX3360	3/0 2,0	75	12	Violet

Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm					
1/2 Round Body, Taper Point  	PX3369	3/0 2,0	75	12	Violet
	PX3858	2/0 3,0	75	12	Violet
	PX4244	0 3,5	75	12	Violet
3/8 Reverse Cutting  	PX3891	2/0 3,0	75	12	Violet



40 mm					
1/2 Round Body, Taper Point  	PX3982	2/0 3,0	75	12	Violet
	PX4368	0 3,5	75	12	Violet
	PX4713	1 4,0	75	12	Violet
	PX4719	1 4,0	90	12	Violet
	PX5006	2 5,0	90	12	Violet
1/2 Round Body, Taper Point LOOP  	PX4388	0 3,5	150	12	Violet
	PX4736	1 4,0	150	12	Violet

45 mm					
1/2 Round Body, Taper Point Heavy  	PX4779	1 4,0	90	12	Violet

50 mm					
1/2 Round Body, Taper Point  	PX4459	0 3,5	90	12	Violet
	PX4818	1 4,0	90	12	Violet
	PX5070	2 5,0	90	12	Violet

Polidiox

Surgical Suture, Polydioxanone (PDO) Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point LOOP 	PX4827	1 4,0	150	12	Violet
1/2 Round Body, Taper Point Heavy, LOOP 	PX4466	0 3,5	150	12	Violet