

Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO) Absorbable,
Synthetic, Monofilament, Violet

Surgical Sutures
Absorbable

Boz®



Definition

Short-term synthetic absorbable, antibacterial, wound support through an extended healing period.

Material

Polydioxanone (PDO)

Antibacterial Agent:

CHA (Chlorhexidine diacetate)

Coating

No coating

Structure

Monofilament

Colour

Violet (dyed)

Size

EP : 1 to 5

USP : 5/0 to 2

Tensile Strenght Retention

2. Week 70%

4. Week 50%

Mass Absorption

180 - 240 Days

Sterilization Method

Ethylene Oxide Gas

Main Indication









Orthopaedics

General Surgery

Abdominal Wall Closure


Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
13 mm							
1/2 Round Body, Taper Point 430 Micron			PXP1897	5/0 1,0	75	12	Violet
16 mm							
1/2 Round Body, Taper Point			PXP1990	5/0 1,0	75	12	Violet
			PXP2504	4/0 1,5	75	12	Violet
			PXP3005	3/0 2,0	75	12	Violet
1/2 Round Body, Taper Point			PXP1999	5/0 1,0	90	12	Violet
			PXP2513	4/0 1,5	90	12	Violet
3/8 Round Body, Taper Point			PXP2005	5/0 1,0	75	12	Violet
			PXP2520	4/0 1,5	75	12	Violet
3/8 Round Body, Taper Point			PXP2523	4/0 1,5	90	12	Violet
3/8 Reverse Cutting			PXP2021	5/0 1,0	75	12	Violet
19 mm							
3/8 Reverse Cutting			PXP2635	4/0 1,5	45	12	Violet
20 mm							
1/2 Round Body, Taper Point			PXP2685	4/0 1,5	75	12	Violet
			PXP3142	3/0 2,0	75	12	Violet
3/8 Reverse Cutting			PXP2719	4/0 1,5	75	12	Violet
22 mm							
1/2 Inside Cutting			PXP3217	3/0 2,0	45	12	Violet






Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	PXP3252	3/0 2,0	75	12	Violet
	PXP3724	2/0 3,0	75	12	Violet
	PXP4190	0 3,5	75	12	Violet
26 mm					
1/2 Round Body, Taper Point 	PXP2828	4/0 1,5	75	12	Violet
	PXP3313	3/0 2,0	75	12	Violet
	PXP3786	2/0 3,0	75	12	Violet
3/8 Reverse Cutting 	PXP3346	3/0 2,0	75	12	Violet
3/8 Inside Cutting 	PXP3360	3/0 2,0	75	12	Violet
30 mm					
1/2 Round Body, Taper Point 	PXP3369	3/0 2,0	75	12	Violet
	PXP3858	2/0 3,0	75	12	Violet
	PXP4244	0 3,5	75	12	Violet
3/8 Reverse Cutting 	PXP3891	2/0 3,0	75	12	Violet
40 mm					
1/2 Round Body, Taper Point 	PXP3982	2/0 3,0	75	12	Violet
	PXP4368	0 3,5	75	12	Violet
	PXP4375	0 3,5	90	12	Violet
	PXP4713	1 4,0	75	12	Violet

Antibacterial Polidiox Plus

Surgical Suture, Polydioxanone (PDO)
Absorbable, Synthetic, Monofilament, Violet

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	PXP4719	1 4,0	90	12	Violet
	PXP5006	2 5,0	90	12	Violet
1/2 Round Body, Taper Point, LOOP 	PXP4388	0 3,5	150	12	Violet
	PXP4736	1 4,0	150	12	Violet
45 mm					
1/2 Round Body, Taper Point Heavy 	PXP4779	1 4,0	90	12	Violet
50 mm					
1/2 Round Body, Taper Point 	PXP4459	0 3,5	90	12	Violet
	PXP4818	1 4,0	90	12	Violet
	PXP5070	2 5,0	90	12	Violet
1/2 Round Body, Taper Point LOOP 	PXP4827	1 4,0	150	12	Violet
1/2 Round Body, Taper Point Heavy LOOP 	PXP4466	0 3,5	150	12	Violet