


Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET)
Nonabsorbable, Synthetic

Surgical Sutures
Non-Absorbable

Boz®



Definition

Synthetic Non- Absorbable

Material

Polyester Polyethylene Terephthalate (PET)

Coating

Silicone, Uncoated

Structure

Braided

Colour

Green, White

Size

EP : 0.5 to 7
USP : 7/0 to 5

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method


Ethylene Oxide Gas

Main Indication

Cardiac Valve Repair
General Surgery
Cardiovascular Procedures
Neurosurgical Procedures




Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/2 Spatula 340 Micron			PS1756	5/0 1,0	45	12	Green
1/4 Spatula 340 Micron			PS1408	6/0 0,7	45	12	Green
			PS1765	5/0 1,0	45	12	White
			PS1766	5/0 1,0	45	12	Green
3/8 Round Body, Taper Point 280 Micron			PS1302	7/0 0,5	45	12	Green
8,1 mm							
1/4 Sabre 380 Micron			PS5531	4/0 1,5	60	12	Green
10 mm							
3/8 Round Body, Taper Point 280 Micron			PS1482	6/0 0,7	75	12	Green
3/8 Spatula 280 Micron			PS1452	6/0 0,7	45	12	Green
13 mm							
3/8 Round Body, Taper Point 280 Micron			PS1573	6/0 0,7	75	12	Green
3/8 Round Body, Taper Point 320 Micron			PS1903	5/0 1,0	75	12	Green
3/8 Round Body, Taper Point 330 Micron			PS1587	6/0 0,7	75	12	Green
3/8 Round Body, Taper Point 330 Micron 6mm x 3mm x 3mm Teflon Pledget			PS1916	5/0 1,0	75	12	Green
3/8 Reverse Cutting 560 Micron			PS2455	4/0 1,5	45	12	White
3/8 Spatula 280 Micron			PS1604	6/0 0,7	45	12	Green











Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
16 mm							
1/2 Round Body,Taper Point			PS3008	3/0 2,0	75	12	Green
1/2 Round Body,Taper Point			PS2515	4/0 1,5	90	12	Green
			PS3011	3/0 2,0	90	12	Green
1/2 Round Body,Taper Point			PS3012	3/0 2,0	90	6	5 Green 5 White
			PS5923	3/0 2,0	90	12	Green
6mm x 3mm x 1.5mm Teflon Pledget							
3/8 Inside Cutting			PS2560	4/0 1,5	75	12	Green
			PS3040	3/0 2,0	75	12	Green
17 mm							
1/2 Round Body,Taper Point			PS2582	4/0 1,5	75	12	Green
			PS3609	2/0 3,0	75	12	Green
1/2 Round Body,Taper Point			PS2587	4/0 1,5	90	12	Green
			PS3066	3/0 2,0	90	12	Green
			PS3613	2/0 3,0	90	12	Green
			PS3615	2/0 3,0	90	6	5 Green 5 White
1/2 Round Body,Taper Point			PS2579	4/0 1,5	75	12	Green
			PS3606	2/0 3,0	75	6	5 Green 5 White
			PS3617	2/0 3,0	90	12	Green
3mm x 3mm x 1.5mm Teflon Pledget							

Polisil



Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
17 mm							
1/2 Round Body, Taper Point 6mm x 3mm x 1.5mm Teflon Pledget	 		PS3067	3/0 2,0	90	6	5 Green 5 White
			PS3616	2/0 3,0	90	6	5 Green 5 White
			PS3618	2/0 3,0	90	12	Green
1/2 Tapercutting	 		PS2570	4/0 1,5	90	12	Green
1/2 Tapercutting 7mm x 3mm x 1.5mm Teflon Pledget	 		PS3592	2/0 3,0	75	6	5 Green 5 White
3/8 Round Body, Taper Point	 		PS5762	3/0 2,0	75	12	Green
20 mm							
1/2 Round Body, Taper Point	 		PS3648	2/0 3,0	75	12	Green
1/2 Round Body, Taper Point	 		PS3650	2/0 3,0	75	12	Green
1/2 Round Body, Taper Point	 		PS2688	4/0 1,5	75	12	Green
			PS2695	4/0 1,5	90	12	Green
			PS5843	3/0 2,0	75	12	Green
			PS3148	3/0 2,0	90	12	Green
1/2 Round Body, Taper Point 6mm x 3mm x 1.5mm Teflon Pledget	 		PS3651	2/0 3,0	75	6	5 Green 5 White
			PS3653	2/0 3,0	90	12	Green
1/2 Round Body, Taper Point 6mm x 3mm x 3mm Teflon Pledget	 		PS2706	4/0 1,5	75	12	Green

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
25 mm					
1/2 Round Body, Taper Point 	PS3727	2/0 3,0	75	12	Green
	PS4191	0 3,5	75	12	Green
1/2 Round Body, Taper Point 	PS3258	3/0 2,0	75	12	Green
	PS3262	3/0 2,0	90	12	Green
	PS3730	2/0 3,0	75	6	5 Green 5 White
	PS3733	2/0 3,0	75	12	Green
	PS3737	2/0 3,0	90	6	5 Green 5 White
	PS3741	2/0 3,0	90	12	Green

25 mm					
1/2 Round Body, Taper Point  3mm x 3mm x 1.5mm Teflon Pledget 	PS3738	2/0 3,0	90	12	Green
1/2 Round Body, Taper Point  6mm x 3mm x 1.5mm Teflon Pledget 	PS3263	3/0 2,0	90	12	Green
	PS3729	2/0 3,0	75	6	5 Green 5 White
	PS3732	2/0 3,0	75	12	Green
1/2 Tapercutting  3mm x 3mm x 1.5mm Teflon Pledget 	PS5826	2/0 3,0	90	6	5 Green 5 White

26 mm					
1/2 Round Body, Taper Point 	PS3788	2/0 3,0	75	12	Green
	PS4217	0 3,5	75	12	Green
1/2 Round Body, Taper Point 	PS3322	3/0 2,0	90	12	Green
	PS4221	0 3,5	90	12	Green




Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
1/2 Round Body,Taper Point		PS3370	3/0 2,0	75	12	Green
		PS3860	2/0 3,0	75	12	Green
		PS4247	0 3,5	75	12	Green
		PS4611	1 4,0	75	12	Green
1/2 Reverse Cutting		PS4268	0 3,5	75	12	Green
		PS4623	1 4,0	75	12	Green
3/8 Reverse Cutting		PS3400	3/0 2,0	75	12	Green
		PS3892	2/0 3,0	75	12	Green
35 mm						
3/8 Reverse Cutting		PS3935	2/0 3,0	75	12	Green
36 mm						
1/2 Round Body,Taper Point		PS5515	1 4,0	75	12	Green
1/2 Reverse Cutting		PS5514	0 3,5	75	12	Green
40 mm						
1/2 Round Body,Taper Point		PS4369	0 3,5	75	12	Green
		PS4714	1 4,0	75	12	Green
		PS5007	2 5,0	90	12	Green
45 mm						
1/2 Round Body,Taper Point		PS5036	2 5,0	75	12	Green

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
50 mm					
1/2 Round Body, Taper Point 	PS4814	1 4,0	75	12	Green
	PS5067	2 5,0	75	12	Green
Straight Reverse Cutting 	PS4049	2/0 3,0	75	12	Green
55 mm					
1/2 Round Body, Taper Point 	PS5181	5 7,0	75	12	Green
1/2 Taper cutting 	PS5173	5 7,0	75	12	Green

Polisil

Surgical Suture, Braided Polyester, Silicone Coated (PET) Nonabsorbable, Synthetic

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
17x45 cm					
POLISIL Non-Needle Suture	PS2965	3/0 2,0	17x45	12	Green
	PS3559	2/0 3,0	17x45	12	Green
	PS4150	0 3,5	17x45	12	Green
	PS4560	1 4,0	17x45	12	Green
13x60 cm					
POLISIL Non-Needle Suture	PS2962	3/0 2,0	13x60	12	Green
	PS3554	2/0 3,0	13x60	12	Green
250 cm					
POLISIL Non-Needle Suture	PS4543	1 4,0	250	12	Green