


Polytef

Surgical Suture, Braided Polyester, PTFE Coated Nonabsorbable, Synthetic

Surgical Sutures
Non-Absorbable

Boz®



Definition

Synthetic Non-Absorbable

Material

Polyester Polyethylene Terephthalate (PET)

Coating

PTFE (polytetrafluoroethylene)

Structure

Braided

Colour

Green

Size

EP : 0.5 to 3

USP : 7/0 to 2/0

Tensile Strength Retention

Permanent

Mass Absorption

Permanent

Sterilization Method

Ethylene Oxide Gas

Main Indication

General Soft tissue approximation or ligation

Cardiovascular, Vascular Surgery


Ophthalmic Surgery

Neurosurgical Surgery


Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
10 mm							
3/8 Reverse Cutting Micropoint 280 Micron			PT5557	7/0 0,5	55	12	Green
3/8 Reverse Cutting Micropoint 330 Micron			PT5556	6/0 0,7	55	12	Green

12 mm							
3/8 Reverse Cutting Micropoint 380 Micron			PT5555	5/0 1,0	55	12	Green

16 mm							
3/8 Reverse Cutting Micropoint			PT5554	5/0 1,0	75	12	Green
3/8 Reverse Cutting Micropoint			PT5553	4/0 1,5	75	12	Green

17 mm							
1/2 Trocar			PT5552	4/0 1,5	75	12	Green

18 mm							
3/8 Reverse Cutting			PT5551	3/0 2,0	75	12	Green

20 mm							
1/2 Trocar Point			PT5550	3/0 2,0	75	12	Green
3/8 Reverse Cutting			PT5565	2/0 3,0	75	12	Green