


Tempowire

Surgical Stainless Steel Wire Nonabsorbable, Multifilament



Boz®

Definition

Non-Absorbable

Material

316L Stainless Steel conforming to
(ISO 5832 Part-1)

Coating

FEP, Polyethylene

Structure

Braided, Twisted

Colour

Blue, Transparent, Orange

Size

EP : 2 to 3.5

USP : 3/0 to 0

Tensile Strenght Retention

Permanent

Mass Absorption

Permanent

Sterilization Method










Ethylene Oxide Gas










Main Indication

Cardiac Surgery

Tempowire

Surgical Stainless Steel Wire Nonabsorbable, Multifilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
17 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle	TW3068	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5918	3/0 2,0	60	12	PE-Blue
22 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Zig - Zag type	TW3194	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5874	3/0 2,0	60	12	PE-Blue
25 mm					
1/2 Round Body, Taper Point  60 mm Straight Reverse Cutting   Breakable needle  Anchor	TW5869	2/0 3,0	60	12	PE-Blue
26 mm					
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   Breakable needle	TW3327	3/0 2,0	60	12	PE-Blue

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm					
1/2 Round Body, Taper Point  85 mm Straight Reverse Cutting   	TW5877	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW3326	3/0 2,0	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW4208	0 3,5	60	12	PE-Blue
1/2 Round Body, Taper Point  63 mm Straight Reverse Cutting   	TW5576	3/0 2,0	60	12	PE-Blue
1/2 Tapercutting  63 mm Straight Reverse Cutting  	TW3296	3/0 2,0	60	12	PE-Blue
1/2 Tapercutting  63 mm Straight Reverse Cutting   	TW4204	0 3,5	60	12	FEP
1/2 Tapercutting  63 mm Straight Reverse Cutting   	TW3297	3/0 2,0	60	12	PE-Blue
3/8 Round Body, Taper Point  63 mm Straight Reverse Cutting 	TW5907	3/0 2,0	60	12	PE-Blue