

TEFLON PLEDGET



Surgical, Polytetrafluoroethylene (PTFE) Felt
Nonabsorbable, Non-woven, White

DESCRIPTION

Boz Pledgets are precut pieces from a sheet of polytetrafluoroethylene (PTFE) fabric that are used to spread the load of the suture against host tissue.

INDICATIONS

Boz Pledgets are indicated as sterile, single use, implantable devices placed in the heart during septal defect repairs, for patch grafting, to repair tissue, and to buttress sutures in cardiovascular surgery.



Code	Dimension	Description
TPFE01	6x3x1.5mm	6 pcs teflon pledget (felt) per pouch
TPFE02	6x3x1.5mm	5 pcs teflon pledget (felt) per pouch
TPFE03	6x3x1.5mm	10 pcs teflon pledget (felt) per pouch
TPFE04	7x3x1.5mm	6 pcs teflon pledget (felt) per pouch
TPFE10	3x3x1.5mm	10 pcs teflon pledget (felt) per pouch
TPFE11	3x3x1.5mm	5 pcs teflon pledget (felt) per pouch
TPFE18	5x5x1,5 mm	10 pcs teflon pledget (felt) per pouch
TPFE19	5x10x1,5 mm	10 pcs teflon pledget (felt) per pouch
TPFE20	6x4x1.5mm	10 pcs teflon pledget (felt) per pouch
TPFE21	6x5x1,5 mm	10 pcs teflon pledget (felt) per pouch
TPFE28	6x8x1,5mm	4 pcs teflon pledget (felt) per pouch
TPFE30	6x10x1,5mm	5 pcs teflon pledget (felt) per pouch
TPFE50	6x6x1,5mm	10 pcs teflon pledget (felt) per pouch
TPFE60	7x3x1,5mm	10 pcs teflon pledget (felt) per pouch
TPFE61	13x8x1,5 mm	10 pcs teflon pledget (felt) per pouch
TPFE62	100x100x1,5 mm	1 pc teflon pledget (felt) per pouch
TPFE63	100x150x1,5 mm	1 pc teflon pledget (felt) per pouch
TPFE64	150x150x1,5 mm	1 pc teflon pledget (felt) per pouch



Sağlık Mahallesi Sağlık 1 Sk. No:33/5 Çankaya/ANKARA/TÜRKİYE
Tel: +90 (312) 254 03 40 (5 hat) - Faks: +90 (312) 254 03 50
web:www.boztibbi.com - email:boz@boztibbi.com

